

HIGH-FLOW Traditional vaporizers



VetFlo™

Conventional anesthesia system with single and multiple anesthesia systems

Kent Scientific's small animal and rodent anesthesia vaporizer is one of the easiest to use and most economical gas anesthetic systems on the market, with wide laboratory application including for in-vivo COVID-19 research.

Perfect for laboratory and veterinary anesthesia delivery, the VetFlo Anesthesia System is completely self-contained and portable.

All systems include:

- Tabletop stand
- Flowmeter
- High pressure oxygen hose · WAG canisters (2)
- One Year Warranty

3 reasons to choose VetFlo

- 1 **New systems** — never reconditioned or re-manufactured
- 2 **Completely self-contained** and portable — ideal for flexible lab setups
- 3 **Simple set-up** and operation



VetFlo systems are available on pole mounted stands for easy mobility.

Starter Kit includes all the accessories needed to begin your procedure.

More information

Toll-Free Phone: 888-5RATTUS (888-572-8887)

Local Phone: 860-626-1172

Fax: 860-626-1179

E-mail: sales@kentscientific.com

kentscientific.com



Providing small animal researchers with over 30 years of innovation.

Anesthesia Products



LOW-FLOW
SomnoFlo®
SomnoFlo® O₂Care
SomnoSuite®

HIGH-FLOW
VetFlo™

LOW-FLOW Electronic vaporizers

Ambient air, less WAG, no calibration.

No more using human amounts of anesthetic for a mouse!

Unlike Traditional Vaporizers, our Low-Flow Vaporizers are engineered with solid state electronics and an integrated electronic vaporizer which uses either ambient air or compressed gas to deliver isoflurane or sevoflurane anesthetic at low flow rates proportionate to the animal's size. This advanced design eliminates the need for annual certification that is required by traditional wick-design vaporizers.

Using less anesthesia benefits the animal during procedures, and significantly reduces exposure of lab personnel to WAG while generating impressive cost savings. Kent Scientific Low-Flow Electronic Vaporizers are used in the world's leading pharmaceutical and academic laboratories to further advance training and research in a host of fields, including COVID-19 related research.

Features & benefits of low-flow

- **Low flow rates** – Saves money by using less than 1 mL/hr of isoflurane
- **Built-in air compressor** – Uses room air or compressed gas
- **Automated anesthetic delivery system** – Ensures precision to 0.1% delivery of anesthetic
- **No servicing or calibration needed** – Cost-effective, reliable equipment
- **Induction chamber purge** – Minimizes WAG exposure

Estimated yearly operational savings

Low use
100 hours/yr

\$3,320

Medium use
250 hours/yr

\$8,160

High use
500 hours/yr

\$16,920

	SomnoSuite	SomnoFlo O ₂ Care*	VetFlo
Flow rate			
Range	25 to 1,000 mL/min	100 to 1,000 mL/min	0.1L/min, 2-4 L/min
Setting	Animal weight (with flow rate displayed)	Flow rate (with animal weight displayed)	Flow rate (no animal weight displayed)
Manual override	Yes	Yes	N/A
Sensor	Electronic flow sensor (monitors and controls anesthetic delivery)	Electronic flow sensor (monitors and controls anesthetic delivery)	None
Operation			
Modes	Ambient air and compressed gas	Ambient air and compressed gas	Compressed gas
Anesthetic	Isoflurane and sevoflurane	Isoflurane and sevoflurane	Isoflurane and sevoflurane
Presets	User-selectable	User-selectable	None
Alarms	Yes	Yes	No
Purge	Yes	Yes	No
Maintenance	No annual calibration	No annual calibration	Annual calibration required
Application			
Multi-animal	Yes	Yes	Yes
Stereotaxic	Yes	Yes	Yes
Filling and delivery			
Method	Glass syringe	Direct from bottle	Funnel or key fill
Capacity	5 mL and 10 mL	100 mL and 250 mL	250 mL
Accuracy	0.1%	0.1%	0.5%
Data			
Output	Yes	No	No
Internal storage	Yes	No	No
Optional modules			
Available	Yes (RightTemp®, MouseSTAT®, RoVent®, CapnoScan®)	No	No
Warranty			
Standard (included)	1-year	1-year	1-year (vaporizer)
Premium (optional)	3-year	3-year	

*Also available: SomnoFlo for everyday low-flow anesthesia applications. Ask your Kent representative which system is the best fit for your protocol.

SomnoFlo® O₂Care

Recent research has highlighted that while pure oxygen can cause unintended off-target effects during anesthesia, 21% oxygen isn't always sufficient to prevent hypoxia. Custom gas blends and external gas mixers, while effective, can often be cumbersome, expensive, and challenging to use.

SomnoFlo O₂Care offers an innovative approach by blending ambient air with 100% compressed oxygen to deliver the ideal anesthetic carrier gas composition for longer procedures.

Need a routine low-flow option? SomnoFlo® is also available for everyday small animal anesthesia workflows.

“Our SomnoFlo O₂Care set up was easy and the digital functions of the vaporizer are a huge plus.”

—Lab Director

“SomnoFlo O₂Care menus are easy to use & understand. I very much like that there are only 2–3 buttons necessary for general use.”

—University Post Doc



SomnoSuite®

The SomnoSuite is designed with a precision syringe pump that delivers precise amounts of anesthesia. Like the SomnoFlo O₂Care, its integrated electronic vaporizer uses either ambient air or compressed gas to deliver precise anesthesia according to the animal's size.

SomnoSuite allows for the addition of modular physiological monitors within its compact design. All SomnoSuites are equipped with a RightTemp® temperature monitoring and homeothermic control module.

SomnoSuite modules

- **MouseSTAT®** Pulse oximeter & heart rate monitor
- **Rovent®** Small animal ventilator
- **CapnoScan®** End-tidal CO₂ monitor
- **RightTemp®** Temperature monitor & homeothermic control

“The SomnoSuite system is phenomenal ... It is one of the best investments I've made for the lab!”

—University Professor

“We love your SomnoSuite system! It's so convenient, so safe and also saves isoflurane!”

—Lab Manager