

TBS-380 Fluorometer

This is a Quick Start Guide to assist in setting up the TBS-380 Fluorometer to run samples and retrieve data as quickly as possible. For more detailed information, see the TBS-380 User's Manual or visit our website at www.turnerbiosystems.com.

1 Connections

1. Assemble the power supply by sliding the appropriate AC clip onto the AC adaptor, and then plug into the wall power outlet.
2. Connect the power supply to the TBS-380.
3. Connect the 9 pin serial cable between the TBS-380 and the computer/printer. The female connector attaches to the TBS-380.
4. Turn on the TBS-380.



2 Calibration

1. Select the right channel by pressing the **A/B** button.
2. Select a standard value from 1 to 999 by pressing the **STD VAL** button and then the **UP** or **DOWN** arrows.
3. Press **ENTER** when finished.
4. Press the **CAL** button and then **ENTER**.
5. Insert a cuvette with a blank solution and press **ENTER**.
6. Insert a cuvette with the calibration standard and press **ENTER**.
7. Press **ENTER** to accept the calibration.

NOTE: If the Minicell Adaptor is used, be sure to insert it into the right position in the sample compartment (see reverse).



3 Sample Analysis

You can read samples in two reading modes: Discrete or Continuous.

Discrete Mode

1. Press the **SET** button.
2. Select "Discrete" using the up or down arrows and press **ENTER**.
3. Insert the sample and press **READ**.

NOTE: For low volume samples use the Minicell Adaptor (see reverse).

Continuous Mode

1. Press the **SET** button.
2. Select "Continuous" using the **UP** or **DOWN** arrows and press **ENTER**.
3. Select a measure interval using the **UP** or **DOWN** arrows and press **ENTER**.
4. Select a number of measurements you want to obtain using the **UP** or **DOWN** arrows and press **ENTER**.
5. Press the **READ** button and the screen will display the setting selected.
6. Press **ENTER** to accept the settings.
7. Insert the sample and press **ENTER**.

4 Installation of Spreadsheet Interface Software

1. Exit all Windows programs.
2. Insert the CD and run the setup program.
3. The setup wizard will install the necessary files.
4. When the setup is complete, an icon named "Spreadsheet Interface Software" will be found in the "Programs" menu.
5. Restart your computer.
6. Open the software and select the right COM port.
7. Click on "Start"



TBS-380

QUICK VIEW DIAGRAM

(See reverse for Quick Start Instructions)



Positioning the Minicell Adaptor for Sample Analysis/ Calibration

Minicell Adaptor in BLUE Position



When "BLUE" channel is selected, insert the Minicell Adaptor with the BLUE label facing you.

Minicell Adaptor in UV Position



When "UV" channel is selected, insert the Minicell Adaptor with the UV label facing you.



TURNER BIOSYSTEMS

645 N Mary Ave, Sunnyvale CA 94085 USA
Phone: 408.636.2400 - Fax: 408.737.7919 - Toll Free: 888.636.2401
www.turnerbiosystems.com