Spectral Cell Analyzer

Safety Guide

Safety Guide

EN

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JP

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FR

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DE

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IT

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ES

US/International

Sony Biotechnology Inc.
1730 North 1st St, 2NW
San Jose, CA, 95112, U.S.A
Voice: +1 800-275-5963
FAX: +1 408-352-4130
sales@sonybiotechnology.com

Japan

Sony Corporation
Medical Business Unit
Life Science Business Division
1-7-1 Konan, Minato-ku,
Tokyo 108-0075, Japan
Voice: 0120-667-010
FAX: 0120-388-060
cytometry@sony.co.jp

Europe

Sony Biotechnology Europe
The Heights, Brooklands,
Weybridge, Surrey, KT13 0XW, UK
Voice: +44 (0) 1932 816000
FAX: +44 (0) 1932 817170
Sales_EU@sonybiotechnology.com

Before operating the instrument, please read this manual thoroughly and retain it for future reference.

WARNING

To reduce the risk of fire or electric shock, do not expose this apparatus to rain or moisture.

To avoid electrical shock, do not open the cabinet. Refer servicing to qualified personnel only.

THIS APPARATUS MUST BE EARTHED.

This symbol is intended to alert the user to the presence of important operating and maintenance (servicing) instructions in the literature accompanying the appliance.

Important Safety Instructions

If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

- Read these instructions.
- Keep these instructions.
- Heed all warnings.
- Follow all instructions.
- When installing, place the unit on a sturdy, level surface.
- Do not block any ventilation openings.
Install in accordance with the manufacturer's instructions.
- Do not install near any heat sources such as radiators, heat registers, stoves, or other apparatus that produce heat.
- Do not defeat the safety purpose of the polarized or grounding-type plug. A polarized plug has two blades with one wider than the other. A grounding-type plug has two blades and a third grounding prong. The wide blade or the third prong are provided for your safety. If the provided plug does not fit into your outlet, consult an electrician for replacement of the obsolete outlet.
- Protect the power cord from being walked on or pinched particularly at plugs, convenience receptacles, and the point where they exit from the apparatus.
- Only use attachments/accessories specified by the manufacturer.
- Unplug this apparatus during lightning storms or when unused for long periods of time.
- Refer all servicing to qualified service personnel. Servicing is required when the apparatus has been damaged in any way, such as power-supply cord or plug is damaged, objects have fallen into the apparatus, does not operate normally, or has been dropped.

When installing the installation space must be secured in consideration of the ventilation and service operation.

- Do not block the ventilation slots at the rear panels, and vents of the fans.
- Leave a space around the unit for ventilation.
- Leave more than 15 cm (6 in.) of space in the rear of the unit to secure the operation area.


When the unit is installed on a bench, leave clearance of at least 40 cm (16 in.) on the left and right, and 15 cm (6 in.) at the front. It is recommended that 40 cm (16 in.) clearance be left above the unit for performance and maintenance considerations.

CAUTION

The supplied power cord is for use with this device only. Do not use the cord with other devices.

Caution

Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous radiation exposure.

This LE-SP6800 Spectral Cell Analyzer is classified as a CLASS 1 LASER PRODUCT. The CLASS 1 LASER PRODUCT label is located on the rear panel.

Laser Properties

Excitation Laser Module

Wavelength : 488 nm
Emission duration : Continuous
Laser output power : 80 mW (max)

Wavelength : 638 nm
Emission duration : Continuous
Laser output power : 110 mW (max)

Wavelength : 405 nm
Emission duration : Continuous
Laser output power : 140 mW (max)

Standard : IEC60825-1 (2007)

Caution

The use of optical instruments with this product will increase eye hazard.

Regulatory Information

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

For the customers in the U.S.A

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules.

These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment.

This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications.

Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense. You are cautioned that any changes or modifications not expressly approved in this manual could void your authority to operate this equipment.

For the customers in Canada

CAN ICES-3 (A)/NMB-3(A)

For the customers in Europe

This product has been manufactured by or on behalf of Sony Corporation, 1-7-1 Konan Minato-ku Tokyo, 108-0075 Japan. Inquiries related to product compliance based on European Union legislation shall be addressed to the authorized representative, Sony Deutschland GmbH, Hedelfinger Strasse 61, 70327 Stuttgart, Germany. For any service or guarantee matters, please refer to the addresses provided in the separate service or guarantee documents. This product may cause interference if used in residential areas. Such use must be avoided unless the user takes special measures to reduce electromagnetic emission to prevent interference to the reception of radio and television broadcasts.

The A-weighted emission sound pressure level at workstations does not exceed 70 dB(A).

For the customers in the U.S.A.

SONY LIMITED WARRANTY - Please visit <http://www.sony.com/psa/warranty> for important information and complete terms and conditions of Sony's limited warranty applicable to this product.

For the customers in Canada

SONY LIMITED WARRANTY - Please visit <http://www.sonybiz.ca/pro/lang/en/ca/article/resources-warranty-product-registration> for important information and complete terms and conditions of Sony's limited warranty applicable to this product.

For the customers in Europe

Sony Professional Solutions Europe - Standard Warranty and Exceptions on Standard Warranty.
Please visit <http://www.pro.sony.eu/warranty> for important information and complete terms and conditions.

For the customers in Korea

SONY LIMITED WARRANTY - Please visit <http://bpeng.sony.co.kr/handler/BPAS-Start> for important information and complete terms and conditions of Sony's limited warranty applicable to this product.

Safety and Warning Labels

This symbol is intended to alert the user to the presence of risk of exposure to biologically transmissible disease.


This symbol is intended to alert the user to the presence of important operating and maintenance (servicing) instructions in the literature accompanying the appliance.


This symbol indicates the protective conductor terminal. Always connect the earth connector.

Main Unit (Front)

- Fluids handled by this unit may contain biological, chemical, or other agents. Always wear gloves, and also wear face mask, goggles, and other protective clothing as required.
- The sample tube holder ② may pinch your fingers if your fingers are in close proximity when the holder is moving. Keep your hands clear at all times, except when loading and unloading the sample tube.
- Touching the safety button ① activates the interlock system, stopping and then lowering the sample tube holder ②.

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
CERTIFICATION
 PRODUCT COMPLIES WITH 21
 CFR 1040. 10 AND 1040. 11
 EXCEPT FOR DEVIATIONS
 PURSUANT TO LASER NOTICE
 No. 50, DATED JUNE 24, 2007
 Sony Corporation
 1-7-1 KONAN, MINATO-KU,
 TOKYO, 108-0075 JAPAN
 MANUFACTURED :
 SKZ
 4-444-974-01

CLASS 1 LASER PRODUCT
LASER KLASSE 1
LUOKAN 1 LASERLAITE
KLASS 1 LASERAPPARAT
 4-444-974-01


- Waste fluid may contain biological, chemical, or other agents. Always wear gloves and other protective clothing, mask, and goggles, as required, when handling waste fluid. The handling of waste fluid should be performed in accordance with biological hazard handling safety procedures.
- Waste fluid may contain virulent or hazardous substances that may be highly contagious or toxic, depending on the sample under test. Always wear gloves and other protective clothing, mask, and goggles, as required, when handling waste fluid.
- Do not touch the surfaces of the waste tank cap and tank nozzle that come into contact with waste fluid. Take care not to spill any fluid.
- When disposing of waste fluid, always follow laboratory rules and procedures for the decontamination to neutralize all biological hazards.
- Make sure that the waste fluid line is securely connected to the connector on the rear panel of main unit.
- Make sure the compressed air supply pressure is in the range 350 kPa to 450 kPa (51 psi to 65 psi).

Main Unit (Front, Internal)

- Sample fluid may contain biological, chemical, or other agents. Always wear gloves and face mask when changing the sample tube.
- Sample fluid may contain biological, chemical, or other agents. Always wear gloves and face mask when changing the flow cell chip.

Main Unit (Right Side, Internal)


CAUTION - CLASS 3B LASER
RADIATION WHEN OPEN
AVOID EXPOSURE TO THE BEAM

4-430-744-01


Class 3B laser radiation may be emitted if the cover is removed. Laser radiation can cause serious injury. Avoid eye and skin exposure to laser beams.

Main Unit (Left Side, Internal)


- Waste fluid may contain biological, chemical, or other agents. Always wear gloves and other protective clothing, mask, and goggles, as required, when emptying or replacing the waste tank. The handling of waste fluid should be performed in accordance with biological hazard handling safety procedures.
- Waste fluid may contain virulent or hazardous substances that may be highly contagious or toxic, depending on the sample under test. Always wear gloves and other protective clothing, mask, and goggles, as required, when handling waste fluid.
- When handling the waste tank, always follow the laboratory rules for the handling of waste fluids.
- Do not touch the surfaces of the waste tank cap and tank nozzle that come into contact with waste fluid. Take care not to spill any fluid.
- Make sure that the cap on the waste tank is correctly and securely fastened at all times.
- Make sure that the waste fluid line is securely connected to the connector on the waste tank.

Fluidics Cart (Rear)

- Waste fluid may contain biological, chemical, or other agents. Always wear gloves and other protective clothing, mask, and goggles, as required, when handling waste fluid. The handling of waste fluid should be performed in accordance with biological hazard handling safety procedures.
- Waste fluid may contain virulent or hazardous substances that may be highly contagious or toxic, depending on the sample under test. Always wear gloves and other protective clothing, mask, and goggles, as required, when handling waste fluid.
- Do not touch the surfaces of the waste tank cap and tank nozzle that come into contact with waste fluid. Take care not to spill any fluid.
- Make sure that the waste fluid line is securely connected to the connector on the rear panel of the fluidics cart.

Waste Tank (Fluidics Cart)

- Waste fluid may contain biological, chemical, or other agents. Always wear gloves and other protective clothing, mask, and goggles, as required, when emptying or replacing the waste tank. The handling of waste fluid should be performed in accordance with biological hazard handling safety procedures.
- Waste fluid may contain virulent or hazardous substances that may be highly contagious or toxic, depending on the sample under test. Always wear gloves and other protective clothing, mask, and goggles, as required, when handling waste fluid.
- When handling the waste tank, always follow the laboratory rules for the handling of waste fluids.
- Do not touch the surfaces of the waste tank cap and tank nozzle that come into contact with waste fluid. Take care not to spill any fluid.
- Make sure that the cap on the waste tank is correctly and securely fastened at all times.
- Make sure that the waste fluid line is securely connected to the connector on the waste tank.

Usage Notes

Use and Storage

Do not subject the unit to severe shocks

The internal mechanism may be damaged or the body warped.

Do not cover the unit while operating

Doing so will cause temperatures to rise inside the unit, possibly resulting in failure.

Do not place objects on top of the unit

Objects placed on top of the unit may fall, possibly causing an accident or injury.

Use only under normal operating conditions

Always use the unit under specified normal operating conditions to ensure correct performance.

Use only genuine parts

Always use only genuine cables and other parts to ensure optimum performance.

Do not relocate the unit

The main unit is quite heavy, and relocating it is fraught with hazard. Do not attempt to relocate the unit by yourself.

Contact your Sony distributor if you wish to relocate the unit.

Before use

Always perform the following before using the unit.

- Check that all tanks are installed and securely connected.
- Check the fluid levels in the tanks.
- Check that the sheath filter is correctly installed.

If not using the unit for an extended period

Perform the following.

- Clean the unit according to the maintenance procedures.
- After shutting down, open the fluidics maintenance door on the left side and turn off the main power switch, and remove the power cord from the power outlet.
- Dispose of any waste fluid in the waste tanks in accordance with laboratory rules and guidelines.
- Remove all sheath fluid from the sheath tank.
- Remove all DI water from the DI water tank.
- Close all doors.
- Store the unit on a level, stable surface.

Transporting the unit

Consult your Sony distributor.

Moving the unit

The unit is a high-precision piece of equipment. Moving the unit can subject the unit to vibration or shock, and may cause the unit to malfunction. If you wish to move the unit, contact your Sony distributor to move and recalibrate the unit.

Care of the unit

Shut down the unit and clean it with a clean, soft, lint-free cloth. If the body of the unit is dirty, spray the cloth with ethanol from a spray bottle. Always wear gloves when cleaning to prevent bacterial contamination.

In the event of operating problems

Consult your Sony distributor.

Regular inspection

It is recommended that your Sony distributor regularly inspect the unit once a year to ensure continued performance.

Product disposal

When disposing of the product, always follow laboratory rules and procedures for the decontamination of the product to neutralize all biological hazards.

Disposal of consumables

Always follow the laboratory rules and guidelines when disposing of used tubes, flow cell chips, gloves, and other consumables.

Data backup

All instrument settings and measurement data are stored on the hard disk in the host computer, and that data may become corrupted or lost if any unforeseen troubles occur. It is recommended that all data be backed up regularly.

SP6800 software

- Microsoft SQL Server 2012 Express is installed on the host computer together with SP6800 software. Refrain from upgrading SQL Server, rewriting the database, or performing other actions. Doing so may cause unexpected problems.
- Refrain from launching other applications that consume significant CPU resources or memory while running SP6800 software. Doing so may cause unexpected problems.

Specifications

General

Configuration	2-laser type (488/638)	2-laser type (405/488)	3-laser type (405/488/638)
Model	LE-SP6800ZB	LE-SP6800ZC	LE-SP6800ZE
Dimensions	Main unit: 60.0 cm × 63.5 cm × 71.3 cm (23 ⁵ / ₈ in. × 25 in. × 28 ¹ / ₈ in.) (w/d/h) Fluidics cart: 78.6 cm × 52.1 cm × 58.0 cm (31 in. × 20 ⁵ / ₈ in. × 22 ⁷ / ₈ in.) (w/d/h)		
Mass	Main unit: Approx. 99 kg (218 lb. 4 oz.) (dry weight), Fluidics cart: Approx. 32 kg (70 lb. 8 oz.) (dry weight)		
Compressed air supply pressure	350 kPa to 450 kPa (51 psi to 65 psi)		
Power supply	100 V to 240 V AC, 50/60 Hz		
Power consumption	220 W (maximum)		
Operating temperature	16 °C to 29 °C (61 °F to 84 °F), room temperature variation: ≤ 5 °C (9 °F)		
Operating humidity	20% to 80% (condensation free)		

Optics System

Configuration		2-laser type (488/638)	2-laser type (405/488)	3-laser type (405/488/638)
Model		LE-SP6800ZB	LE-SP6800ZC	LE-SP6800ZE
Laser specifications				
	Wavelength	488 nm, 638 nm semiconductor lasers	405 nm, 488 nm semiconductor lasers	405 nm, 488 nm, 638 nm semiconductor lasers
	Maximum output (irradiation power at flow cell chip)	488 nm: 40 mW 638 nm: 60 mW	405 nm: 60 mW 488 nm: 40 mW	405 nm: 60 mW 488 nm: 40 mW 638 nm: 60 mW
	Irradiation method	2-spot dual-axis irradiation		
Detection system				
	Scatter measurement	Forward scatter (FSC), Side scatter (SSC)		
	Spectroscopic method	Prism array spectroscopy		
	Fluorescence measurement	32-channel PMT (500 nm to 800 nm detection wavelength)	32-channel PMT (500 nm to 800 nm detection wavelength) PMT × 2 (420 nm to 440 nm, 450 nm to 470 nm)	
	Cell passage position measurement	4-segment photodiode (4 optical sensors)		

Measurement system

Configuration	2-laser type (488/638)	2-laser type (405/488)	3-laser type (405/488/638)
Model	LE-SP6800ZB	LE-SP6800ZC	LE-SP6800ZE
Signal resolution	Height: 20 bits, Area: 32 bits (equivalent to 4.29 million channels), Sampling frequency: 50 MHz		
Measurement parameters	64 color parameters, FSC, SSC, cell passage position coordinates (X,Y)	66 color parameters, FSC, SSC, cell passage position coordinates (X, Y)	
Pulse processing	Height, Area, Width (all measured)		
Optical axis adjustment method	Automatic adjustment using CoreFinder™ technology		
Event rate	10,000 eps (nominal), 20,000 eps (maximum)		

Fluidics system

Configuration	2-laser type (488/638)	2-laser type (405/488)	3-laser type (405/488/638)
Model	LE-SP6800ZB	LE-SP6800ZC	LE-SP6800ZE
Sample loader type	Single auto-loader		
Supported sample tube	5.0 ml polystyrene round tubes		
Flow rate	Low (3 m/s), Mid (5 m/s), High (10 m/s)		
Optical detection channel size	200 μm $\times$ 200 μm		

Performance

Configuration	2-laser type (488/638)	2-laser type (405/488)	3-laser type (405/488/638)
Model	LE-SP6800ZB	LE-SP6800ZC	LE-SP6800ZE
Fluorescence sensitivity	FITC: 120 MESF, PE: 70 MESF		
Fluorescence detection linearity	FITC $R^2 \geq 0.995$, PE $R^2 \geq 0.995$		
Detectable cell size	0.5 μm to 40 μm (beads)		
Fluorescence detection resolution	488 nm laser: CH16, 638 nm laser: CH27 ¹⁾ $\leq 2.5\%$ (HPCV)		

1) Excluding the LE-SP6800ZC.

Design and specifications are subject to change without notice.

Note

Always verify that the unit is operating properly before use. SONY WILL NOT BE LIABLE FOR DAMAGES OF ANY KIND INCLUDING, BUT NOT LIMITED TO, COMPENSATION OR REIMBURSEMENT ON ACCOUNT OF THE LOSS OF PRESENT OR PROSPECTIVE PROFITS DUE TO FAILURE OF THIS UNIT, EITHER DURING THE WARRANTY PERIOD OR AFTER EXPIRATION OF THE WARRANTY, OR FOR ANY OTHER REASON WHATSOEVER.

Environmental Notices

Disposal of Old Electrical & Electronic Equipment (Applicable in the European Union and other European countries with separate collection systems)

This symbol on the product or on its packaging indicates that this product shall not be treated as household waste. Instead it shall be handed over to the applicable collection point for the recycling of electrical and electronic equipment. By ensuring this product is disposed of correctly, you will help prevent potential negative consequences for the environment and human health, which could otherwise be caused by inappropriate waste handling of this product. The recycling of materials will help to conserve natural resources. For more detailed information about recycling of this product, please contact your local Civic Office, your household waste disposal service or the shop where you purchased the product.

Trademarks

- Microsoft and SQL Server are registered trademarks of Microsoft Corporation in the United States and/or other countries.

The products or system names appearing in this document are trademarks or registered trademarks of their respective owners. Further, the ® or ™ symbols are not used in the text.

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