

Set solvents

NOTE

The selection of the solvents itself must be done in the standard pump method user interface.

- Open the pump method dialog using the button **Advanced 2D pump settings...** and change the selection of the solvents there.
- After closing the dialog, the solvent settings should be updated immediately.

- 1 Set the percentage of solvent B to any value from 0 – 100 % in steps of 0.01 %.

Solvents

Percentage settings

A: %

B: %

Solvent information

A1: 100.0 % Water V.03

B1: 100.0 % Acetonitrile V.03

Figure 64 2D-LC solvent settings


Solvent A always delivers the remaining percentage of volume. If the rate of solvent B is for example set to 20 %, solvent A, following the calculation $\%A = 100 - \%B$, automatically is set to 80 %.

The name of the selected solvents and their solvent channels (**A1:...** or **A2:...** and **B1:...** or **B2:...**) are shown in the corresponding text fields.

NOTE

The corresponding Percent B value in the Standard Pump user interface will be ignored as long as the 2D-LC functionality is enabled (see [“Configuration”](#) on page 148).

Set flow

A screenshot of a software dialog box titled "Flow settings". It contains two rows of controls. The first row has the label "2D Flow" followed by a text input field containing "1.00" and the unit "ml/min". The second row has a checked checkbox labeled "use idle flow" followed by a text input field containing "1.00" and the unit "ml/min".

Flow settings		
2D Flow	1.00	ml/min
☒ use idle flow	1.00	ml/min

Figure 65 Flow settings

- 1 Set the **2D Flow** (range 0 – 5.0 mL/min).

This defines the flow in the 2nd dimension being used while 2D-LC is active (within 2D time segments where mode is not equal to OFF)

- 2 To set and use **Idle Flow** select check box **use idle flow**.

This defines the flow in the 2nd dimension that is used while the 2D-LC mode is OFF (range 0 – 5.0 mL/min).

NOTE

If **use idle flow** is not selected, the **2D Flow** is also used while 2D-LC mode is OFF.

Set Solvent Composition Gradient

Set Solvent Composition Gradient

The timetable in the **²D Gradient** group allows changing the solvent composition.

Percent B ranges from 0 – 100 %.

Change the solvent composition at a specified time

- 1 To change the solvent composition (%B) at the specified time apply a percent B range from 0 – 100 %

NOTE

Different start conditions in the first row may cause step gradients and RI-effects (density differences of the different liquid phases may cause different DAD detection through baseline disturbances).

²D gradient		
	Time [min]	% B
▶ ⊕	0.00	10.00
⊕	1.25	60.00

The time axis relates to the Stoptime of the 2nd dimension pump. **Time [min] = 0.00** marks the start of the maybe repetitive gradient cycles, a time greater than **Stoptime ²D** will be ignored.

Define shifted gradients

- 1 To modify an entry in the timetable over the runtime of the 1st dimension (shifted gradient), click **+**-sign at the beginning of the line and add one or more lines.

2D Gradient

Time [min]	% B
0.00	20.00
0.10	50.00

+ - ✕ ◀ ▶

NOTE

A gray colored **+** indicates that the associated nested table has no entries, otherwise the **+** is black.

- To add a new entry, use the add button (**+** sign at the table bottom) in the timetable or in the shifted gradient table depending on the current focus. A new empty line is added at the end of the table, after editing the new line the table will be sorted automatically ascending by time.

OR

To delete the currently selected entry, use the delete button (**-** at the table bottom).

2D Gradient

Time [min]	% B
0.00	20.00
0.10	50.00

Shifted Gradient

Time [min]	% B
5.00	60.00

+ - ✕ ◀ ▶

The time column specifies time values relative to the runtime of the 1st dimension. In the example above, the original timetable entry 20 %B at time = 0.0 min will be changed to 42 %B at time = 4.0 min doing linear interpolation in between. Between 4.0 min and 8.0 min, the value will change from 42 %B to 25 %B.

Setup ²D Gradient graphically

The user can graphically setup the ²D gradient including the initial composition (%B) value, the ²D-stoptime and the modulation (repetition) time.




Figure 66 ²D Gradient window in edit mode

- 1 Click  to enable the graphical editing capabilities.
- 2 To add a new gradient point, move the cursor within the drawing area close to a new gradient point until the cursor changes to  and click.
- 3 To delete a gradient point, move the cursor close to the gradient point to be deleted until the cursor changes to , select the right segment and click.
- 4 To move a gradient point, move the cursor close to the gradient point to be moved until the cursor changes to , select the left segment and drag.
- 5 To change the stop time, move the cursor close to the red dotted vertical line until the cursor changes to  and drag.
- 6 To change the modulation time, move the cursor close to the green dotted vertical line until the cursor changes to  and drag.
- 7 To change the initial composition, move the cursor close to the filled circle most left near the y-axis until the cursor changes to  and drag the point.

Set ²D Time Segments

The content of the **²D Time Segments** table specifies when (within the runtime of the 1st dimension) the selected 2D-LC mode is active.

Table 8 Definitions ²D Time Segments

Column name	Description
Time	Specifies when a new segment starts (or ends)
Mode	Following options exist: - Time based The specified time defines the beginning of a time segment. - Peak based The peak detector is enabled at the specified time. - Off The time segments ends at the specified time.
Maximum peak duration (Comprehensive mode only)	Only valid in case of trigger mode = peak-based. After that time the 2D-gradient repetition ends regardless of the peak detector state.
Sampling time (Heartcutting mode only)	Set the time the loop remains in the flow path of the 1 st dimension.
Add transfer volume	- Checked: Valve is switched at the specified time plus the time to deliver the delay volume - Unchecked Valve is switched at the specified time (This check box is available only for Time based mode)


NOTE

If the **²D Time Segments** table is empty, no 2D-LC operation will be executed at all.

Set ²D Time Segments for Comprehensive mode

- To specify when a new segment starts, fill in required time in the time column. After the time defined, the 2D-gradient repetition ends regardless of the peak detector state.

Trigger table (Comprehensive)



- 2 To specify the mode and time, select **Time based**, **Peak based** or **Off** from the drop-down list in the **Mode** column fill the **Time** field.

- **Time based**

The specified time defines the beginning of a time segment where comprehensive 2D-LC is active. The ²D-gradient repetition starts immediately and ends when the ¹D-Stoptime is reached or at the time specified in the next timetable entry. The actual gradient cycle is always completed except the ¹D stoptime is reached

- **Peak based**

The peak detector is enabled at the specified time. The ²D-gradient repetition is started when a begin peak is detected and ends either with peak-end or when max. peak duration is reached. The time segment ends when the ¹D-Stoptime is reached or at the time specified in the next timetable entry. It is possible to collect multiple peaks within one time segment.

- **Off**

The time segments ends at the specified time.

- 3 In case of trigger-mode, define the the time in the **Max. peak duration** column.
- 4 To add or delete table rows, use the + and - icons below the table.

The **²D Time Segments** now are defined.

Define Peak Detector Parameter

This section allows parameterizing the peak detector to be used for peak-triggered 2D-LC operation (comprehensive or heart cutting).

Peak detector

Mode	Threshold ▼		
Upslope	5.00	mAU/s	Threshold 5.000
			mAU
Downslope	5.00	mAU/s	Upper threshold 3000.000
			mAU

Figure 67 Overview on Peak detector parameters

The *stop time* for a 2D-LC measurement must be set for the 2D pump, which can be accessed through the **advanced** settings. It must be at least the 1D run time and applies to the entire measurement, not to partial 2D-only runs/gradients for parked peaks.

Multiple Heart-Cutting *automatically extends this run time*, if required, as analyzing parked peaks takes usually longer than the 1st dimension run only.

If you define a 1D stop time, it will be applied unchanged, for example the analysis will stop after that time without processing any parked peaks. This is not recommended and will lead to a warning in the gradient preview.

⚠ Stop time mismatch

The estimated 2D runtime [12.30 min] exceeds the current instrument stop time [10.00 min]. Please set stop time in the 2D pump only and disable in the first dimension.

NOTE

If no peak detector is configured (see "[Configuration](#)" on page 148) this section is disabled. The currently configured peak detector (name & serial number of the detector) is shown in the section header.

NOTE

To facilitate the determination of parameters, it is possible to preview **Threshold** and **Slope** in the reference chromatogram.

- 1 Go to **Instrument> Setup 2DLC** and tab **Advanced**.
- 2 Select **Peak detection mode** from the drop-down list.

The following options are available:

- **Off**

The peak detector is not used.

- **Threshold only**
Detects peaks based on threshold values only.
 - **Threshold/Slope** values
Detects peaks based on both - threshold and slope.
 - **Slope only**
Detects peaks based on slope values only.
- 3 To define **Upslope** (slope of the rising peak), add the required values to the corresponding field.
 - 4 To define **Downslope** (slope of the falling peak), add the required values to the corresponding field.
 - 5 To define **Threshold** (height of the peak that triggers collection), add the required values to the corresponding field.
 - 6 To define **Upper threshold** (height of the peak that ensures that collection is not switched off even for a saturated signal that might be expected to do so), add the required values to the corresponding field.

Gradient Preview Functionality

The gradient preview provides the following functions:

- Displays the gradient (%B) of the 1st dimension pump
- Displays the gradient (%B) of the 2nd dimension during the runtime of the first dimension, depending on the selected 2D-LC mode (comprehensive OR heart-cutting)
- Allows to graphically edit the gradient shifting
- Displays a reference signal by which a user can easily setup the trigger table or optimize his peak detector settings

Reference Signal

The user can load a chromatographic signal from an LC detector - a so-called reference signal. The signal will be shown in the gradient preview. This signal is automatically shown in the gradient preview of the setup dialog as long as the reference signal is part of the method. The signal can also be removed from the method or replaced by another signal.

Loading / removing a reference signal is triggered by toolbar buttons (or the context menu) of the gradient preview.

When a reference signal is loaded, a reference signal related y-axis is shown on the right side of the gradient preview window. The grid of the graphic window is either adjusted to the absorbance axis (right axis) or the %B axis (left) and can be changed by clicking on the corresponding axis. The signal name (e.g. DAD1 A, Sig= 280, 190 Ref=550,100) is shown in the legend of the graphic window.




Figure 68 Gradient preview


The reference signal offers the following:

- Simplified set-up of (time-based) heart cuts in case the chromatogram of the sample is known in advance
- Preview of peaks, which would be analyzed in the 2nd dimension that is based on the current peak detector settings (threshold, slope)


Setup Second Dimension Gradient with the Graphical User Interface

The gradient preview allows to edit gradient shifting graphically. This replaces the editing of large timetables by a few mouse operations.


- 1 To enter the editing mode, use the context menu.




Timetable entries are marked with circles.




- 2 Drag the mouse to a new %B value at a specified runtime of the 1st dimension.




This draws a straight line. When releasing the mouse, a new timetable entry is made and the gradient rollout is automatically updated.



- 3 Repeat step 2 on page 139 with other timetable entries.



- 4 Move the mouse cursor near to a shift line to change menu context and insert or delete shift points as needed.



Method Development of Active Solvent Modulation (ASM)

ASM method development helps finding the optimal dilution of ¹D solvents in the sample loop for best ²D resolution at lowest cycle time.

After switching on the ASM functionality (see ["Method parameters"](#) on page 142), execute the steps in the following order:

- 1 ["Optimizing the dilution by using ASM capillaries"](#) on page 143
- 2 ["Optimizing the sample loop flush"](#) on page 143
- 3 ["Including the ASM phase to the 2D gradient"](#) on page 144
- 4 ["Optimizing dilution through method settings"](#) on page 145

Method parameters

The screenshot shows the 'Advanced settings' tab of a software interface. It contains two main sections: 'Peak detector' and 'Active solvent modulation (ASM)'. The 'Peak detector' section has a 'Mode' dropdown set to 'Threshold', and input fields for 'Upslope' (5.00 mAU/s), 'Downslope' (5.00 mAU/s), 'Threshold' (5.000 mAU), and 'Upper threshold' (3000.000 mAU). The 'Active solvent modulation (ASM)' section is highlighted with a red box and contains a checked 'use ASM' checkbox (with '(ASM factor 1.7)' next to it) and a 'Flush sample loop' field set to '3.0' times, with '(0.41 min)' in parentheses.

Section	Parameter	Value	Unit
Peak detector	Mode	Threshold	
	Upslope	5.00	mAU/s
	Downslope	5.00	mAU/s
	Threshold	5.000	mAU
Peak detector	Upper threshold	3000.000	mAU
	Active solvent modulation (ASM)		
Active solvent modulation (ASM)	use ASM	checked	(ASM factor 1.7)
	Flush sample loop	3.0	times (0.41 min)

Figure 69 Method parameters for the ASM Valve (example)

Advanced settings of 2D-LC method parameters allow switching on and off the use of the ASM functionality.

- If this option is off, it works as a standard 2D-LC valve without dilution.
- If this option is on, the user can set how often he wants to flush the sample loop during the ASM phase.

Optimizing the dilution by using ASM capillaries

A choice of four different ASM capillaries is available for achieving best results. Longer capillaries reduce, shorter capillaries increase the dilution of ¹D solvent in the sample loop.

Install and configure different ASM capillaries (see “Configure the ASM Valve” on page 190) for optimizing the results.

Capillary p/n	Length (mm)	Inner diameter (mm)	Volume (μl)	ASM factor	Split ratio (loop:ASM)
5500-1300	85	0.12	0.96	5	1:4
5500-1301	170	0.12	1.9	3	1:2
5500-1302	340	0.12	3.8	2	1:1
5500-1303	680	0.12	7.7	1.5	1:0.5

flow through ASM capillary
ASM factor
ASM back pressure

Optimizing the sample loop flush

Activate ASM in the software and set Flush sample loop to 3.0 times.

NOTE

Flushing the sample loop 3 times is typically enough and the recommended default. Less time may be sufficient and can be verified during optimization. The user interface displays how long this will take.

Active solvent modulation (ASM)

☒ use ASM (ASM factor 1.7)

Flush sample loop times (0.41 min)

Figure 70 Set Flush sample loop (example)

Including the ASM phase to the ^2D gradient




Figure 71 Programming the ^2D gradient table (example)

Gradients that were programmed for the second dimension originally without ASM Valve must be shifted by the delay caused by this dilution during the ASM phase such that the analytical gradient starts after the ASM phase.

If the ASM phase takes for example 0.41 min (based on selected ASM capillary, flush factor and ^2D flow rate), all times are shifted compared to a ^2D gradient without ASM.

- Gradient ends later and the gradient stop time is increased by 0.41 min
- ^2D Cycle time is increased accordingly
- One line is added to the gradient table for the ASM phase
- All times for the analytical gradient are shifted by 0.41 min.
This is true for shifted gradient steps as well (if applicable).

Optimizing dilution through method settings




Figure 72 Optimizing separation by using a lower percentage of B for the ASM and column equilibration phase (example)

For optimizing separation, you may use a lower percentage of B for the ASM phase and column equilibration phase compared to the original gradient for increasing dilution before the ²D column.

If for example the original analytical gradient started at 20 % B, you may use an ASM phase of for example 2 % B for diluting ¹D solvent more strongly during the ASM phase by changing the gradient start condition and adding a line to the ²D gradient table for the ASM phase. The starting point for the analytical gradient does not change. The solvent composition of the equilibration phase is automatically reduced to the start condition.

Apply high-resolution sampling with small cut sizes. Small cut sizes reduce the transfer of solvent volume from ¹D to ²D, which can further improve solvent compatibility and ²D resolution.

6

Run the System

Connect the capillaries to the 2D-LC valve	147
Standard Heart-Cutting 2D-LC	148
Configuration	148
Checkout/Familiarization Procedure	150
Multiple Heart-Cutting and High Resolution Sampling 2D-LC	159
Configuration	160
Checkout/Familiarization Procedure	163
Comprehensive 2D-LC	175
Configuration	176
Checkout/Familiarization Procedure	179
Active Solvent Modulation (ASM)	189
Configuration	189

This chapter describes how to run the Agilent 1290 Infinity II 2D-LC Solution in the modes standard heart-cutting, multiple heart-cutting, high resolution sampling and comprehensive 2D-LC.

Connect the capillaries to the 2D-LC valve

NOTE

Plumbing of 2D-LC valve

The correct plumbing of the 2D-LC valve differs between comprehensive versus heart-cutting mode and cocurrent versus countercurrent mode.

Use the 2D-LC software to find out the correct plumbing of the 2D-LC valve ports.

1 Install all capillaries to the 2D-LC valve.

Use the 2D-LC software to find out the correct plumbing of the 2D-LC valve ports:

- Standard heart-cutting:
 [“Configure Valve and Loop”](#) on page 149
- Multiple heart-cutting
 [“Configure Valve and Loop”](#) on page 161
- High resolution sampling
- [“Configure Valve and Loop”](#) on page 161
- Comprehensive
 [“Configure Valve and Loop”](#) on page 149

Standard Heart-Cutting 2D-LC

This section describes how to run the Agilent 1290 Infinity II 2D-LC Solution in standard heart-cutting 2D-LC.

Configuration

Check box

☒ Enable 2D-LC

Pumps

1D Pump: 2D - Binary Pump (G7120A)

2D Pump: 2D - Binary Pump (G7120A)

Detectors

Module name	Usage	Peak trigger	Transfer volume [µl]
1D-DAD (G4212B)	1D Detector	☒	10.00
2D-DAD (G7117B)	1D Detector	☐	

Columns

1D Column: Eclipse XDB-C18 (auto1D-G)

2D Column: SB-C18 (auto1D-7)

Valve topology

Select topology: 2pos/4port duo 2 loops (concurrent)

2D-LC Valve

2D-LC Valve (G1170A) 8Port2Positions42442...

Multiple Heart-Cutting Valves

Deck A (Ports 1 / 8)

Deck B (Ports 5 / 4)

Diverter Valve

Valve (G1170A) 6Port2Positions1300Bar

Waste: Port 1 -> 6 MSD; Port 1 -> 2

Capillaries

Loop: 50675926 Capillary 0.35x420 (40 µl)

Transfer: 5500-1270 Capillary ST 0.12x170 5/µl

ASM: Generic Capillary 100x0.12mm 1.1µl

Valves, loops, and capillaries

Ok Cancel

Figure 73 Overview 2D-LC configuration graphical user interface

The configuration of the 2D-LC-system is done via the configuration dialog in the software. The order of configuration is mandatory. The following configuration parameters are available:

- **Pumps**
Section to define which pump is in the first and which one in the second dimension.
- **Detectors**
Section to define which detector is in the second dimension and which detector should be used for peak detection (optional).
- **Columns**
Section to define the columns being used in 1st and 2nd dimension.
- **Valve topology, 2D-LC Valve, Multiple Heart-Cutting Valves, Diverter Valve, Capillaries**
Section to identify the modulation valve(s) used for toggling the loop(s) and section to define the volume of the sampling loop(s).

Configure Valve and Loop

To run 2D-LC, it must be defined, which valve is used for 1st and 2nd dimension.

- 1st dimension:

The **2D-LC Valve** drop-down list contains all configured valves which can be used for 2D-LC functionality.

- 2nd dimension(only relevant for multiple heart-cutting 2D-LC):

If more than one valve matches the current valve/loop configuration, the user can select from a drop-down list, which valve is used to connect 1st and 2nd dimension.

- **Diverter Valve**

If more than one valve matches the **Diverter Valve** configuration, a list-box is shown where the user can select the diverter valve.

The **Identify** button triggers the blinking of the status LED of the corresponding column compartment or valve drive. The button is only enabled in the Online version of the ChemStation.

All possible loop configurations depending on the selected valves are listed separately and illustrated on screen.




Figure 74 2D-LC valve and loop configuration (concurrent)

Software required 1290 Infinity 2D-LC Acquisition Software

Preparations Check box **Enable 2D-LC** selected

NOTE

Valves may be part of a column compartment (G7116A/B, G1316C) or a valve drive (G1170A).

- 1 Select **2D-LC Valve**.
- 2 Select **topology** under **Valve topology**.
- 3 To save settings click **OK**.

Valves and loops are configured for 2D-LC.

Checkout Familiarization Procedure

Checkout runs - heart-cutting: 1290 Infinity Binary or Quaternary LC in ¹D

The familiarization procedure illustrates the system's 2D-LC capabilities and supports the user to start the method for a specific analytical task. The familiarization procedure will guide the user through the most important setups and analysis function.

The sample provided with the familiarization procedure can be determined with a UV-detector and a mass spectrometer. The methods to analyze the starter sample are delivered together with the full package to ensure a smooth familiarization and checkout procedure. With the given method, peaks will overlap in the first dimension and will be separated in the second dimension.

The Agilent 1290 Infinity II 2D-LC Solution is delivered together with all required parts for a complete familiarization procedure for (multiple) heart-cutting and comprehensive 2D-LC.

Parts required	p/n	Description
	5190-6895	2D-LC starter sample, 1 x 2 mL Includes 2 mL
	858700-902	RRHD SB-C18, 2.1x100 mm, 1.8 µm, 1200 bar ¹ D
	857768-901	RRHD Bonus-RP, 2.1x50 mm, 1.8 µm, 1200 bar ² D, Heart-cutting
	G2453-85060	Formic Acid-Reagent Grade 5 mL (5 cc)

Hardware required

- “Single Heart-Cutting Configuration” on page 60
- Capillary ST 0.35 x 420 mm M/M 40 µl (5067-5926)

Software required CD

Preparations

Solvents needed:

- ¹D: mobile phase A = water with 0.2 % Formic Acid-Reagent Grade 5 mL (5 cc) (G2453-85060), B = methanol
- ²D: mobile phase A = water with 0.2 % Formic Acid-Reagent Grade 5 mL (5 cc) (G2453-85060), B = acetonitrile

Preparations:

- 1 Prepare dilution solvent (20 MeOH in mobile phase A): Add 300 µL MeOH to 1200 µL Mobile Phase A.
- 2 Prepare 400 µL sample: Add 40 µL 2D-LC starter sample to 360 µL dilution solvent.
- 3 Load method Checkout_MHC_Comp_1290BinX1290Bin.M from the 2D-LC Addon SW CD and modify the settings for your single heart or multiple heart cutting configuration.

1 Apply method parameters for ¹D

Table 9 Checkout method parameter settings ¹D

Module	Menu Path	Parameter	Value
Column			RRHD SB-C18, 2.1x 100 mm, 1.8 µm, 1200 bar (858700-902)
¹ D Pump	Set up Instrument Method...> Setup Method 1D Pump	Solvent A	H ₂ O + 0.2 % formic acid
		Solvent B	Methanol
		Timetable [1/100 events] (Gradient)	Time [min]:0.0 min, B[%]20
			Time [min]:50 min, B[%]100
			Stoptime: 40 min
			Posttime: 10 min
		Flow rate	0.300 mL/min
		Posttime	6 min
¹ D Column Compartment	Instrument> Setup Method 1D Column Compartment	Temperature	40 °C
¹ D Detector	Instrument> Setup Method 1D Detector	Signal A	Wavelength: 254 nm Bandwidth: 4 nm Reference Wavelength: 360 nm Reference Bandwidth: 100 nm
		Peakwidth:	5 Hz
¹ D Sampler	Instrument> Setup Method 1D Sampler	Injection volume:	2.0 µL

2 Apply method parameters for ²D.

Module	Menu Path	Parameter	Value
Column			RRHD Bonus-RP, 2.1x 50 mm, 1.8 µm, 1200 bar (857768-901)
² D Pump	Setup ² D-Pump	Solvent A	H ₂ O + 0.2 % formic acid
		Solvent B	Acetonitrile
		²D gradient	Time [min]: 0.0 min, B[%] 10
			Time [min]: 1.25 min, B[%] 60
		²D Gradient stoptime:	1.25 min
		²D cycletime	1.75 min
		²D Flow	1.0 mL/min
	Setup ² D-Pump Pre-viewEdit mode	Gradient shift:	0->20 min from 10->30 % B (only downslope)
	Setup ² D-Pump Advanced ²D pump settings...	Stoptime	40 min (will be automatically prolonged, if peaks in ² D are not worked off)
¹ D Column Compartment		Temperature	40 °C
² D Detector	Instrument> Setup Method 2D Detector	Signal A	Wavelength: 254 nm Bandwidth: 4 nm Reference Wavelength: 360 nm Reference Bandwidth: 100 nm
		Peakwidth:	≥40 Hz

3 Program and/or find the following cuts in the predefined method:

Cut-#	Cut-Time [min] 1290 Binary LC	Cut-Time [min] 1290 Quaternary LC ¹
1	4.25	4.35
2	6.58	6.86
3	10.05	10.4
4	13.3	13.7
5	16.8	17.15
6	23.9	24.6

¹ The Cut-Time can vary slightly depending on the configuration

4 Run the method with the 2D-LC starter sample, 1 x 2 mL (5190-6895), 1:10 diluted with Methanol/Water (20/80; v/v) with 0.1 % formic acid.

Checkout runs - heart-cutting: 1260 Infinity Binary in ¹D

The familiarization procedure illustrates the system's 2D-LC capabilities and supports the user to start the method for a specific analytical task. The familiarization procedure will guide the user through the most important setups and analysis function.

The sample provided with the familiarization procedure can be determined with a UV-detector and a mass spectrometer. The methods to analyze the starter sample are delivered together with the full package to ensure a smooth familiarization and checkout procedure. With the given method, peaks will overlap in the first dimension and will be separated in the second dimension.

The Agilent 1290 Infinity II 2D-LC Solution is delivered together with all required parts for a complete familiarization procedure for (multiple) heart-cutting and comprehensive 2D-LC.

Parts required	p/n	Description
	5190-6895	2D-LC starter sample, 1 x 2 mL Includes 2 mL
	858700-902	RRHD SB-C18, 2.1x100 mm, 1.8 µm, 1200 bar ¹ D
	857768-901	RRHD Bonus-RP, 2.1x50 mm, 1.8 µm, 1200 bar ² D, Heart-cutting
	G2453-85060	Formic Acid-Reagent Grade 5 mL (5 cc)

Hardware required

- "Single Heart-Cutting Configuration" on page 60
- Capillary ST 0.35 x 420 mm M/M 40 µl (5067-5926)

Software required CD

Preparations

Solvents needed:

- ¹D: mobile phase A = water with 0.2 % Formic Acid-Reagent Grade 5 mL (5 cc) (G2453-85060), B = methanol
- ²D: mobile phase A = water with 0.2 % Formic Acid-Reagent Grade 5 mL (5 cc) (G2453-85060), B = acetonitrile

Preparations:

- 1 Prepare dilution solvent (20 MeOH in mobile phase A): Add 300 µL MeOH to 1200 µL Mobile Phase A.
- 2 Prepare 400 µL sample: Add 40 µL 2D-LC starter sample to 360 µL dilution solvent.
- 3 Load method Checkout_MHC_Comp_1290BinX1290Bin.M from the 2D-LC Addon SW CD and modify the settings for your single heart or multiple heart cutting configuration.

1 Apply method parameters for ¹D

Table 10 Checkout method parameter settings ¹D

Module	Menu Path	Parameter	Value
Column			RRHD SB-C18, 2.1x 100 mm, 1.8 µm, 1200 bar (858700-902)
¹ D Pump	Set up Instrument Method...> Setup Method 1D Pump	Solvent A	H ₂ O + 0.2 % formic acid
		Solvent B	Methanol
		Timetable [1/100 events] (Gradient)	Time [min]:0.0 min, B[%]20
			Time [min]:50 min, B[%]100
			Stoptime: 40 min
			Posttime: 10 min
		Flow rate	0.300 mL/min
		Posttime	6 min
¹ D Column Compartment	Instrument> Setup Method 1D Column Compartment	Temperature	40 °C
¹ D Detector	Instrument> Setup Method 1D Detector	Signal A	Wavelength: 254 nm Bandwidth: 4 nm Reference Wavelength: 360 nm Reference Bandwidth: 100 nm
		Peakwidth:	5 Hz
¹ D Sampler	Instrument> Setup Method 1D Sampler	Injection volume:	2.0 µL

2 Apply method parameters for ²D.

Module	Menu Path	Parameter	Value
Column			RRHD Bonus-RP, 2.1x 50 mm, 1.8 µm, 1200 bar (857768-901)
² D Pump	Setup ² D-Pump	Solvent A	H ₂ O + 0.2 % formic acid
		Solvent B	Acetonitrile
		²D gradient	Time [min]: 0.0 min, B[%] 10
			Time [min]: 1.25 min, B[%] 60
		²D Gradient stoptime:	1.25 min
		²D cycletime	1.75 min
		²D Flow	1.0 mL/min
	Setup ² D-Pump Pre-viewEdit mode	Gradient shift:	0->20 min from 10->30 % B (only downslope)
	Setup ² D-Pump Advanced ²D pump settings...	Stoptime	40 min (will be automatically prolonged, if peaks in ² D are not worked off)
¹ D Column Compartment		Temperature	40 °C
² D Detector	Instrument> Setup Method 2D Detector	Signal A	Wavelength: 254 nm Bandwidth: 4 nm Reference Wavelength: 360 nm Reference Bandwidth: 100 nm
		Peakwidth:	≥40 Hz

3 Program and/or find the following cuts in the predefined method:

Cut-#	Cut-Time [min] 1260 Binary LC ¹
1	9.5
2	13.13
3	17.6
4	21.2
5	25.25
6	31.55

¹ The Cut-Time can vary slightly depending on the configuration

4 Run the method with the 2D-LC starter sample, 1 x 2 mL (5190-6895), 1:10 diluted with Methanol/Water (20/80; v/v) with 0.1 % formic acid.

Multiple Heart-Cutting and High Resolution Sampling 2D-LC

This section describes how to run the Agilent 1290 Infinity II 2D-LC Solution in multiple heart-cutting and high resolution sampling 2D-LC.

Configuration

Overview Configuration Dialog

Module name	Usage	Peak trigger	Transfer volume [µl]
DAD 1D (G7117B)	1D Detec...	☒	0.00
DAD 2D (G1315A)	2D Detector		
G6110A MSD (G6110A)	2D Detector		

Figure 75 Overview 2D-LC configuration graphical user interface

The configuration of the 2D-LC-system is done via the configuration dialog in the software. The order of configuration is mandatory. The following configuration parameters are available:

- **Pumps**

Section to define which pump is in the first and which one in the second dimension.

- **Detectors**

Section to define which detector is in the second dimension and which detector should be used for peak detection (optional).

- **Columns**

Section to define the columns being used in 1st and 2nd dimension.

- **Valve topology, 2D-LC Valve, Multiple Heart-Cutting Valves, Diverter Valve, Capillaries**

Section to identify the modulation valve(s) used for toggling the loop(s) and section to define the volume of the sampling loop(s).

Configure Valve and Loop

To run Multiple Heart-cutting/High-resolution sampling 2D-LC, the following parameters must be defined:

- **Select topology:**
The drop-down list contains all possible valve/loop topologies.
- **2D-LC valve(s):**
 - **#1:**
If more than one valve matches the current valve / loop configuration, a list-box is shown where the user can select the 2D-LC valve.
 - **#2:**
If more than one valve matches the current valve / loop configuration, a list-box is shown where the user can select the 2D-LC valve.

Shows the 2nd valve to be used for the injection on the 2nd dimension. This field is only visible if a **Dual 2pos/6port** configuration is selected.
- **Multiple Heart-Cutting Valve(s)**
 - **Deck A (Ports 1 / 8):**
Shows the 6port/14port valve used for the parking deck 1. The field is only enabled in case of a valve / loop configuration with one or two Multiple Heart-Cutting Valves (6 or 12 loops). If more than one 6/14 valve is available, a list-box is shown where the user can select the appropriate valve.
 - **Deck B (Ports 5 / 4):**
Select the 6port/14port valve used for the parking deck 2. The list-box is only enabled in case of a valve / loop configuration with two Multiple Heart-Cutting Valves (6 or 12 loops).
- **Diverter Valve**
If more than one valve matches the **Diverter Valve** configuration, a list-box is shown where the user can select the diverter valve.
- **Capillaries**
A list-box is shown where the user can select the capillaries.
- The **Identify** button triggers the blinking of the status LED of the corresponding valve or column compartment. The button is only enabled in the Online version of the ChemStation.

All possible loop configurations depending on the selected valves are listed separately and illustrated on screen.




Figure 76 2D-LC valve and loop configuration

Preparations

- OpenLAB ChemStation Edition C.01.10 (or higher)
- 1290 Infinity 2D-LC Acquisition Software
- Check box **Enable 2D-LC** selected

NOTE

Valves must be part of the 1290 Infinity Valve Drive (G1170A).

- 1 Select the valve / loop combination for the injection on the 1st dimension (**Select topology**).
- 2 Specify the volume of the loop (**Loop size**).
- 3 Select the valve for the injection on the 2nd dimension (**2D-LC valve(s)**).
- 4 Select the valve for parking peaks (**Multiple Heart-Cutting Valve(s)**).
- 5 Select the capillaries used in your system.
- 6 To save settings click **OK**.

Checkout/Familiarization Procedure

Checkout runs - MHC: 1290 Infinity Binary or Quaternary LC in ¹D

Parts required	p/n	Description
	5190-6895	2D-LC starter sample, 1 x 2 mL Includes 2 mL
	858700-902	RRHD SB-C18, 2.1x100 mm, 1.8 µm, 1200 bar ¹ D
	857768-901	RRHD Bonus-RP, 2.1x50 mm, 1.8 µm, 1200 bar ² D, Heart-cutting
	G2453-85060	Formic Acid-Reagent Grade 5 mL (5 cc)

Hardware required ["Multiple Heart-Cutting Configuration"](#) on page 61

Software required CD

- Preparations**
- Solvents needed:
- ¹D: mobile phase A = water with 0.2 % Formic Acid-Reagent Grade 5 mL (5 cc) (G2453-85060), B = methanol
 - ²D: mobile phase A = water with 0.2 % Formic Acid-Reagent Grade 5 mL (5 cc) (G2453-85060), B = acetonitrile
- Preparations:
- 1 Prepare dilution solvent (20 MeOH in mobile phase A): Add 300 µL MeOH to 1200 µL Mobile Phase A.
 - 2 Prepare 400 µL sample: Add 40 µL 2D-LC starter sample to 360 µL dilution solvent.
 - 3 Load method Checkout_MHC_Comp_1290BinX1290Bin.M from the 2D-LC Addon SW CD and modify the settings for your single heart or multiple heart cutting configuration.

1 Apply method parameters for ¹DTable 11 Checkout method parameter settings ¹D

Module	Menu Path	Parameter	Value
Column			RRHD SB-C18, 2.1x 100 mm, 1.8 µm, 1200 bar (858700-902)
¹ D Pump	Set up Instrument Method...> Setup Method 1D Pump	Solvent A	H ₂ O + 0.2 % formic acid
		Solvent B	Methanol
		Timetable [1/100 events] (Gradient)	Time [min]:0.0 min, B[%]20
			Time [min]:50 min, B[%]100
			Stoptime: 40 min
			Posttime: 10 min
		Flow rate	0.300 mL/min
		Posttime	6 min
¹ D Column Compartment	Instrument> Setup Method 1D Column Compartment	Temperature	40 °C
¹ D Detector	Instrument> Setup Method 1D Detector	Signal A	Wavelength: 254 nm Bandwidth: 4 nm Reference Wavelength: 360 nm Reference Bandwidth: 100 nm
		Peakwidth:	5 Hz
¹ D Sampler	Instrument> Setup Method 1D Sampler	Injection volume:	2.0 µL

2 Apply method parameters for ²D.

Module	Menu Path	Parameter	Value
Instrument	Instrument> Setup 2D-LC>>		Verify that the 2D-LC mode is set to Multiple Heart-Cutting
Column			RRHD Bonus-RP, 2.1x 50 mm, 1.8 µm, 1200 bar (857768-901)
² D Pump	Setup ² D-Pump	Solvent A	H ₂ O + 0.2 % formic acid
		Solvent B	Acetonitrile
		²D gradient	Time [min]: 0.0 min, B[%] 10
			Time [min]: 1.25 min, B[%] 60
		²D Gradient stoptime:	1.25 min
		²D cycletime	1.75 min
		²D Flow	1.0 mL/min
	Setup ² D-Pump Pre-viewEdit mode	Gradient shift:	0->20 min from 10->30 % B (only downslope)
	Setup ² D-Pump Advanced ²D pump settings...	Stoptime	40 min (will be automatically prolonged, if peaks in ² D are not worked off)
¹ D Column Compartment		Temperature	40 °C
² D Detector	Instrument> Setup Method 2D Detector	Signal A	Wavelength: 254 nm Bandwidth: 4 nm Reference Wavelength: 360 nm Reference Bandwidth: 100 nm
		Peakwidth:	≥40 Hz

3 Program and/or find the following cuts in the predefined method:

Cut-#	Cut-Time [min] 1290 Binary LC	Cut-Time [min] 1290 Quaternary LC ¹
1	4.25	4.35
2	6.58	6.86
3	10.05	10.4
4	12.19	12.58
5	13.3	13.7
6	13.44	13.85
7	13.58	14
8	13.72	14.15
9	13.86	14.3
10	16.8	17.15
11	16.94	17.3
12	17.08	17.45
13	17.22	17.6
14	17.36	17.75
15	23.9	24.6

¹ The Cut-Time can vary slightly depending on the configuration

4 Run the method with the 2D-LC starter sample, 1 x 2 mL (5190-6895), 1:10 diluted with Methanol/Water (20/80; v/v) with 0.1 % formic acid.

Checkout runs - MHC: 1260 Infinity Binary in ¹D

Parts required	p/n	Description
	5190-6895	2D-LC starter sample, 1 x 2 mL Includes 2 mL
	858700-902	RRHD SB-C18, 2.1x100 mm, 1.8 µm, 1200 bar ¹ D
	857768-901	RRHD Bonus-RP, 2.1x50 mm, 1.8 µm, 1200 bar ² D, Heart-cutting
	G2453-85060	Formic Acid-Reagent Grade 5 mL (5 cc)

Hardware required "Multiple Heart-Cutting Configuration" on page 61

Software required CD

Preparations

Solvents needed:

- ¹D: mobile phase A = water with 0.2 % Formic Acid-Reagent Grade 5 mL (5 cc) (G2453-85060), B = methanol
- ²D: mobile phase A = water with 0.2 % Formic Acid-Reagent Grade 5 mL (5 cc) (G2453-85060), B = acetonitrile

Preparations:

- 1 Prepare dilution solvent (20 MeOH in mobile phase A): Add 300 µL MeOH to 1200 µL Mobile Phase A.
- 2 Prepare 400 µL sample: Add 40 µL 2D-LC starter sample to 360 µL dilution solvent.
- 3 Load method Checkout_MHC_Comp_1290BinX1290Bin.M from the 2D-LC Addon SW CD and modify the settings for your single heart or multiple heart cutting configuration.

1 Apply method parameters for ¹DTable 12 Checkout method parameter settings ¹D

Module	Menu Path	Parameter	Value
Column			RRHD SB-C18, 2.1x 100 mm, 1.8 µm, 1200 bar (858700-902)
¹ D Pump	Set up Instrument Method...> Setup Method 1D Pump	Solvent A	H ₂ O + 0.2 % formic acid
		Solvent B	Methanol
		Timetable [1/100 events] (Gradient)	Time [min]:0.0 min, B[%]20
			Time [min]:50 min, B[%]100
			Stoptime: 40 min
			Posttime: 10 min
		Flow rate	0.300 mL/min
		Posttime	6 min
¹ D Column Compartment	Instrument> Setup Method 1D Column Compartment	Temperature	40 °C
¹ D Detector	Instrument> Setup Method 1D Detector	Signal A	Wavelength: 254 nm Bandwidth: 4 nm Reference Wavelength: 360 nm Reference Bandwidth: 100 nm
		Peakwidth:	5 Hz
¹ D Sampler	Instrument> Setup Method 1D Sampler	Injection volume:	2.0 µL

2 Apply method parameters for ²D.

Module	Menu Path	Parameter	Value
Instrument	Instrument> Setup 2D-LC>>		Verify that the 2D-LC mode is set to Multiple Heart-Cutting
Column			RRHD Bonus-RP, 2.1x 50 mm, 1.8 µm, 1200 bar (857768-901)
² D Pump	Setup ² D-Pump	Solvent A	H ₂ O + 0.2 % formic acid
		Solvent B	Acetonitrile
		²D gradient	Time [min]: 0.0 min, B[%] 10
	Setup ² D-Pump Pre-viewEdit mode		Time [min]: 1.25 min, B[%] 60
		²D Gradient stoptime:	1.25 min
		²D cycletime	1.75 min
		²D Flow	1.0 mL/min
		Gradient shift:	0->20 min from 10->30 % B (only downslope)
		Stoptime	40 min (will be automatically prolonged, if peaks in ² D are not worked off)
	Setup ² D-Pump Advanced ²D pump settings...		
¹ D Column Compartment		Temperature	40 °C
² D Detector	Instrument> Setup Method 2D Detector	Signal A	Wavelength: 254 nm Bandwidth: 4 nm Reference Wavelength: 360 nm Reference Bandwidth: 100 nm
		Peakwidth:	≥40 Hz

3 Program and/or find the following cuts in the predefined method:

Cut-#	Cut-Time [min] 1260 Binary LC ¹
1	9.5
2	13.13
3	17.6
4	20.2
5	20.45
6	20.7
7	21.2
8	21.45
9	21.7
10	22.2
11	24.6
12	25.25
13	25.5
14	25.75
15	27

¹ The Cut-Time can vary slightly depending on the configuration

4 Run the method with the 2D-LC starter sample, 1 x 2 mL (5190-6895), 1:10 diluted with Methanol/Water (20/80; v/v) with 0.1 % formic acid.

Checkout run - high-resolution sampling

Parts required	p/n	Description
	5190-6895	2D-LC starter sample, 1 x 2 mL Includes 2 mL
	858700-902	RRHD SB-C18, 2.1x100 mm, 1.8 µm, 1200 bar ¹ D
	857768-901	RRHD Bonus-RP, 2.1x50 mm, 1.8 µm, 1200 bar ² D, Heart-cutting
	G2453-85060	Formic Acid-Reagent Grade 5 mL (5 cc)

Hardware required "Multiple Heart-Cutting Configuration" on page 61

Software required CD

Preparations Solvents needed:

- ¹D: mobile phase A = water with 0.2 % Formic Acid-Reagent Grade 5 mL (5 cc) (G2453-85060), B = methanol
- ²D: mobile phase A = water with 0.2 % Formic Acid-Reagent Grade 5 mL (5 cc) (G2453-85060), B = acetonitrile

Preparations:

- 1 Prepare dilution solvent (20 MeOH in mobile phase A): Add 300 µL MeOH to 1200 µL Mobile Phase A.
- 2 Prepare 400 µL sample: Add 40 µL 2D-LC starter sample to 360 µL dilution solvent.
- 3 Load method Checkout_MHC_Comp_1290BinX1290Bin.M from the 2D-LC Addon SW CD and modify the settings for your single heart or multiple heart cutting configuration.

1 Apply method parameters for ¹D.Table 13 Checkout method parameter settings ¹D

Module	Tool	Parameter	Value
Column			RRHD SB-C18, 2.1x 100 mm, 1.8 µm, 1200 bar (858700-902)
¹ D Pump	Instrument> Setup Method 1D Pump	Solvent A	H ₂ O + 0.2 % formic acid
		Solvent B	Methanol
		Timetable [1/100 events] (Gradient)	Time [min]: 0.0 min, B[%] 20
			Time [min]: 20 min, B[%] 55
			Time [min]: 21 min, B[%] 95
			Stoptime: select As Injector/No Limit
			Posttime: select Off
		Flow rate	0.600 mL/min
¹ D Column Compartment	Instrument> Setup Method 1D Column Compartment	Posttime	Off
		Temperature	40 °C
¹ D Detector	Instrument> Setup Method 1D Detector	Signal A	Wavelength: 254 nm Bandwidth: 4 nm Reference Wavelength: 360 nm Reference Bandwidth: 100 nm
		Signal B	Wavelength: 260 nm Bandwidth: 4 nm Reference Wavelength: 360 nm Reference Bandwidth: 100 nm
		Peakwidth:	5 Hz
¹ D Sampler	Instrument> Setup Method 1D Sampler	Injection volume:	1.0 µL , 2.0 µL , 5.0 µL , and 3.0 µL

2 Apply method parameters for ²D.

Module	Tool	Parameter	Value
Column			RRHD Bonus-RP, 2.1x 50 mm, 1.8 µm, 1200 bar (857768-901)
² D Pump	Setup ² D-Pump	Solvent A	H ₂ O + 0.2 % formic acid
		Solvent B	Acetonitrile
		²D gradient	Time [min]: 0.0 min, B[%] 10
			Time [min]: 1.25 min, B[%] 60
		²D Gradient stoptime:	1.25 min
		²D cycletime	1.75 min
		²D Flow	1.0 mL/min
	Setup ² D-Pump Pre-viewEdit mode	Gradient shift:	0->20 min from 10->30 % B (only downslope)
	Setup ² D-Pump Advanced ²D pump settings...	Stoptime	23 min
		Posttime	3 min
¹ D Column Compartment		Temperature	40 °C
² D Detector	Instrument> Setup Method 2D Detector	Signal A	Wavelength: 254 nm Bandwidth: 4 nm Reference Wavelength: 360 nm Reference Bandwidth: 100 nm
		Signal B	Wavelength: 260 nm Bandwidth: 4 nm Reference Wavelength: 360 nm Reference Bandwidth: 100 nm
		Peakwidth:	≥40 Hz

- 3 In **2D-LC Mode** select **HiRes sampling**.
- 4 Run a survey run with the 2D-LC starter sample, 1 x 2 mL (5190-6895), 1:10 diluted with Methanol/Water (20/80; v/v) with 0.1 % formic acid of the first dimension to find out the exact retention time of the peak that elutes at around 11.7 min.
- 5 Run HiRes sampling with the sampling time adjusted to the ¹D retention time of the peak of interest in your system.
Settings for High-Resolution sampling:
 - Sampling time: 3 s
 - Number of cuts: 10
 - Injection volume: 1, 2, 5, and 3 μ L,

Comprehensive 2D-LC

This section describes in detail the installation, configuration, method parameters, data analysis and checkout/familiarization of comprehensive two dimensional liquid chromatography with the Agilent 1290 Infinity II 2D-LC Solution.

Configuration

Overview Configuration Dialog

Configure 2D-LC: 2D-LC Testsystem

Check box ☒ Enable 2D-LC

Pumps

1D Pump: 1D - Binary Pump (G4220B) ?

2D Pump: 2D - Binary Pump (G7120A) ?

Detectors

Module name	Usage	Peak trigger	Transfer volume (µl)
1D-DAD (G4212B)	1D Detector	☒	10.00
2D-DAD (G7117B)	2D Detector	☐	


Columns

1D Column: Eclipse XDB-C18 (autoID-G) ?

2D Column: SB-C18 (autoID-7) ?

Valve topology

Select topology: 2pos/4port diva 2 loops (concurrent) ?



2D-LC Valve

2D-LC Valve (G1170A) 8Port2Positions42442... ?

Multiple Heart-Cutting Valves

Deck A (Ports 1 / 8) ?

Deck B (Ports 5 / 4) ?

Diverter Valve

Valve (G1170A) 6Port2Positions13008ar ?

Waste: Port 1 -> 6 MSD; Port 1 -> 2

Capillaries ...

Loop: 5067-5926 Capillary 0.35x420 (40 µl)

Transfer: 5500-1270 Capillary ST 0.12x170 5/µl

ASM: Generic Capillary 100x0.12mm 1.1µl

Ok Cancel

Valves, loops, and capillaries

Figure 77 Overview 2D-LC configuration graphical user interface

The configuration of the 2D-LC-system is done via the configuration dialog in the software. The order of configuration is mandatory. The following configuration parameters are available:

- **Pumps**
Section to define which pump is in the first and which one in the second dimension.
- **Detectors**
Section to define which detector is in the second dimension and which detector should be used for peak detection (optional).
- **Columns**
Section to define the columns being used in 1st and 2nd dimension.
- **Valve topology, 2D-LC Valve, Multiple Heart-Cutting Valves, Diverter Valve, Capillaries**
Section to identify the modulation valve(s) used for toggling the loop(s) and section to define the volume of the sampling loop(s).

Configure Valve and Loop

To run 2D-LC, it must be defined, which valve is used for 1st and 2nd dimension.

- 1st dimension:

The **2D-LC Valve** drop-down list contains all configured valves which can be used for 2D-LC functionality.

- 2nd dimension(only relevant for multiple heart-cutting 2D-LC):

If more than one valve matches the current valve/loop configuration, the user can select from a drop-down list, which valve is used to connect 1st and 2nd dimension.

- **Diverter Valve**

If more than one valve matches the **Diverter Valve** configuration, a list-box is shown where the user can select the diverter valve.

The **Identify** button triggers the blinking of the status LED of the corresponding column compartment or valve drive. The button is only enabled in the Online version of the ChemStation.

All possible loop configurations depending on the selected valves are listed separately and illustrated on screen.




Figure 78 2D-LC valve and loop configuration (concurrent)

Software required 1290 Infinity 2D-LC Acquisition Software

Preparations Check box **Enable 2D-LC** selected

NOTE

Valves may be part of a column compartment (G7116A/B, G1316C) or a valve drive (G1170A).

- 1 Select **2D-LC Valve**.
- 2 Select **topology** under **Valve topology**.
- 3 To save settings click **OK**.

Valves and loops are configured for 2D-LC.

Checkout/Familiarization Procedure

Checkout run - comprehensive (standard setup)

The familiarization procedure illustrates the system's 2D-LC capabilities and supports the user to start the method for a specific analytical task. The familiarization procedure will guide the user through the most important setups and analysis function, described in the chapters before.

The sample provided with the familiarization procedure can be determined with a UV-detector and a mass spectrometer. The methods to analyze the starter sample are delivered together with the full package to ensure a smooth familiarization and checkout procedure. With the given method, peaks will overlap in the first dimension and will be separated in the second dimension.

The Agilent 1290 Infinity II 2D-LC Solution is delivered together with all required parts for a complete familiarization procedure for (multiple) heart-cutting and comprehensive 2D-LC.

Parts required	p/n	Description
	5190-6895	2D-LC starter sample, 1 x 2 mL
	858700-902	RRHD SB-C18, 2.1x100 mm, 1.8 µm, 1200 bar 1D
	959757-302	RRHD Eclipse Plus C18, 3.0x50 mm, 1.8 µm
	G2453-85060	Formic Acid-Reagent Grade 5 mL (5 cc)

Software required CD

- Preparations**
- Solvents needed:
- 1D: mobile phase A = water with 0.2 % Formic Acid-Reagent Grade 5 mL (5 cc) (G2453-85060), B = methanol
 - 2D: mobile phase A = water with 0.2 % Formic Acid-Reagent Grade 5 mL (5 cc) (G2453-85060), B = acetonitrile
- Preparations:
- 1 Prepare dilution solvent (20 % MeOH in mobile phase A): Add 500 µL MeOH to 2000 µL Mobile Phase A.
 - 2 Prepare 400 µL sample: Add 40 µL 2D LC starter sample to 360 µL dilution solvent.
 - 3 Load method Checkout_MHC_FullComp_1290BinX1290Bin.MCheckout_MHC_FullComp_1290QuatX1290Bin.M, or Checkout_MHC_FullComp_1260BinX1290Bin.M(depends on your pump setup) from the CD.


- Run the full comprehensive checkout run using single storage loops. The method for the checkout run is available on the 2D-LC Addon SW DVD under checkout methods.
 - Detection:

UV Detection at 254 nm, BW 4 nm; reference at 360 nm, BW 100 nm
 - Acquisition rate:

20 Hz
 - Sample:

2D-LC starter sample, 1 x 2 mL (5190-6895), 1:10 diluted with Methanol/Water (20/80; v/v) with 0.1 % formic acid.
 - Injection volume:

2 μ L
- Verify that the Valve & Loop configuration is set up as shown below:



- 3 Verify the ¹D pump method is set up as shown below:

Method of G4220A (DE92900288)

1st Dim Binary Pump (G4220A)

Flow
0.100 mL/min

Solvents

A: 60.00 %
1: 100.0 % Water V.03
2: Organic
0.2 % FA

B: 40.00 %
1: 100.0 % Methanol V.03
2: Aqueous

Pressure Limits
Min: 0.00 bar Max: 1,200.00 bar

Stop time
☐ As Injector/No Limit
☒ 40.00 min


Post time
☐ Off
☒ 10.00 min

Advanced
☒ Timetable (2/100 events)


☐ function centric view

Time [min]	A [%]	B [%]	Flow [mL/min]	Max. Pressure Limit [bar]
0.00	60.00	40.00	0.100	1200.00
34.00	40.00	60.00	---	---
34.50	10.00	90.00	---	---

- 4 In **Instrument> Setup 2D-LC**, verify that the 2D-LC mode is set to **Comprehensive** and that the 2D pump and modulation method are set up as shown below:



- 5 Run the full comprehensive checkout run using single storage loops and review the obtained data with the GC Image software. The resulting separation should look similar to the one shown below:



Checkout run - comprehensive (multiple heart-cutting setup)


The familiarization procedure illustrates the system's 2D-LC capabilities and supports the user to start the method for a specific analytical task. The familiarization procedure will guide the user through the most important setups and analysis function, described in the chapters before.

The sample provided with the familiarization procedure can be determined with a UV-detector and a mass spectrometer. The methods to analyze the starter sample are delivered together with the full package to ensure a smooth familiarization and checkout procedure. With the given method, peaks will overlap in the first dimension and will be separated in the second dimension.

The Agilent 1290 Infinity II 2D-LC Solution is delivered together with all required parts for a complete familiarization procedure for (multiple) heart-cutting and comprehensive 2D-LC.

Parts required	p/n	Description
	5190-6895	2D-LC starter sample, 1 x 2 mL
	858700-902	RRHD SB-C18, 2.1x100 mm, 1.8 µm, 1200 bar 1D
	959757-302	RRHD Eclipse Plus C18, 3.0x50 mm, 1.8 µm
	G2453-85060	Formic Acid-Reagent Grade 5 mL (5 cc)
Software required	CD	
Preparations	Solvents needed: 1D: mobile phase A = water with 0.2 % Formic Acid-Reagent Grade 5 mL (5 cc) (G2453-85060), B = methanol 2D: mobile phase A = water with 0.2 % Formic Acid-Reagent Grade 5 mL (5 cc) (G2453-85060), B = acetonitrile Preparations: 1 Prepare dilution solvent (20 % MeOH in mobile phase A): Add 500 µL MeOH to 2000 µL Mobile Phase A. 2 Prepare 400 µL sample: Add 40 µL 2D LC starter sample to 360 µL dilution solvent. 3 Load method Checkout_MHC_FullComp_1290BinX1290Bin.MCheckout_MHC_FullComp_1290QuatX1290Bin.M, or Checkout_MHC_FullComp_1260BinX1290Bin.M(depends on your pump setup) from the CD.	

- 1 Repeat the full comprehensive checkout run using MHC valves instead of single storage loops. For this purpose, disconnect the transfer capillaries from the 2D-LC valve to the storage loops and install MHC valves between ports 4 and 5, respectively ports 1 and 8, of the 2D-LC valve (cocurrent configuration).
- 2 In **Instrument> Configure 2D-LC**, change the Valve & Loop configuration to 40 μ L storage loops as shown below:



- Detection:
UV Detection at 254 nm, BW 4 nm; reference at 360 nm, BW 100 nm
- Acquisition rate:
20 Hz
- Sample:
2D-LC starter sample, 1 x 2 mL (5190-6895), 1:10 diluted with Methanol/Water (20/80; v/v) with 0.1 % formic acid.
- Injection volume:
2 μ L

- 3 Verify the ¹D pump method is set up as shown below:

Method of G4220A (DE92900288)

1st Dim Binary Pump (G4220A)

Flow: 0.100 mL/min

Solvents

A: 60.00 % 1: 100.0 % Water V.03 0.2 % FA
2: Organic

B: 40.00 % 1: 100.0 % Methanol V.03
2: Aqueous

Pressure Limits
Min: 0.00 bar Max: 1,200.00 bar

Stop time
☐ As Injector/No Limit
☒ 40.00 min

Post time
☐ Off
☒ 10.00 min

Advanced
☒ Timetable (2/100 events)


☐ function centric view

Time [min]	A [%]	B [%]	Flow [mL/min]	Max. Pressure Limit [bar]
0.00	60.00	40.00	0.100	1200.00
34.00	40.00	60.00	---	---
34.50	10.00	90.00	---	---


Add Remove Clear All Clear Empty
Cut Copy Paste Shift Times 0.00 min

Ok Apply Cancel

- 4 In **Instrument> Setup 2D-LC**, verify that the 2D-LC mode is set to **Comprehensive** and that the 2D pump and modulation method are set up as shown below:



- 5 Run the full comprehensive checkout run using single storage loops. The resulting separation should look similar to the one shown below:



- 6 Compare the result to that obtained using the MHC valves.

Active Solvent Modulation (ASM)

Configuration

Adjusting the split ratio

Different ASM capillaries are available for adjusting the split ratio and therefore the dilution.

The method can therefore be optimized either for optimum resolution (strong dilution) or lowest cycle time (weak dilution).

Configure the ASM Valve

Configure 2D-LC: 2DLC-MSD

☒ Enable 2D-LC

Pumps

¹D Pump: 1D Quat. Pump (G7104A)

²D Pump: 2D Binary Pump (G7120A)

Detectors

Module name	Usage	Peak trigger	Transfer volume [μl]
1D DAD (G1315C)	¹ D Detector	☒	10.00
2D DAD (G7117B)	² D Detector		
G6110A MSD (G6110A)	¹ D Detector		

Columns

¹D Column: SB-C18 (autoID-7)

²D Column: Eclipse Plus C18 (autoID-10)

Valve topology

1 Select topology: 2D-LC Valve ASM 2x6 loops (countercurrent)

2D-LC Valve

2 2D-LC Valve ASM (G1170A:) Generic

Multiple Heart-Cutting Valves

Deck A (Ports 1 / 8)
Deck A (G1170A:) 14Port6Positions1300BarNpl

Deck B (Ports 5 / 4)
Deck B (G1170A:) 14Port6Positions1300Bar

Diverter Valve
Not configured

Capillaries

Loop: 5067-5926 Capillary 0.35x420 (40 μl)

Transfer: 5500-1270 Capillary 0.12x170 (1.9 μl)

3 ASM: Generic Capillary 500x0.12mm 5.7ul

Ok Cancel

Figure 79 ASM Valve configuration (overview)

- 1 Select a topology for using the ASM Valve.
- 2 Choose the ASM Valve as 2D-LC Valve. This is usually done automatically based on installed valves.
- 3 Choose an ASM capillary. This defines the split ratio.

Preparations

All modules including the ASM Valve are configured in OpenLAB CDS ChemStation Edition.

- 1 Select a valve topology using an ASM Valve.

HINT

For minimum carry-over, please use counter-current installation for the ASM Valve.

- 2 Select the ASM Valve as 2D-LC Valve (which is usually pre-selected).

- 3 Define the ASM capillary.
 - a To configure capillaries, click on **Capillaries...** (see [Figure 79](#) on page 190).
 - b Select any of the pre-defined ASM capillaries.

The 'Setup Capillaries' dialog box displays a table of predefined capillaries. The 'ASM capillary' row is highlighted with a blue selection bar.

	Capillary Name (P/N)	Length [mm]	Diameter [mm]	Volume [μ l]
Sample loop capillary	5067-5926 Capillary 0.35x420 (40 μ l)	420	0.35	40.4
Transfer capillary between 2D-LC valve and MHC-valve	5500-1270 Capillary 0.12x170 (1.9 μ l)	170	0.12	1.9
ASM capillary	5500-1300 Capillary 0.12x85 (1.0 μ l)	85	0.12	1.0

Below the table, the 'ASM factor' is set to 5.1. 'Ok' and 'Cancel' buttons are at the bottom right.

Figure 80 Configuration of the ASM valve with predefined capillaries

OR

If you are using a different capillary, you can choose **Generic Capillary**

The 'Setup Capillaries' dialog box shows the 'Generic Capillary' selected in the 'ASM capillary' dropdown. The length is set to 500 mm and the volume is 5.7 μ l.

	Capillary Name (P/N)	Length [mm]	Diameter [mm]	Volume [μ l]
Sample loop capillary	5067-5926 Capillary 0.35x420 (40 μ l)	420	0.35	40.4
Transfer capillary between 2D-LC valve and MHC-valve	5500-1270 Capillary 0.12x170 (1.9 μ l)	170	0.12	1.9
ASM capillary	Generic Capillary	500	0.12	5.7

The 'ASM factor' is set to 1.7. 'Ok' and 'Cancel' buttons are at the bottom right.

Figure 81 ASM valve configuration (overview)

In this case, you need to enter two of following three parameters: length, diameter or volume. These parameters are required for calculating the flush volume and back pressure, see ["Understanding the ASM factor"](#) on page 41

The ASM factor is calculated and displayed based on selected capillaries.

- 4 Install capillary connections as displayed in figure **Valve topology** in the UI, see [Figure 79](#) on page 190.

NOTE

Please note that ASM capillaries are labeled with ASM (in contrast to transfer and other capillaries).

NOTE

Please note that port positions given in the MHC valve configuration in 2D-LC Software A.01.04 refer to standard 2D-LC valves, not ASM. Please use correct ports displayed figure **Valve topology**. This will be corrected in 2D-LC Software A.01.04 SR1.



7

Investigate the effects of using different gradients in the ²Dimension


Investigate the effects of using different gradients in the 2Dimension 193

This chapter describes, how shifted gradients in the second dimension can be used to enlarge the accessible two-dimensional separation space.

Investigate the effects of using different gradients in the 2nd Dimension

When combining separation systems with related separation mechanisms in the first and second dimension (as in RPxRP), orthogonality is limited. As a result, only a part of the available two-dimensional separation space will be occupied. In such a case, shifted gradients in the second dimension can be used to enlarge the accessible two-dimensional separation space.


- 1 To investigate the effects of using different gradients in the second dimension, firstly run a comprehensive 2D-LC separation with the same second dimension gradient from 5-95 % B repeated during the whole run. The 1D pump method should be set up as during the checkout runs (see below):




Investigate the effects of using different gradients in the 2Dimension

Investigate the effects of using different gradients in the 2Dimension

- 2 In **Instrument> Setup 2D-LC**, set up a 2D pump and modulation method with repeating gradients from 5 – 95 % B as shown below:



- 3 Run the comprehensive 2D-LC analysis. The resulting separation should look similar to the one shown below:




NOTE

Notice how the peaks are distributed around a diagonal line, indicating related separation mechanisms in the first and second dimension.


Investigate the effects of using different gradients in the 2Dimension

Investigate the effects of using different gradients in the 2Dimension

- 4 To improve the separation in the second dimension, a shallower second dimension gradient (e.g. from 25 – 75 % B) could be used. The setup of this 2D method is shown below (this is just shown for explanation purpose; you do not need to run this method!):



The separation resulting from using repeating gradients from 25 – 75 % B in the second dimension is shown below:



NOTE

Notice how the peaks are slightly further separated in the second dimension compared to using repeating gradients from 5 – 95 % B. Also notice that the last peaks eluting from the first dimension column are not eluted in one modulation cycle from the second dimension column (wrap-around; see marked area).

To be able to use even shallower gradients in the second dimension to further improve the separation and to also avoid the occurrence of wrap-around, continuously shifted gradients can be used in the second dimension (as was done during the checkout runs).

- 5 Compare the separations resulting from using the same second dimension gradient (from 5 – 95 % B and also from 25 – 75 % B) repeating during the whole run to the separation obtained using continuously shifted second dimension gradients in the checkout run.


NOTE

Notice how the peaks are spread more widely across the two-dimensional separation space (the accessible two-dimensional separation space is enlarged) when shifted gradients are used. Also, notice the effect that using continuously shifted second dimension gradients has on the second dimension retention times of consecutive fractions of the same first dimension peak.

Investigate the effects of using different gradients in the 2Dimension


Investigate the effects of using different gradients in the 2Dimension

- 6 Apart from using continuously shifted gradients in the second dimension, as was done during the checkout runs, it is also possible to stepwise shift the second dimension gradients. For this purpose, keep the valve & loop configuration as well as the 1D pump method the same. In **Instrument> Setup 2D-LC**, set up a 2D pump and modulation method with stepwise shifted gradients as shown below:



- 7 Run the comprehensive 2D-LC analysis with stepwise shifted gradients in the second dimension.

The resulting separation should look similar to the one shown below:



NOTE

Notice how consecutive fractions of the same first dimension peak have exactly the same retention time in the second dimension, as they experienced exactly the same second dimension gradient (in contrast to using continuously shifted gradients in the second dimension, which leads to consecutive fractions of one first dimension peak experiencing slightly different second dimension gradients). But be careful! This is only true if the stepwise shifting of the second dimension gradients is performed at times, when no peaks are eluting from the first dimension column.

In case your resulting separation looks different from the one shown above: Your peaks might show a different first dimension retention time due to the use of another first dimension pump (in the separation shown above, a 1290Bin Pump was used in the first dimension). Check whether the stepwise shifting of the second dimension gradients was performed at times when peaks eluted from the first dimension column in your separation and understand the effect this can have on the second dimension retention times of consecutive fractions of the same first dimension peak!

8

Data Analysis

Data Analysis for Heartcutting 2D-LC (LC-LC) 201

2D-LC Viewer 202

High-resolution sampling - Results 207

Quantitation 211

Data Analysis for Comprehensive 2D-LC (LCxLC) 219

Overview 222

Installation 222

Use GCImage Software 223

This chapter describes the analyzation of data in 2D-LC and is separated in a section heart-cutting 2D-LC and a section comprehensive 2D-LC.

Data Analysis for Heartcutting 2D-LC (LC-LC)

For data-analysis of heart-cutting 2D-LC data OpenLAB CDS ChemStation edition is usually fully sufficient.

Again, the data will be stored in one data-file. If more than one peak was analyzed in the second dimension they will simply follow one after the other in a distance of the second dimension run-time.

If a detector right after the first dimension column was used, e.g. as peak detector for peak triggered operation, these data will be available as a second data-trace.




Figure 82 Example of a heartcutting 2D-LC experiment. Data analysis and display done by OpenLAB CDS ChemStation Edition, the data of a detector placed after the first dimension column (blue) and of the 2D data detector (green). The additional impurities (marked red) could not be detected by 1D-LC only.

2D-LC Viewer

Using the Multiple Heart-Cutting Upgrade kit, the Agilent 1290 Infinity II 2D-LC Solution offers the possibility to store multiple peaks in several sample loops. These stored samples are then injected to the second dimension one by one. Thus long ²D gradients are possible without loss of ¹D peaks. But it is quite difficult to review the ²D results using the standard ChemStation Data Analysis. Especially as parked peaks are analyzed in a different order as they have been parked in (this is necessary to avoid carry-over). The **2D-LC Viewer** offers the opportunity to view and analyze second dimension chromatograms comfortably. The viewer can also be used for the analysis of standard heart-cutting 2D-LC data.

Overview 2D-LC Viewer



Figure 83 Overview of the 2D-LC Viewer graphical user interface

1	Tab 2D-LC Viewer
2	Separate displays tabs
3	¹ D Chromatogram(s)
4	Sampling table (¹ D)
5	Peak table (² D)
6	² D Chromatogram(s)
7	MS Spectra

The **2D-LC Viewer** provides the following functions:

- Tab pages enable the user to switch between
 - **2D-LC Viewer**, and
 - **Data Analysis**
- All panes are connected. Highlighting a cut or a chromatogram in one of the fields, will automatically highlight it in the other fields.
- **Sampling table (¹D)**
- **Peak table (²D)**
- Toolbar with the elements:
 - Print to printer
Prints the report according to the options set in the **Report Options** dialog using the standard print dialog.
 - Print preview
Shows the rendered report in a preview window. It is possible to print the report directly from the preview window.
 - Report options
Shows up the Report Options Dialog
 - Auto scale
Resets all chromatogram windows to their default scaling
 - **²D chromatogram** (checkbox)
Hide / unhide the full **²D Signal** window
 - **¹D Signal** list box
Used to select a signal from the ¹D detector
 - **²D Signal** list box
Used to select a signal from the ²D detector
- **¹D Chromatogram**
- **²D Chromatogram** (hidden, if **²D chromatogram** checkbox unchecked)
- **²D Chromatogram(s)**
- **MS Spectrum**

Sampling table (¹D)

The table lists all heart-cuts which have been analyzed in the 2nd dimension.

Cut #	¹ D Cut start [min]	Sampling time [s]	Mode	² D Run start [min]
Cut group: 1				
1	6.37	3.20	Time	6.45
Cut group: 2				
2	11.80	3.20	Time	20.09
3	11.85	3.20	Time	18.09
4	11.91	3.20	Time	16.09
5	11.96	3.20	Time	14.09
6	12.01	3.20	Time	12.09
7	12.07	3.20	Time	29.00
8	12.12	3.20	Time	27.00
9	12.17	3.20	Time	25.00
10	12.23	3.20	Time	23.00

Figure 84 Sampling table (¹D) (example)

The different columns can be selected or deselected by right mouseclicking the headline of the table.

Cut #	¹ D Cut start [min]	Sampling time [s]	Mode	² D Run start
Cut group: 1				
1	6.37	3.20	Time	6.45
Cut group: 2				
2	11.80	3.20	Time	20.09
3	11.85	3.20	Time	18.09
4	11.91	3.20	Time	16.09
5	11.96	3.20	Time	14.09
6	12.01	3.20	Time	12.09
7	12.07	3.20	Time	29.00
8	12.12	3.20	Time	27.00
9	12.17	3.20	Time	25.00
10	12.23	3.20	Time	23.00

Table 14 Legend for Sampling Table

PosNr	Description
Cut #	The current number of the heart-cut
¹ D Cut start [min]	Time when the heart-cut starts (peak begin or time value in trigger table)
¹ D Sampling time [min]	The duration (in minutes) of the heart-cut in the 1st dimension. The duration is determined either by the loop fill time, the end-of-peak detection or the max peak duration
Trigger	Indicates whether heart-cut was taken based on a peak-trigger (Peak) or based on a time given in the trigger table
² D run start [min]	Time when the analyses of this heart-cut in the 2nd dimension starts (gradient start)
¹ D Ret. Time [min]	The retention time (as given by the integrator) of the highest peak in the 1st dim. signal within heart-cut time range. The table cell is empty if no peak found or the signal isn't integrated. (Column not visible by default)
Deck	Number of the deck (1 or 2) where the cut (peak) has been parked / analyzed (Column not visible by default)
Loop	Number of the loop (1... 6) where the cut (peak) has been parked / analyzed (Column not visible by default)

¹D Chromatogram




Figure 85 ¹D Chromatogram (example)

Heart-cuts can be selected using left mouse button, multiple selection using **Ctrl-key** + left mouse button is also supported.

The selected signal (see toolbar) from the ¹D detector is shown.

- Heart cuts are indicated by a grey rectangle area
- Selected heart-cut(s) is (are) marked in a blue rectangle
- Heart-cuts are annotated using the retention time if available.
Otherwise the heart-cuts are annotated using their current number.
- Peaks (cuts) that couldn't be taken during acquisition will be marked by a warning icon on the x-axis at the time the heart-cut should have been taken.
- A tooltip provides more information about time and reason why the peak couldn't be cut.

2D Chromatogram




Figure 86 Full ²D Signal (example - only visible, if Full ²D Signal (checkbox) is checked)

This window shows the selected signal of the ²D detector containing the individual analyses of the heart-cuts.

- The selected heart-cut(s) is marked as a blue area.
- The area of a heart-cut is marked with a gray rectangle when hovering with the mouse over the chromatogram window.
- A heart-cut can be selected by clicking in such a rectangle.
- Multiple selections are supported using **Ctrl**-click.
- All heart-cuts are annotated using the heart-cut number (see also **Sampling table**).
- **F** indicates a bypass (or flush) gradient, which was used to flush the transfer capillaries after switching the 2D-LC valve.




Figure 87 2D Chromatogram (example from 2D Chromatogram above)

The 2D chromatogram of the selected heart-cut is shown as an individual run (x-axis starting at time 0). Chromatograms are overlaid if multiple heart-cuts are selected

For further details, refer to the online help.

High-resolution sampling - Results




Figure 88 Results peak group 1




Figure 89 Results peak group 2

The **1D chromatogram** shows 2 groups of peaks, which have been defined in the setup of a High-Resolution measurement. This information of compounds in **1D** belonging together will later be used for data analysis. The heart-cut table lists all cuts created in the first dimension, structured by signals and peak groups.

The **2D Chromatogram(s)** shows an overlay of all cuts, which have been selected in the heart-cut table. The Peak Table lists all peaks (peak areas), which have been found in the second dimension by applying ChemStation integrator settings.

Areas of peaks which belong together are summed up automatically.

Structure of Sampling and Peak Tables

Sampling table (1D)				
Cut #	1D Cut start [min]	Sampling time [min]	Mode	1D Run start [min]
1	9.14	0.08	Time	26.50
2	9.24	0.08	Time	22.30
3	9.35	0.08	Time	18.10
4	9.46	0.08	Time	13.90
5	9.56	0.08	Time	9.70
6	13.14	0.08	Time	49.36
7	13.24	0.08	Time	45.16
8	13.34	0.08	Time	40.96
9	13.44	0.08	Time	36.76
10	13.54	0.08	Time	32.56

Figure 90 Sampling table

1	First level: All signals for all runs
2	Second level: Peak groups
3	Third level: Cuts numbered continuously per signal

Peak table (2D)				
Cut # /	2D Retention time [min]	Area	Height	Width
Signal: MSD1 213, EIC=212.8:213.8				
Component:5 100.000 Area%				
3	2.781	1826339.0	968824.7	0.028
		1826339.0		
Signal: MSD1 259, EIC=258.7:259.7				
Component:8 100.000 Area%				
4	2.916	380586.9	219993.6	0.028
		380586.9		
Signal: MSD1 TIC, MS File				
Component:3 49.514 Area%				
3	2.781	2915092.3	1537154.3	0.029
		2915092.3		
Component:4 2.960 Area%				
3	2.919	174246.7	92484.4	0.028
		174246.7		
Component:6 6.068 Area%				
4	2.780	357227.6	175817.6	0.030
		357227.6		
Component:7 16.949 Area%				
4	2.917	997879.5	535319.3	0.028

Figure 91 Peak Table

1	First level: All signals for all runs
2	Second level: Components with matching 2D retention times
3	Third level: Cuts and area sum per compound

Summed signal

The optional summed signal (black line) shows the sum of all selected cuts.

This allows easy navigation: one area of interest in ¹D corresponds to one summed 2D chromatogram (for trying, choose option **Show summed signal** only).

For a single wavelength, this corresponds to the entire ²D area used for quantitation.




Figure 92 Summed signal

Quantitation

Compared to a standard quantitation, a 2D-LC quantitation requires 2 additional steps (identification of compounds, summation of areas) carried out by the 2D-LC software.





Figure 93 Necessary steps for quantitation in 2D-LC

Identify compounds

This procedure exemplarily shows a quantitative measurement of the 2D-LC checkout sample using a 3-level calibration of two compounds at two wavelengths.



- 1 Select signal.

NOTE

Use a signal with a low concentration.

Sequence: ZDLC_Quant_2016-04-29

Ready/Reprocess Data Mode



Overlay	Type	Line	Inj	Vial	Sample Name	Acq. Method	Sequence Method	Sample Type	Data File
▶		1	1	P2-E1	ZDLC df 500	Checkout_HiRes_12...	Checkout_HiRes_12...	Sample	001-P2-E1-...
+		2	1	P2-E2	ZDLC df 200	Checkout_HiRes_12908...	Checkout_HiRes_12908...	Sample	002-P2-E2-2...
+		3	1	P2-E3	ZDLC df 100	Checkout_HiRes_12908...	Checkout_HiRes_12908...	Sample	003-P2-E3-2...
+		4	1	P2-E4	ZDLC df 50	Checkout_HiRes_12908...	Checkout_HiRes_12908...	Sample	004-P2-E4-2...
+		5	1	P2-E6	ZDLC df 10	Checkout_HiRes_12908...	Checkout_HiRes_12908...	Sample	005-P2-E6-2...
+		7	1	P2-C3	Sample 20	Checkout_HiRes_12908...	Checkout_HiRes_12908...	Sample	007-P2-C3-5...

- 2 Integrate the cuts.


NOTE

Integration is done using the ChemStation integrator.

- a Click on  to **Set Integration Range**.




The 2D-LC Viewer offers simplified access for changing the integration of one or multiple limits.



3 Verify small peaks in 2D. Use Full 2D Signal or 2D run start time for navigation.

NOTE


Peak integration is based on ChemStation integrator settings. Please verify, if all peaks of interest have been found and integrated. If not, ChemStation settings need to be adjusted. For finding the corresponding peak in the 2D signal, use the peak table and heart-cut table.




4 Adjust integrator settings in ChemStation.

NOTE

ChemStation has no information about cuts and individual ²D chromatograms. So it is important to know, where to look at. In the example above, Cut #2 starts at 16.97 minutes. At 17.939 minutes is compound 2 of cut #2, which is not integrated. In this example, decreasing the area reject parameter from 5 to 1 enables integrating the peak.



5 Apply the adjusted integrator settings to other signals.

**NOTE**

Specific Events For Signal:
DAD2 B Specific
DAD Default
DAD2 C Specific

In integration events, change to a different wavelength for applying settings here. Consider using same limits for other signals, apply to sequence method.


Calibration

Create a calibration table


Sequence: 2D Training USA 2015-08-04 11-16-08

Overlay	Type	Line	Inj	Vial	Sample Name	Acq. Method	Sample Type	Data File
+		1	1	P1-A1	2D-LC Checkout	Checkout_MHC_Qua...	Sample	sLCxLC Quant 1.D
+		2	1	P1-A1	2D-LC Checkout	Checkout_MHC_Quant.M	Sample	sLCxLC Quant 2.D
+		3	1	P1-A1	2D-LC Checkout	Checkout_MHC_Quant.M	Sample	sLCxLC Quant 5.D
+		4	1	P1-A1	2D-LC Checkout	Checkout_MHC_Quant.M	Sample	sLCxLC Quant 3.D

- 1 Create a new calibration table from the MHC Viewer.




- 2 Enter amount/concentration.



Add levels

Overlay	Type	Line	Inj	Vial	Sample Name	Acq. Method	Sample Type	Data File
		1	1	P1-A1	2D-LC Checkout	Checkout_MHC_Quant.M	Sample	sLCxLC Quant 1.D
+		2	1	P1-A1	2D-LC Checkout	Checkout_MHC_Qua...	Sample	sLCxLC Quant 2.D
+		3	1	P1-A1	2D-LC Checkout	Checkout_MHC_Quant.M	Sample	sLCxLC Quant 5.D
+		4	1	P1-A1	2D-LC Checkout	Checkout_MHC_Quant.M	Sample	sLCxLC Quant 3.D

- 1 Select and add more signals and levels as needed.



- 2 Set concentration for each level.

Add new level



Default amount ng/μl

OK Cancel

Calibration table

Calibration Table						
Enter Delete Insert... Print OK Help						
#	RT	Signal	Compound	Lvl	Amt[ng/ul]	Area
1	12.617	DAD2 A	Compound1	1	1.000	228.780
				2	2.000	464.650
				3	5.000	1172.200
	12.617	DAD2 B		1	1.000	69.095
				2	2.000	141.200
				3	5.000	355.910
2	12.686	DAD2 A	Compound2	1	1.000	365.560
				2	2.000	744.530
				3	5.000	1901.500
	12.686	DAD2 B		1	1.000	206.250
				2	2.000	419.610
				3	5.000	1065.300

Peak Table						
Cut #	/	2D Retention time (rel.) [min]	Area	Height	Width	Symmetry
5		0.901	320.121	336.944	0.015	0.906
6		0.901	210.633	226.007	0.014	0.907
7		0.900	138.931	147.254	0.015	0.897
8		0.897	67.435	71.419	0.015	0.902
9		0.896	28.679	30.527	0.015	0.902
10		0.896	11.465	12.105	0.015	0.890
			1172.219			






Quantification

Quantifying samples with unknown concentration based on the calibration table can be done using ChemStation reporting as usual.

Quantification

Overlay	Type	Line	Inj	Vial	Sample Name	Acq. Method	Sample Type	Data File
		1	1	P1-A1	2D-LC Checkout	Checkout_MHC_Quant.M	Sample	sLCxLC Quant 1.D
		2	1	P1-A1	2D-LC Checkout	Checkout_MHC_Quant.M	Sample	sLCxLC Quant 2.D
		3	1	P1-A1	2D-LC Checkout	Checkout_MHC_Quant.M	Sample	sLCxLC Quant 5.D
		4	1	P1-A1	2D-LC Checkout	Checkout_MHC_Qua...	Sample	sLCxLC Quant 3.D



Signal 1: DAD2 A, Sig=254,4 Ref=360,100

RetTime	Type	Area	Amt/Area	Amount	Grp	Name
[min]		[mAU*s]		[ng/ul]		
12.617	BB	705.28630	4.27790e-3	3.01715		Compound1
12.686	BBA	1134.78678	2.64492e-3	3.00142		Compound2

Totals : 6.01857

Data Analysis for Comprehensive 2D-LC (LCxLC)

Typically very complex samples are analyzed by comprehensive 2-dimensional liquid chromatography. The compounds which are often co-eluting from the first dimension are further separated in the second dimension. With the Agilent 1290 Infinity II 2D-LC Solution always one large data-file spanning the run-time of the two-dimensional analysis will be acquired. As an example, a 2-dimensional analysis of a mixture of 26 polyphenolic standard compounds is shown in a one dimensional data analysis display (Figure 94 on page 220). Theoretically, the data can be analyzed with OpenLAB CDS ChemStation edition software.

But for easier data-analysis and better visualization of the comprehensive 2D-LC data special software is recommended. Agilent recommends GC Image LCxLC edition Software from GC Image LLC, Nebraska, USA. A trial download can be found on www.GCImage.com as well as an online manual. Agilent 2D-LC data files also including UV spectra and mass spectra data can be directly imported. This software, with the information of the modulation time, is capable to extract the data and isolate each second dimension run. Data will be reconstructed in a two-dimensional display of the retention times. This can be displayed as a colored 2-dimensional map of compound peaks (Figure 95 on page 220). After baseline correction the peaks can be automatically detected by a peak detection algorithm inherent in the 2D-LC data analysis software (Figure 96 on page 221). Since the third dimension is the intensity of the peaks a 3-dimensional plot of the data is possible (Figure 97 on page 221). With the given data set further qualitative and quantitative data analysis is possible.




Figure 94 Display of two-dimensional LC data with a one-dimensional data analysis software




Figure 95 2D-LC plot of the optimized separation of 26 polyphenolic compounds




Figure 96 2DLC plot after baseline correction and with software detected peak annotation




Figure 97 3-Dimensional display of the separation of the 26 compound standard mixture. The first dimension separation takes 40 minutes and each second dimension separation takes 39 seconds. The back side shows a generated first dimension chromatogram and gives the impression which peaks are coeluting and separated in the second dimensions.

Overview

GC Image LC x LC Edition (short GC Image) is a software for for visualization and data analysis of full comprehensive two-dimensional liquid chromatograms:

- M8700AA GC Image LCxLC Edition for UV and Single Quad measurements
- M8710AA GC Image LCxLC-HRMS Edition for UV and/or High Resolution MS measurements (Q-TOF)

Installation

Parts required

Description

CD with software

License dongle (Wibu Key)

Activation code

- 1 The CD contains two executables: LCxLC2.8r3-HRMS.exe (or higher), LCxLC2.8r3-HRMS-64.exe (or higher). Choose the appropriate version for your operating system. Corresponding versions are available for the UV only detection.
- 2 Double-click the chosen executable and follow the instructions on the screen.
- 3 Activate the software with the USB key. Insert the USB dongle and wait. The driver will install automatically.
- 4 Activate R2.8 (or higher) in the Windows Start Menu.
- 5 Enter the activation code, which is shipped with the software.

Use GCImage Software

GCImage is a powerful expert software with many sophisticated features for display, data analysis, compound identification, library search, workflow automation, reporting etc.

The basic knowledges to successfully use the software are the following:

- Import 2D ChemStation data files
- Setting the modulation period
- Choosing a color mapping
- Navigate in the display
- Navigate in the display
- Detect peaks (Blobs)

Preparations


The USB dongle needs always to be inserted when working with GCImage software. If not, you will be asked to insert it.

Basic knowledges


1 Start up LCImage

LCImage offers optionally a password protected user management system. If you don't need it, simply click „Login with system“, which is based on Windows user account.



- 2 Import the UV signal from the second dimension detector.




3 Import parameters




4 Fit view



5 Correct Baseline




6 Shift phase




7 Zoom into an interesting region by using the right mouse button and dragging over the display


- 8 Adjust colors: LC Image offers refined possibilities for optimizing the color scales. Play around with settings for improving the contrast.




- 9 Select a data range.




- 10 By clicking the „Show 3D perspective“ button or the corresponding menu item, you can easily create a customizable 3D plot.



11 View single 2D chromatograms




12 Select blobs



MS Data

- 1 Import MS data: The import functionality of MS data is very similar to those of UV measurements. Additionally, you can for example filter to a certain mass range („range limit“), that you are interested in.



- 2 By clicking on „Show 1D view“, you can display the TIC for that 2D slice.
- 3 By clicking on data points or blobs in the 2D view, you can display MS spectra of corresponding plots.

9

Troubleshooting and Diagnostics

Overview of the Module's Indicators and Test Functions	232
User Interfaces	233
Agilent Lab Advisor Software	234
The Basic Principle of Troubleshooting	235
Pressure too high	238
Pressure too low	239
Peak area and peak height related	240
Retention time related	241
Missing signal linearity	242
Drifting signal	243
Signal noisy	244
Recommended Tests to Conclude Troubleshooting	245

This chapter gives an overview about the troubleshooting and diagnostic features and the different user interfaces.

Overview of the Module's Indicators and Test Functions

For an overview of the module's indicators and test functions, refer to the manuals of the modules installed in your system.

User Interfaces

- Depending on the user interface, the available tests and the screens/reports may vary.
- Preferred tool should be Agilent Lab Advisor Software, see [“Agilent Lab Advisor Software”](#) on page 234.
- The Agilent OpenLAB ChemStation C.01.03 and above do not include any maintenance/test functions.
- Screenshots used within these procedures are based on the Agilent Lab Advisor Software.

Agilent Lab Advisor Software

The Agilent Lab Advisor Software (basic license, shipped with an Agilent LC pump) is a standalone product that can be used with or without a chromatographic data system. Agilent Lab Advisor helps to manage the lab for high-quality chromatographic results by providing a detailed system overview of all connected analytical instruments with instrument status, Early Maintenance Feedback counters (EMF), instrument configuration information, and diagnostic tests. By the push of a button, a detailed diagnostic report can be generated. Upon request, the user can send this report to Agilent for a significantly improved troubleshooting and repair process.

The Agilent Lab Advisor software is available in two versions:

- Lab Advisor Basic
- Lab Advisor Advanced

Lab Advisor Basic is included with every Agilent 1200 Infinity Series and Agilent InfinityLab LC Series instrument.

The Lab Advisor Advanced features can be unlocked by purchasing a license key, and include real-time monitoring of instrument actuals, all various instrument signals, and state machines. In addition, all diagnostic test results, calibration results, and acquired signal data can be uploaded to a shared network folder. The Review Client included in Lab Advisor Advanced allows to load and examine the uploaded data no matter on which instrument it was generated. This makes Data Sharing an ideal tool for internal support groups and users who want to track the instrument history of their analytical systems.

The optional Agilent Maintenance Wizard Add-on provides an easy-to-use, step-by-step multimedia guide for performing preventive maintenance on Agilent 1200 Infinity and Agilent InfinityLab LC Series instrument.

The tests and diagnostic features that are provided by the Agilent Lab Advisor software may differ from the descriptions in this manual. For details, refer to the Agilent Lab Advisor software help files.

The Basic Principle of Troubleshooting

Troubleshooting key Concept – Divide and Conquer

The following troubleshooting concept, shows exemplarily how to approach problems in 2D-LC chromatography.


Most of the following explanations can also be used to isolate and detect standard LC issues.

The basic principle of troubleshooting should always be a step by step approach to the 2D-LC problem. As a first step, find out whether the cause of the error is either:

- The application method, or
- The 2D-LC instrument

For a recommended approach to isolate the cause of the issue, see the graphic below. All examples use symbols as described in the following table.

Table 15 Description for symbols as used in troubleshooting decision trees

Symbol	Description
	Shows and describes a problem in the 2D-LC system. Indicates the starting point for a series of actions and decisions leading to a solution for the problem.
	Illustrates, that the user must identify what an observation means. Then the user must take a decision, which further way of troubleshooting to follow.
	Shows, the user must act to proceed and come to the next decision or solution.




Figure 98 Example for a strategy to identify the application method or instrument as root cause for issues in 2D-LC chromatography

After ruling out the application method as the cause of the issue, one can start to search for the problem's root cause within the 2D-LC Instrument hardware.

Common HPLC hardware issues, along with the location of each problem's respective troubleshooting procedure are listed below:

- "Pressure too high" on page 238
- "Pressure too low" on page 239
- "Peak area and peak height related" on page 240
- "Retention time related" on page 241
- "Missing signal linearity" on page 242
- "Drifting signal" on page 243
- "Signal noisy" on page 244

Pressure too high




Figure 99 Example for a strategy to eliminate issues related to too high pressure in 2D-LC instruments

Pressure too low




Figure 100 Example for a strategy to eliminate issues related to too low pressure in 2D-LC instruments

Peak area and peak height related




Figure 101 Example for a strategy to eliminate issues related to peak problems in 2D-LC instruments

Retention time related




Figure 102 Example for a strategy to eliminate issues related to retention time in 2D-LC instruments

Missing signal linearity




Figure 103 Example for a strategy to eliminate issues related to missing signal linearity in 2D-LC instruments

Drifting signal




Figure 104 Example for a strategy to eliminate issues related to drifting signal in 2D-LC instruments

Signal noisy




Figure 105 Example for a strategy to eliminate issues related to signal noise in 2D-LC instruments

Recommended Tests to Conclude Troubleshooting

The following table shows the most important tests to conclude troubleshooting.

- For further detailed information, see:
 - Maintenance information in the specific manual of each module.
 - *Troubleshooting Guide* poster 5994-0709EN.
 - *Best Practice for Using an Agilent LC System* 01200-90090.
- For additional help, contact your local Agilent Technologies service representative.

Table 16 Recommended Tests for 2D-LC System Troubleshooting

Pump	Column Com- partment	Autosampler	Valve	Detector	2D-LC Instru- ment
Pressure Test Leak Test	Thermostat Test Pressure Test (if column valve is pres- ent)	Pressure Test Inject stan- dards or inject different volu- mina or blanks	Switching valve posi- tion/Check pressure read- ing Pressure Test	Lamp Intensity Test Wavelength calibration In addition there are detector spe- cific tests.	Run Checkout For 2D-LC Instruments • Pressure test of the 1D-LC Part • Pressure Test of the 2D-LC Part
Pump charac- teristic • Pump Rip- ple (1260 Pump) • Tuning (1290 pump)					

What Are Error Messages	247
General Error Messages	248
Timeout	248
Shutdown	249
Remote Timeout	250
Lost CAN Partner	250
Leak Sensor Short	251
Leak Sensor Open	251
Compensation Sensor Open	252
Compensation Sensor Short	252
Fan Failed	253
Leak	253
Module-Specific Error Messages	254
Initialization of Valve Failed	254
Valve Switching Failed	255
Valve Tag Violation	255
Pressure Cluster Partner Missing	256
Position Cluster Partner Missing	256
External Valve falls into resident mode	257

This chapter describes the meaning of error messages, and provides information on probable causes and suggested actions how to recover from error conditions.

What Are Error Messages

Error messages are displayed in the user interface when an electronic, mechanical, or hydraulic (flow path) failure occurs which requires attention before the analysis can be continued (for example, repair, or exchange of consumables is necessary). In the event of such a failure, the red status indicator at the front of the module is switched on, and an entry is written into the module logbook.

If an error occurs outside a method run, other modules will not be informed about this error. If it occurs within a method run, all connected modules will get a notification, all LEDs get red and the run will be stopped. Depending on the module type, this stop is implemented differently. For example, for a pump the flow will be stopped for safety reasons. For a detector, the lamp will stay on in order to avoid equilibration time. Depending on the error type, the next run can only be started, if the error has been resolved, for example liquid from a leak has been dried. Errors for presumably single time events can be recovered by switching on the system in the user interface.

Special handling is done in case of a leak. As a leak is a potential safety issue and may have occurred at a different module from where it has been observed, a leak always causes a shutdown of all modules, even outside a method run.

In all cases, error propagation is done via the CAN bus or via an APG/ERI remote cable (see documentation for the APG/ERI interface).

General Error Messages

General error messages are generic to all Agilent series HPLC modules and may show up on other modules as well.

Timeout

Error ID: 0062

The timeout threshold was exceeded.

Probable cause	Suggested actions
1 The analysis was completed successfully, and the timeout function switched off the module as requested.	Check the logbook for the occurrence and source of a not-ready condition. Restart the analysis where required.
2 A not-ready condition was present during a sequence or multiple-injection run for a period longer than the timeout threshold.	Check the logbook for the occurrence and source of a not-ready condition. Restart the analysis where required.

Shutdown

Error ID: 0063

An external instrument has generated a shutdown signal on the remote line.

The module continually monitors the remote input connectors for status signals. A LOW signal input on pin 4 of the remote connector generates the error message.

Probable cause	Suggested actions
1 Leak detected in another module with a CAN connection to the system.	Fix the leak in the external instrument before restarting the module.
2 Leak detected in an external instrument with a remote connection to the system.	Fix the leak in the external instrument before restarting the module.
3 Shut-down in an external instrument with a remote connection to the system.	Check external instruments for a shut-down condition.
4 The degasser failed to generate sufficient vacuum for solvent degassing.	Check the vacuum degasser for an error condition. Refer to the <i>Service Manual</i> for the degasser or the pump that has the degasser built-in.

Remote Timeout

Error ID: 0070

A not-ready condition is still present on the remote input. When an analysis is started, the system expects all not-ready conditions (for example, a not-ready condition during detector balance) to switch to run conditions within one minute of starting the analysis. If a not-ready condition is still present on the remote line after one minute the error message is generated.

Probable cause	Suggested actions
1 Not-ready condition in one of the instruments connected to the remote line.	Ensure the instrument showing the not-ready condition is installed correctly, and is set up correctly for analysis.
2 Defective remote cable.	Exchange the remote cable.
3 Defective components in the instrument showing the not-ready condition.	Check the instrument for defects (refer to the instrument's documentation).

Lost CAN Partner

Error ID: 0071

During an analysis, the internal synchronization or communication between one or more of the modules in the system has failed.

The system processors continually monitor the system configuration. If one or more of the modules is no longer recognized as being connected to the system, the error message is generated.

Probable cause	Suggested actions
1 CAN cable disconnected.	• Ensure all the CAN cables are connected correctly. • Ensure all CAN cables are installed correctly.
2 Defective CAN cable.	Exchange the CAN cable.
3 Defective mainboard in another module.	Switch off the system. Restart the system, and determine which module or modules are not recognized by the system.

Leak Sensor Short

Error ID: 0082

The leak sensor in the module has failed (short circuit).

The current through the leak sensor is dependent on temperature. A leak is detected when solvent cools the leak sensor, causing the leak sensor current to change within defined limits. If the current increases above the upper limit, the error message is generated.

Probable cause	Suggested actions
1 Defective leak sensor.	Please contact your Agilent service representative.
2 Leak sensor incorrectly routed, being pinched by a metal component.	Please contact your Agilent service representative.

Leak Sensor Open

Error ID: 0083

The leak sensor in the module has failed (open circuit).

The current through the leak sensor is dependent on temperature. A leak is detected when solvent cools the leak sensor, causing the leak-sensor current to change within defined limits. If the current falls outside the lower limit, the error message is generated.

Probable cause	Suggested actions
1 Leak sensor not connected to the main board.	Please contact your Agilent service representative.
2 Defective leak sensor.	Please contact your Agilent service representative.
3 Leak sensor incorrectly routed, being pinched by a metal component.	Please contact your Agilent service representative.

Compensation Sensor Open

Error ID: 0081

The ambient-compensation sensor (NTC) on the main board in the module has failed (open circuit).

The resistance across the temperature compensation sensor (NTC) on the main board is dependent on ambient temperature. The change in resistance is used by the leak circuit to compensate for ambient temperature changes. If the resistance across the sensor increases above the upper limit, the error message is generated.

Probable cause	Suggested actions
1 Defective main board.	Please contact your Agilent service representative.

Compensation Sensor Short

Error ID: 0080

The ambient-compensation sensor (NTC) on the main board in the module has failed (open circuit).

The resistance across the temperature compensation sensor (NTC) on the main board is dependent on ambient temperature. The change in resistance is used by the leak circuit to compensate for ambient temperature changes. If the resistance across the sensor falls below the lower limit, the error message is generated.

Probable cause	Suggested actions
1 Defective main board.	Please contact your Agilent service representative.

Fan Failed

Error ID: 0068

The cooling fan in the module has failed.

The hall sensor on the fan shaft is used by the main board to monitor the fan speed. If the fan speed falls below a certain limit for a certain length of time, the error message is generated.

Depending on the module, assemblies (e.g. the lamp in the detector) are turned off to assure that the module does not overheat inside.

Probable cause	Suggested actions
1 Fan cable disconnected.	Please contact your Agilent service representative.
2 Defective fan.	Please contact your Agilent service representative.
3 Defective main board.	Please contact your Agilent service representative.

Leak

Error ID: 0064

A leak was detected in the module.

The signals from the two temperature sensors (leak sensor and board-mounted temperature-compensation sensor) are used by the leak algorithm to determine whether a leak is present. When a leak occurs, the leak sensor is cooled by the solvent. This changes the resistance of the leak sensor which is sensed by the leak-sensor circuit on the mainboard.

Probable cause	Suggested actions
1 Loose fittings.	Ensure all fittings are tight.
2 Broken capillary.	Exchange defective capillaries.

Module-Specific Error Messages

For further module-specific errors, please see the manual of the module in question.

Initialization of Valve Failed

Error ID: 24000

During the initialization process the motor of the valve drive moves to some special positions depending on the installed valve head. A failure in this process means either that the movement couldn't be performed properly or it was not noticed correctly by the sensor.

Probable cause	Suggested actions
1 Mechanical problems. Friction too high or blockages on the valve drive's motor or on the valve head.	• Check valve head for correct installation • Try to identify the source of trouble by installing a different valve head if possible. • Contact your Agilent Service representative.
2 Defect Sensor on the Valve Drive Motor	• Check valve head for correct installation • Try to identify the source of trouble by installing a different valve head if possible. • Contact your Agilent Service representative.

Valve Switching Failed

Error ID: 24001

The valve drive was not able to operate the valve head correctly. Either due to mechanical reasons or the movement couldn't be detected correctly.

Probable cause	Suggested actions
1 Mechanical problems. Friction too high or blockages on the valve drive's motor or on the valve head.	• Check valve head for correct installation • Try to identify the source of trouble by installing a different valve head if possible. • Contact your Agilent Service representative.
2 Defect Sensor on the Valve Drive Motor	• Check valve head for correct installation • Try to identify the source of trouble by installing a different valve head if possible. • Contact your Agilent Service representative.

Valve Tag Violation

Error ID: 24006

The valve drive identified a different valve head than it had identified during the last initialization.

Probable cause	Suggested actions
1 A valve head has been exchanged (hot-plugged) while the valve drive was still powered on.	Change the valve head. It is important to have the valve switched off for at least 10 s after or before a new valve head has been installed.

NOTE

Soft power-down power supply of the valve drive.

Whenever you want to power cycle the valve drive for a re-boot, it needs to be powered off for at least 10 seconds.

Pressure Cluster Partner Missing

The connection from the valve drive to a defined pressure cluster partner is lost.

Probable cause	Suggested actions
1 Communication issues	Check the CAN cable connections of the modules.
2 Configuration mismatch	Check and correct if necessary the valve configuration and presence of defined pressure cluster partner.

Position Cluster Partner Missing

Probable cause	Suggested actions
1 Communication issues	Check the CAN cable connections of the modules.
2 Configuration mismatch	Check and correct if necessary the valve configuration and presence of defined position cluster partner.

External Valve falls into resident mode

Error ID: Flashing status indicator

The valve drive was not able to operate correctly

Probable cause	Suggested actions
1 Communication issues	• Check the CAN cable connections of the modules. • Check if the hosted module is present.
2 Configuration mismatch	• Check if the firmware on the entire stack is out of the same firmware set. • Check if the limit of 3 hosted modules for each host module is not exceeded. • Check if the dipswitch settings are correct. • Check if the firmware on the entire stack has to be the latest version.

Introduction to Maintenance	259
Warnings and Cautions	260
Overview of Maintenance	262
Cleaning the Module	263
Correcting Leaks	264
Correcting Leaks (G7116B)	264
Correcting Leaks (G1170A)	265
Replace Valve Heads	266
Replace Valve Heads (G7116B)	266
Replace Valve Heads (G1170A)	269
Replacing Parts of the Valve Head	272
Replacing the Fuses of the Infinity Valve Drive	274
Replace the Module Firmware	276

This chapter describes the maintenance of the 2D-LC Solution.

Introduction to Maintenance

The 2D-LC solution is designed for easy maintenance. The most frequent maintenance can be done from the front with the modules in place in the system stack. Examples are maintenance of the needle, needle seats, rotor seals, valve heads, or replacing heat exchangers.

Warnings and Cautions

WARNING**Personal injury or damage to the product**

Agilent is not responsible for any damages caused, in whole or in part, by improper use of the products, unauthorized alterations, adjustments or modifications to the products, failure to comply with procedures in Agilent product user guides, or use of the products in violation of applicable laws, rules or regulations.

- ✓ Use your Agilent products only in the manner described in the Agilent product user guides.
-

WARNING**Electrical shock**

Repair work at the module can lead to personal injuries, e.g. shock hazard, when the cover is opened.

- ✓ Do not remove the cover of the module.
 - ✓ Only certified persons are authorized to carry out repairs inside the module.
-

WARNING**Sharp metal edges**

Sharp-edged parts of the equipment may cause injuries.

- ✓ To prevent personal injury, be careful when getting in contact with sharp metal areas.
-

WARNING

Toxic, flammable and hazardous solvents, samples and reagents

The handling of solvents, samples and reagents can hold health and safety risks.

- ✓ When working with these substances observe appropriate safety procedures (for example by wearing goggles, safety gloves and protective clothing) as described in the material handling and safety data sheet supplied by the vendor, and follow good laboratory practice.
- ✓ The volume of substances should be reduced to the minimum required for the analysis.
- ✓ Do not operate the instrument in an explosive atmosphere.

CAUTION

Hot heat exchangers 

The column compartment has two heat exchanger assemblies that might be hot.

- ✓ Allow them to cool down before starting repairs.

CAUTION

Safety standards for external equipment

- ✓ If you connect external equipment to the instrument, make sure that you only use accessory units tested and approved according to the safety standards appropriate for the type of external equipment.

Overview of Maintenance

The following pages describe maintenance procedures (simple repairs) that can be done without opening the main cover.

Table 17 Maintenance procedures

Procedure	Typical Frequency	Notes
Cleaning the Module	If required	
Correcting Leaks	If a leak has occurred	Check for leaks
Maintain the Column Switching Valve	If valve leaks	
Replace Valve Heads	If the valve performance shows indication of leakage or wear	
Replacing Parts of the Valve Head	If leak sensor is defective	
Replacing the Fuses of the Infinity Valve Drive	When a fuse is defective	
Replace the Module Firmware	If required	

Cleaning the Module

To keep the module case clean, use a soft cloth slightly dampened with water, or a solution of water and mild detergent.

WARNING

Liquid dripping into the electronic compartment of your module can cause shock hazard and damage the module

- ✓ Do not use an excessively damp cloth during cleaning.
 - ✓ Drain all solvent lines before opening any connections in the flow path.
-

Correcting Leaks

Correcting Leaks (G7116B)

When	If a leakage has occurred at the heat exchanger or at the capillary connections or at the column switching valve.				
Tools required	<table><tr><td>Description</td></tr><tr><td>Tissue</td></tr><tr><td>Pipette</td></tr><tr><td>Wrench, 1/4 – 5/16 inch (for capillary connections)</td></tr></table> 1 Remove the door. 2 Use a pipette and tissue to dry the leak sensor area. 3 Observe the capillary connections and the column switching valve for leaks and correct, if required. 4 Re-install the door.	Description	Tissue	Pipette	Wrench, 1/4 – 5/16 inch (for capillary connections)
Description					
Tissue					
Pipette					
Wrench, 1/4 – 5/16 inch (for capillary connections)					

Correcting Leaks (G1170A)

When If leakage has occurred at the capillary connections or at the valve.

Tools required **Description**

Tissue

Pipette

Wrench, 1/4 – 5/16 inch
(for capillary connections)

- 1** Use a pipette and tissue to dry the leak sensor area.
- 2** Observe the capillary connections and the valve for leaks and correct, if required.

Replace Valve Heads

Replace Valve Heads (G7116B)

Several optional valve heads are available, which can be installed and exchanged easily.

Parts required**Description**

Any Agilent Quick Change Valve Head.

CAUTION

The valve actuator contains sensitive optical parts, which need to be protected from dust and other pollution. Pollution of these parts can impair the accurate selection of valve ports and therefore bias measurement results.

- ✓ Always install a valve head for operation and storage. For protecting the actuator, a dummy valve head (part of Transportation Lock Kit (G1316-67001)) can be used instead of a functional valve. Do not touch parts inside the actuator.

CAUTION**Column Damage or Bias Measurement Results**

Switching the valve to a wrong position can damage the column or bias measurement results.

- ✓ Fit the lobe to the groove to make sure the valve is switched to the correct position.

CAUTION**Valve Damage**

Using a low pressure valve on the high pressure side can damage the valve.

- ✓ When using multiple column compartments as part of a method development solution, make sure that the high pressure valve head is connected to the autosampler and the low pressure valve head is connected to the detector.

WARNING

Toxic, flammable and hazardous solvents, samples and reagents


The handling of solvents, samples and reagents can hold health and safety risks.

- ✓ Be sure that no solvent can drop out of the solvent connections when removing them from your valve head.
- ✓ When working with these substances observe appropriate safety procedures (for example by wearing goggles, safety gloves and protective clothing) as described in the material handling and safety data sheet supplied by the vendor, and follow good laboratory practice.


1 Switch off the module.




2 Push the valve head to bring it to its outer position.



3 Remove all capillary connections from the valve head.




4 Unscrew the valve head.




Maintenance

Replace Valve Heads


- 5 Put the new valve head onto the valve drive such that the lobe fits to the groove.




- 6 Screw the valve head onto the valve drive using the union nut.




- 7 Install all required capillary connections to the valve.



- 8 Push the valve head until it snaps in and stays in the rear position.



- 9 Switch on the module.



Replace Valve Heads (G1170A)

The following procedure shows installation only. To remove the valve, follow the instructions in reverse order.

NOTE

The following procedure exemplarily shows a valve head installation. For correct capillary connections see **Valve topology** in the GUI.

CAUTION

The valve actuator contains sensitive optical parts, which need to be protected from dust and other pollution. Pollution of these parts can impair the accurate selection of valve ports and therefore bias measurement results.

- ✓ **Always install a valve head for operation and storage. For protecting the actuator, a dummy valve head can be used instead of a functional valve. Do not touch parts inside the actuator.**

NOTE

For a correct installation of the valve head, the outside pin (red) must completely fit into the outside groove on the valve drive's shaft (red). A correct installation is only possible if the two pins (green and blue) on the valve head fit into their corresponding grooves on the valve drive's actuator axis. Their match depends on the diameter of the pin and groove.


NOTE

The tag reader reads the valve head properties from the valve head RFID tag during initialization of the module. Valve properties will not be updated, if the valve head is replaced while the module is on. Selection of valve port positions can fail, if the instrument does not know the properties of the installed valve.

NOTE


To allow correct valve identification, power off the module for at least 10 s.

- 1 Insert the valve head into the valve shaft.




OR

If the outside pin does not fit into the outside groove, you have to turn the valve head until you feel that the two pins snap into the grooves. Now you should feel additional resistance from the valve drive while continuously turning the valve head until the pin fits into the groove.




- 2 When the outer pin is locked into the groove, manually screw the nut onto the valve head.



NOTE

Fasten the nut with the 5043-1767 Valve Removal tool.

3 Install all required capillary connections to the valve.



4 Power on or power-cycle your module, so the valve head gets recognized during module initialization.

Replacing Parts of the Valve Head

When If valve leaks.

Tools required **Description**

Wrench, 1/4 inch
OR
Hexagonal key, 9/64 inch

- 1 Remove capillaries from ports.
- 2 Loosen each fixing stator screw two turns at a time. Remove bolts from head.
- 3 Remove the stator head (and stator face if applicable).
- 4 Remove the stator ring.
- 5 Remove the rotor seal (and isolation seal if damaged or contaminated).
- 6 Install the new isolation seal (if required). Ensure the metal spring inside the ring faces towards the valve body.
- 7 Install the new rotor seal.
- 8 Replace the stator ring. Ensure the stator ring is flush with the valve body.
- 9 Place the new (if required) stator face in place on the stator head. Reinstall the stator head.
- 10 Insert the stator screws in the stator head. Tighten the screws alternately two turns at a time until the stator head is secure.
- 11 Reconnect the pump capillaries to the valve ports.

CAUTION

Wrong use of Pressure Test may damage valve.

The current implementation of the Pressure Test automatically uses the maximum pressure generated by the pump used in the system.

- ✓ **Do not use the test for modules having a lower maximum pressure than the pump as this will damage the valve. For example do not use 400 bar valve in a TCC or Flexible Cube in combination with a 600 bar pump.**

- 12 Perform a **Pressure Test** to ensure the valve is pressure tight.




Figure 106 Valve Head Parts (example)

1	Stator screws
2	Stator head assembly
3	Stator ring screws (not available)
4	Stator ring (available for service only)
5	Rotor seal
6	Bearing ring
7	Spanner nut (available for service only)

NOTE

Figure 106 on page 273 illustrates replacement parts for the valve heads, with the 6-column Selector valve as an example. The valves can vary in their appearance and do not necessarily include all of the illustrated parts. Neither, every spare part is available for each flavor of the valve.

Replacing the Fuses of the Infinity Valve Drive

When If the flow module shows no reaction.

Tools required **Description**
Screwdriver


Parts required	#	p/n	Description
	2	2110-1486	Fuse 2 AT250 V

WARNING


Electrical shock

- ✓ Disconnect the module from line power before changing a fuse or trying to open the hatch of the power input socket.
- ✓ Never reconnect the line power before having the power input socket closed.

1 Switch off the instrument. The line switch is located at the rear of the module.




2 Disconnect the power cable from the power input socket.




Replacing the Fuses of the Infinity Valve Drive

- 3** To access the fuse drawer, gently lift the outer plastic housing of the power inlet socket using a flat screwdriver.




- 4** Pull out the fuse drawer as shown.




- 5** Replace the defect fuse(s).

- 6** Slide in the fuse drawer and push till it fits tightly.



- 7** Finally, close the fuse drawer housing, reconnect the instrument to the power line and switch it on.



Replace the Module Firmware

When

The installation of newer firmware might be necessary

- if a newer version solves problems of older versions or
- to keep all systems on the same (validated) revision.

The installation of older firmware might be necessary

- to keep all systems on the same (validated) revision or
- if a new module with newer firmware is added to a system or
- if third party control software requires a special version.

Tools required**Description**

Agilent Lab Advisor software

Parts required**# Description**

- | | |
|---|---|
| 1 | Firmware, tools and documentation from Agilent web site |
|---|---|

Preparations

Read update documentation provided with the Firmware Update Tool.

To upgrade/downgrade the module's firmware carry out the following steps:

- 1 Download the required module firmware, the latest FW Update Tool and the documentation from the Agilent web.
<http://www.agilent.com/en-us/firmwareDownload?whid=69761>
- 2 For loading the firmware into the module follow the instructions in the documentation.

Module Specific Information

There is no specific information for this module.

2D-LC Loops	278
2D-LC Capillaries	279
ASM Capillaries	280
Pressure Release Kit	281
2D-LC Easy Starter Kit	282
Valve Drive Parts	283
Valve Driver Parts Infinity II	284
Valve Head Parts	285
Technical specifications	286
Tools	286
Vale Options Overview (for 2D-LC)	287
Obsolete Valve Heads	290
MS Diverter Valve	291
Additional Parts for the MS Diverter Valve Setup	292
Valve Options Overview (for G7116B)	293
Additional Heater Devices	295
Additional Heater Devices (for G7116B)	296
Accessories and Consumables (for G7116B)	297
InfinityLab Quick Connect and Quick Turn Fittings	299
InfinityLab Quick Connect Fittings	300
InfinityLab Quick Connect Fitting Replacement Capillaries	301
InfinityLab Quick Turn Fitting	302
Capillary Kits	303
Important Customer Web Links	304

This chapter provides information on parts material required for the solution.

2D-LC Loops

2D-LC Loops for Standard 2D-LC

p/n	Description
5067-5440	Calibrated loop kit for 2D-LC
5067-5446	Loop housing kit
5067-5424	20 µL Loop 2D-LC
5067-5425	40 µL Loop 2D-LC
5067-5437	60 µL Loop 2D-LC
5067-5426	80 µL Loop 2D-LC

2D-LC Loops for MHC valve Fitting M4

p/n	Description
5067-6643	Capillary ST 0.35x104 mm, M/M, 10
5067-6644	Capillary ST 0.35x208 mm, M/M, 20 µL
5067-5926	Capillary ST 0.35x 420 mm M/M 40 µl
5067-6645	Capillary ST 0.35x831 mm, M/M, 80 µL
5067-6646	Capillary ST 0.35x1247 mm, M/M, 120 µL
5067-6647	Capillary ST 0.35x1870 mm, M/M, 180 µL
5067-6141	M4 Blank nut
5023-2504	Hex driver SW-4 slitted

2D-LC Capillaries

1200 Infinity Series 2D-LC Capillary Kit

p/n	Description
5021-1820	Flex capillary, 0.12 mm x 105 mm, no fittings
G1316-87321	Capillary column-heat exchanger 105 mm lg, 0.17 mm i.d.
5021-1822	Capillary, 0.12 mm x 280 mm
5021-1823	Capillary column – detector SST 400 mm lg, 0.12 mm i.d.
5021-1819	Capillary ST 0.17 mm x 400 mm S/S
5065-9964	Capillary ST 0.12 mm x 500 mm
5067-4609	Capillary ST 0.17 mm x 500 mm SX/-
5067-4669	Capillary ST 0.12 mm x 600 mm S/SL
01078-87305	Capillary, 0.17 mm x 80 cm, male fit
5065-4454	Fitting screw long, front and back ferrules 10/pk
G1316-60005	Low Dispersion Heat Exchanger Double Assy
G7116-60015	Quick Connect Heat Exchanger Standard
5500-1188	Quick Turn Capillary ST 0.12 mm x 105 mm, long socket

InfinityLab 2D-LC Capillary Kit legacy

p/n	Description
5067-4651	Capillary ST 0.12 mm x 280 mm SL/SX
5067-4669	Capillary ST 0.12 mm x 600 mm S/SL
5500-1245	Capillary ST 0.17 mm x 400 mm SI/SI
5500-1251	Capillary ST 0.12 mmX 400 mm SL/SL
5500-1240	Capillary ST 0.17 mm x 105 mm SL/SL
5500-1227	Capillary ST 0.17 mm x 150 mm SL-SL
5500-1217	Capillary, ST, 0.17 mm x 900 mm SI/SX
5067-4608	Capillary ST 0.17 mm x 280 mm SX/S
5067-4670	Capillary ST 0.17 mm ID 600 mm pre-swaged

ASM Capillaries

ASM Valve Capillary Replacement Kit

p/n	Description
5500-1300	Capillary ST 0.12x85M/M ASM
5500-1301	Capillary ST 0.12x170M/M ASM
5500-1302	Capillary ST 0.12x340M/M ASM
5500-1303	Capillary ST 0.12x680M/M ASM
5500-1376	Capillary ST 0.12x170M/M transfer

Pressure Release Kit




Figure 107 Pressure release kit, parts

Item	p/n	Description
	G4236-60010	2D-LC Pressure Release Kit
	0100-0969	TEE, ST, 1/16 inch, Low Dead Volume Not shown
1	5021-1816	Capillary i.d. 0.17 mm, 105 mm lg
2	5022-2184	Union, stand LC flow, no fitting
3	G7167-87307	500 µL Loop extension
4	G4212-60022	Pressure Relief Valve
5	5067-5939	Splitter-Capillary 0.05-ID L-1000 mm

2D-LC Easy Starter Kit

2D-LC Easy Starter Kit

p/n	Description
5190-6895	2D-LC starter sample, 1 x 2 mL
G2453-85060	Formic Acid-Reagent Grade 5 mL (5 cc)
858700-902	RRHD SB-C18, 2.1x100 mm, 1.8 µm, 1200 bar
857768-901	RRHD Bonus-RP, 2.1x50 mm, 1.8 µm, 1200 bar
959757-302	RRHD Eclipse Plus C18, 3.0x50 mm, 1.8 µm

Valve Drive Parts

Item	p/n	Description
1	5043-0275	Clamp guide For attaching the valve to a rail assembly
2	5067-4792	Leak sensor assembly External leak sensor
3	5043-0271	Holder leak plane
4	5043-0270	Leak plane
5	5068-0106	Spanner nut
	2110-1486	Fuse 2 AT250 V
6	5067-4634	Valve rail assembly
7	5067-1510	Rail assy for column organizer




Figure 108 Valve Drive Parts

Valve Driver Parts Infinity II

Item	p/n	Description
1	5067-6138	Valve Holder Kit Right-IF-II-G For G7116A/B
	5067-6139	Valve Holder Kit Left-IF-II-G For G7116A/B (Not shown)
2	5067-5685	Clamp Guide Kit-IF-II
3	5067-4792	Leak sensor assembly External leak plane
4	5043-0271	Holder leak plane
5	5043-0270	Leak plane
6	2110-1486	Fuse 2 AT250 V
7	5063-6527	Tubing, Silicon Rubber, 1.2 m, ID/OD 6/9 mm
8	5181-1519	CAN cable, Agilent module to module, 1 m
9	5500-1156	T-Tube Connector ID6.4
10	5043-0269	Adapter-profile For G1170A (Multiple valve drives can be connected with adapter profiles)




Figure 109 Parts for 1260/1290 Infinity II Valve Drive

Valve Head Parts

NOTE

The figure below illustrates replacement parts for the valve heads, with the 12ps/13pt Selector valve as an example. The valves can vary in their appearance and do not necessarily include all of the illustrated parts. Neither, every spare part is available for each flavor of the valve.




Figure 110 Valve Head Parts (example)

1	Stator screws
2	Stator head assembly
3	Stator ring screws (not available)
4	Stator ring (available for service only)
5	Rotor seal
6	Bearing ring
7	Spanner nut (available for service only)

Technical specifications

Table 18 Technical specifications

Max. Pressure:	1300 bar
Liquid Contacts:	Stainless Steel, PEEK
Connections:	Accepts 10-32 male threaded and M4 fittings

Tools

Tool for extra fittings

p/n	Description
8710-2462	Hex Key Driver 3/32 inch
5023-2504	Hex driver SW-4 slitted For M4 fittings
5067-6141	M4 Blank nut For plugging unused valve ports
5067-6127	Blank Nut SL

Vale Options Overview (for 2D-LC)

The 1300 bar InfnitLab Quick Change Valves are backward compatible to the 1200 bar Valves.

NOTE

The service life of a stator depends on the stress to which the 2D-LC valve is subjected. Therefore, a visual inspection of the surface during maintenance is very important. If scratches or heavy wear is visible during the inspection or the valve has passed 30000 switches, the stator must be replaced.

G4136A 2D-LC Valve Kit, Standard

#	p/n	Description
1	5067-4244	2D-LC Valve Head, 1300 bar
1	5067-5440	Calibrated loop kit for 2D-LC
1	5067-6171	Capillary Kit 2D-LC, Infinity Classic (optional) Internal part, not orderable
1	G4236-68000	2D-LC Easy Starter Kit Internal part, not orderable
2	G4242-64000	Multiple Heart-Cutting Valve
1	5067-6585	Capillary Kit 2D-LC, 1290 Infinity II Internal part, not orderable
1	G1680-63721	Network LAN Switch

5067-4244 2D-LC Valve Head, 1300 bar

#	p/n	Description
3	1535-4857	Stator screws, 10/pk
1	1535-4045	Bearing ring
1	5068-0214	Rotor Seal (VHP)
1	5068-0120	Stator ring
1	5068-0115	Stator

G4243A 2D-LC Valve Kit, ASM

#	p/n	Description
1	5067-4266	2D-LC ASM Valve Head, 1300 bar
1	G4236-68000	2D-LC Easy Starter Kit Internal part, not orderable
1	G1680-63721	Network LAN Switch
1	5500-1300	Capillary ST 0.12x85M/M ASM
1	5500-1301	Capillary ST 0.12x170M/M ASM
1	5500-1302	Capillary ST 0.12x340M/M ASM
1	5500-1303	Capillary ST 0.12x680M/M ASM
4	5500-1376	Capillary ST 0.12x170M/M transfer

5067-4266 2D-LC ASM Valve Head, 1300 bar

p/n	Description
5068-0019	Stator screws
5068-0257	Bearing Ring
5068-0240	Rotor Seal (VHP)
5068-0239	Stator

G4242-64000 Multiple Heart-Cutting Valve

#	p/n	Description
1	5067-4273	6-column selector valve head, 1300 bar
6	5067-5926	Capillary ST 0.35x 420 mm M/M 40 µl
2	5500-1270	Capillary ST 0.12 mm x 170 mm S/M
1	5043-0269	Adapter-profile

5067-4273 6-column selector valve head, 1300 bar

#	p/n	Description
5	5068-0089	Stator screws
1	1535-4045	Bearing ring
1	5068-0242	Rotor Seal (PEEK)
1	5068-0241	Stator Head

Obsolete Valve Heads

The following 1200 bar valve heads are no longer orderable:

5067-4214 2D-LC Valve 1200 bar legacy

p/n	Description
5068-0186	Rotor Seal (Vespel)
5068-0115	Stator
1535-4857	Stator screws, 10/pk
1535-4045	Bearing ring

Multiple Heart-Cutting Valve legacy

#	p/n	Description
1	5067-4142	Valve head 6 column selector (1200 bar)
6	5067-5926	Capillary ST 0.35x 420 mm M/M 40 µl
1	5974-0197	Capillary Cover Label
2	5067-5113	Capillary ST 0.17 mm x 250 mm SL/M
2	5067-6188	Capillary ST 0.17 mm x 500 mm SL-M

5067-4142 6-Column selector valve 1200 bar legacy

p/n	Description
5068-0077	Stator Head
5068-0067	Rotor Seal (Vespel)
5068-0089	Stator screws
1535-4045	Bearing ring

MS Diverter Valve

G4231A 2 position/6 port valve head, 800 bar

p/n	Description
5067-4282	2 position/6 port valve head, 800 bar
5067-4730	2/10 Cap kit 0.17 mm
5067-4249	2/6 Cap Kit 0.12 mm, incl. QC-HEX
5067-4250	2/6 Cap Kit 0.12 mm, incl. LD-HEX
5067-6597	2/6 Cap Kit 0.17 mm, incl. QC-HEX

Alternative diverter valves (2 position / 6 port, PEEK Rotor Seal)

p/n	Description
5067-4137	Valve Head 2 Postion / 6 Port, 600 bar
5067-4282	2 position/6 port valve head, 800 bar
0101-1409	Rotor Seal (PEEK)

Alternative diverter valves (2 position / 10 port, PEEK Rotor Seal)

p/n	Description
5067-4145	2ps/10pt valve head, 600 bar
5067-4283	2ps/10pt valve head, 800 bar
0101-1415	Rotor Seal (PEEK)

Alternative diverter valves (2 position / 6 port, Vespel Rotor Seal)

p/n	Description
5067-4117	2ps/6pt ultra high pressure valve head, 1200 bar
5068-0008	Rotor Seal (Vespel)

Alternative diverter valves (2 position / 10 port, Vespel Rotor Seal)

p/n	Description
5067-4118	2ps/10pt ultra high pressure valve head, 1200 bar
5068-0012	Rotor Seal (Vespel)

Additional Parts for the MS Diverter Valve Setup

p/n	Description
G4212-60022	Pressure Relief Valve
5067-4606	Capillary ST 0.12 mm x 400 mm SX/-
0890-1915	Capillary PK 0.13 mm x 150 cm
5500-1228	Capillary ST 0.3 mm x 80 mm SL-SL
5063-6591	PEEK Fittings 10/PK
0100-0969	TEE, ST, 1/16 inch, Low Dead Volume
5067-6127	Blank Nut SL
5062-2462	Tube PTFE 0.7 mm x 5 m, 1.6 mm od

Valve Options Overview (for G7116B)

Valve Options Overview

NOTE

The 1300 bar InfinitLab Quick Change Valves are backward compatible to the 1200 bar Valves.

Table 19 Replacement parts standard valve heads for G7116B

Valve Head	Rotor Seal	Stator Head	Stator Screws	Stator Ring
5067-4233 8 Column Selector 1300 bar	5068-0200 (PEEK)	5068-0199	5068-0089	n.a.
5067-4241 2ps/6pt Valve 1300 bar	5068-0207 (PEEK)	5068-0006	1535-4857	n.a.
5067-4240 2ps/10pt Valve 1300 bar	5068-0205 (PEEK)	5068-0011	5068-0019	n.a.
5067-4273 6 Column Selector 1300 bar	5068-0242 (PEEK)	5068-0241	5068-0089	n.a.
5067-4284 6 Column Selector 800 bar	5068-0298	5068-0241	5068-0089	n.a.

Obsolete Valve Heads

The following 1200 bar valve heads are no longer orderable:

Table 20 Replacement parts obsolete valve heads for G7116B

Valve Head	Rotor Seal	Stator Head	Stator Screws	Stator Ring
5067-4121 8ps/9pt, 1200 bar	5068-0002 (Vespel)	5068-0001	1535-4857	5068-0127
5067-4117 2ps/6pt, 1200 bar	5068-0008 (Vespel)	5068-0006	1535-4857	5068-0127
5067-4118 2ps/10pt, 1200 bar	5068-0012 (Vespel)	5068-0011	5068-0019	n.a
5067-4142 6 Column Selector, 1200 bar	5068-0067 (Vespel)	5067-0077	5068-0089	n.a

Additional Heater Devices

Table 21 Heat Exchanger Overview

Flow rate	0.075 mm i.d. capillary	0.12 mm i.d. capillary
< 2 mL/min	<i>Ultra-low Dispersion</i>	<i>Standard Flow</i>
	G7116-60021 (Internal volume: 1.0 µL)	G7116-60015 (Internal volume: 1.6 µL)
> 2 mL/min		<i>High Flow</i>
		G7116-60031 (Internal volume: 3.0 µL)

For details, see [Table 22](#) on page 296.

Additional Heater Devices (for G7116B)

Blank heater assemblies without capillaries and fittings:

Table 22 InfinityLab Quick Connect Heat Exchanger

Item	Description
A long, thin, cylindrical metal heat exchanger with two threaded ports at each end. It is shown with two grey quick-connect fittings attached to the ports. The device has a label with the InfinityLab logo and part number G7116-60015.	Quick Connect Heat Exchanger Standard (G7116-60015)
A long, thin, cylindrical metal heat exchanger with two threaded ports at each end. It is shown with two grey quick-connect fittings attached to the ports. The device has a label with the InfinityLab logo and part number G7116-60021.	Quick Connect Heat Exchanger Ultra Low Dispersion (G7116-60021) NOTE Use InfinityLab Quick Turn Fittings to connect to the Quick Connect Heat Exchanger Ultra Low Dispersion.
A long, thin, cylindrical metal heat exchanger with two threaded ports at each end. It is shown with two grey quick-connect fittings attached to the ports. The device has a label with the InfinityLab logo and part number G7116-60031.	Quick Connect Heat Exchanger High Flow (G7116-60031)

Accessories and Consumables (for G7116B)

Accessory Kit (for G7116B)

The Accessory Kit (for G7116B) contains accessories and tools needed for the installation and maintenance.

p/n	Description
5181-1516	CAN cable, Agilent module to module, 0.5 m
5063-6527	Tubing, Silicon Rubber, 1.2 m, ID/OD 6/9 mm
5500-1191	Capillary ST 0.12 mm x 280 mm, long socket
5067-5966	InfinityLab Quick Turn Fitting
5067-5957	InfinityLab Quick Connect Assy ST 0.12 mm x 105 mm
G7116-60015	Quick Connect Heat Exchanger Standard
G7116-68003	Column Holder Lamella, 2/pk
5043-0915	Fitting mounting tool
G7116-60006	Divider Assembly MCT
5022-2184	Union, stand LC flow, no fitting
	Double Drain Connector

Available Consumables (for G7116B)

p/n	Description
G7116-68003	Column Holder Lamella, 2/pk
G7116-68004	Column Holder Clamp (2/pk)
5500-1191	Capillary ST 0.12 mm x 280 mm, long socket (Capillary from column outlet to detector, no fittings.)
G7116-60006	Divider Assembly MCT For separating different temperature zones between left and right heater elements.
5067-5917	InfinityLab Column Identification Tag Blank column ID tag (column ID tag reader kit is required)

Number Kit

p/n	Description
5067-6654	Number Kit 1-8 colored Column Info in red, blue, green, cyan, yellow, black, white, and gray

InfinityLab Quick Connect and Quick Turn Fittings

For further info check either the consumables catalog or [“Important Customer Web Links”](#) on page 304.

InfinityLab Quick Connect Fittings




Figure 111 InfinityLab Quick Connect Fitting

p/n	Description
5067-5965	InfinityLab Quick Connect LC fitting (fitting without pre-installed capillary)
5043-0924	Front Ferrule for Quick Connect/Turn Fitting
5067-5961	InfinityLab Quick Connect Assy ST 0.075 mm x 105 mm
5067-6163	InfinityLab Quick Connect Assy ST 0.075 mm x 150 mm
5067-6164	InfinityLab Quick Connect Assy ST 0.075 mm x 220 mm
5067-6165	InfinityLab Quick Connect Assy ST 0.075 mm x 280 mm
5067-5957	InfinityLab Quick Connect Assy ST 0.12 mm x 105 mm
5067-5958	InfinityLab Quick Connect Assy ST 0.12 mm x 150 mm
5067-5959	InfinityLab Quick Connect Assy ST 0.12 mm x 220 mm
5067-5960	InfinityLab Quick Connect Assy ST 0.12 mm x 280 mm
5067-6166	InfinityLab Quick Connect Assy ST 0.17 mm x 105 mm
5067-6167	InfinityLab Quick Connect Assy ST 0.17 mm x 150 mm
5067-6168	InfinityLab Quick Connect Assy ST 0.17 mm x 220 mm
5067-6169	InfinityLab Quick Connect Assy ST 0.17 mm x 280 mm

InfinityLab Quick Connect Fitting Replacement Capillaries

p/n	Description
5500-1174	InfinityLab Capillary ST 0.075 mm x 105 mm
5500-1175	InfinityLab Capillary ST 0.075 mm x 150 mm
5500-1176	InfinityLab Capillary ST 0.075 mm x 220 mm
5500-1177	InfinityLab Capillary ST 0.075 mm x 250 mm
5500-1178	InfinityLab Capillary ST 0.075 mm x 280 mm
5500-1173	InfinityLab Capillary ST 0.12 mm x 105 mm
5500-1172	InfinityLab Capillary ST 0.12 mm x 150 mm
5500-1171	InfinityLab Capillary ST 0.12 mm x 220 mm
5500-1170	InfinityLab Capillary ST 0.12 mm x 280 mm
5500-1179	InfinityLab Capillary ST 0.12 mm x 400 mm
5500-1180	InfinityLab Capillary ST 0.12 mm x 500 mm
5500-1181	InfinityLab Capillary ST 0.17 mm x 105 mm
5500-1182	InfinityLab Capillary ST 0.17 mm x 150 mm
5500-1183	InfinityLab Capillary ST 0.17 mm x 220 mm
5500-1230	InfinityLab Capillary ST 0.17 mm x 280 mm
5500-1231	InfinityLab Capillary ST 0.17 mm x 500 mm
5500-1259	InfinityLab Capillary ST 0.25 mm x 150 mm
5500-1260	InfinityLab Capillary ST 0.25 mm x 400 mm

InfinityLab Quick Turn Fitting




Figure 112 InfinityLab Quick Turn Fitting

p/n	Description
5067-5966	InfinityLab Quick Turn Fitting
5043-0924	Front Ferrule for Quick Connect/Turn Fitting

Capillary Kits

NOTE

Further capillary kits can be found in the *Agilent 1290 Infinity Valve Drive and Valve Heads User Manual* or on the webpage.

Table 23 Common capillary kit

Part Number	Connection	Amount
Capillary ST 0.12 mm x 340 mm S/SX (5067-4647)	Autosampler to valve	1
Capillary ST 0.17 mm x 700 mm S/SX (5067-4648)	2 nd pump to valve	1
Capillary ST 0.12 mm x 90 mm S/SX (5067-4649)	Valve to heat exchanger	2
Capillary ST 0.12 mm x 150 mm SL/SX (5067-4650)	Short column to valve	2
Capillary ST 0.12 mm x 280 mm SL/SX (5067-4651)	Long column to valve	2
Capillary ST 0.12 mm x 120 mm SX/SX (5067-4652)	Valve to valve	1
Capillary ST 0.12 mm x 200 mm S/SX (5067-4653)	Valve to detector	1
Waste tubing, 2 m (0890-1713)	Valve to waste	2 m
Waste tube (G1375-87326) (includes M4 PEEK fitting)	Valve to waste	1
Plastic fitting (0100-1259)		4
Bag - plastics (9222-0518)		1

Important Customer Web Links

- Videos about specific preparation requirements for your instrument can be found by searching the *Agilent YouTube* channel at <https://www.youtube.com/user/agilent>
- To access *Agilent University*, visit <http://www.agilent.com/crosslab/university/> to learn about training options, which include online, classroom and onsite delivery.
A training specialist can work directly with you to help determine your best options.
- A useful *Agilent Resource Center* web page is available, which includes short videos on maintenance, quick lists of consumables for new instruments, and other valuable information.
- Need technical support, FAQs, supplies? – visit our *Support Home page* at <http://www.agilent.com/search/support>
- Get answers. Share insights. Build connections:
Join the *Agilent Community* at <https://community.agilent.com/welcome>



13

Alternative ways to install the System

Legacy Stack Configuration 306

1D/2D Switching 309

This chapter describes alternative ways to install and setup the system.

Legacy Stack Configuration

The following configurations optimize the system flow path, ensuring minimum delay volume in an Agilent Infinity I configuration.

NOTE

The capillary connections should be as short as possible, to ensure optimum performance of the system.




Figure 113 Legacy stack configuration for Multiple Heart-cutting 2D-LC




Figure 114 Previous stack configuration for Multiple Heart-cutting 2D-LC based on the column organized/valve holder

Table 24 1290 Infinity Binary LC in first dimension

	Partnumber	Description	Comment
1st Dim	G4220A	1290 Infinity Binary Pump	
	G4226A	1290 Infinity Autosampler	
	G1330B	1290 Thermostat	
	G1316C	1290 Thermostatted Column Compartment	
	G4212A	1290 Infinity Diode-Array Detector	
2nd Dim	G2198AA	2D-LC Acquisition Software	
	G4220A	1290 Infinity Binary Pump	
	G1170A	1290 Infinity Valve Drive	For 2D-LC valve
	G4236A	2D-LC Valve Kit, 1200 bar	
	G4212A	1290 Infinity Diode-Array Detector	

NOTE

There are a few other configurations when using multiple detectors at different positions. Examples are:

- After first dimension column, or
- At the waste line in addition to the standard 2D-LC detector after the second dimension columns.

If quantitative information is required, it is recommended to use same detector types and flow-cells.

NOTE

¹D flow cells require a minimum pressure stability of 60 bar (~870 psi) (which excludes FLD and RID detectors).

NOTE

²D detectors should be connected to CAN bus for automatic run time extension. To be considered for ELSD or ³rd party detectors.

NOTE

Avoid high data rates in the first dimension for peak-based sampling. Otherwise, RI effects may be recognized as peaks, which can cause unwanted sampling. 5 Hz give enough data points for correct peak detection of long ¹D runs.

1D/2D Switching

Usually, it is necessary to reconnect capillaries to switch between performing 1D-LC methods and comprehensive 2D-LC methods using the same LC system. The 1D/2D switching setup is a flexible Agilent 1290 Infinity 2D-LC Solution with one detector. It contains an additional 2-position/6 port valve to allow switching between 1D-LC and comprehensive 2D-LC analysis with a few mouse clicks, and without any hardware change. For details, see Automated Switching Between 1D-LC and Comprehensive 2D-LC Analysis (5991-4843EN).




Figure 115 Setup of the Agilent 1290 Infinity 2D-LC Solution with an additional 2-position/6-port valve positioned for 1D-LC analysis.




Figure 116 Setup of the Agilent 1290 Infinity 2D-LC Solution with an additional 2-position/6-port valve positioned for comprehensive 2D-LC analysis; (A) Collection of the effluent from the first dimension column in loop 1 (between ports 1 and 8); (B) Collection of the effluent from the first dimension column in loop 2 (between ports 4 and 5), analysis of the content of loop 1 on the second dimension column.

Agilent 2D-LC Solution with the selectivity of mass selective detection (MSD)

Agilent 2D-LC Solution with an Agilent Single Quadrupole Detector 312

Using the Single Quadrupole Detector with 2D-LC 313

Agilent 2D-LC Solution with high-end MS detection 318

User Interface/Features 320

This chapter describes the different options to use the Agilent 2D-LC Solution with mass selective detection (MSD).

Today the 2D-LC software is an OpenLab CDS ChemStation Edition plug in, which comprises three functions: system control, data acquisition, and data analysis/reporting. The current 2D-LC software requires configuring at least one UV detector in either or both dimensions. A single quadrupole detector may instead be configured for the second dimension. Only UV detectors may be configured in 1D as peak triggering detectors. MS detectors can be used in the second dimension for acquiring 2D-LC data either in OpenLab CDS ChemStation Edition (single quadrupole) or MassHunter workstation (QQQ/Q-TOF).

Agilent 2D-LC Solution with an Agilent Single Quadrupole Detector


Agilent OpenLab 2D-LC Software is a fully integrated solution from method setup to data analysis for 2D-LC/MS within OpenLab CDS ChemStation Edition. It combines the separation power of two-dimensional chromatography with the selectivity of mass selective detection (MSD) in the second dimension. It can be used with all Agilent Single Quadrupole Detectors supported by ChemStation including latest LC/MSD (G6125B) and LC/MSD XT (G6135B). Only UV detectors may be configured in ¹D as peak triggering detectors.

NOTE

The 2D-LC system has been optimized for detectors connected to the CAN bus interface. For non-CAN detectors (e.g. third party detectors), suitable run time settings need to be considered in order to get complete signal data, see KPR 25 in the software status bulletin for details.


Using the Single Quadrupole Detector with 2D-LC

As a destructive detector, the single quadrupole detector is usually installed as a ²D detector. It can be configured as described in section "Detector configuration".



When setting up the MSD (**Instrument> Set up MSD signals> ...**), set **noLimit** as the **Stop Time**. This is required as the run time needs to be extended for unparking cuts.

2D-LC Viewer for Single Quadrupole Data





Select signals

☒ ☒
☐ DAD2 A, Sig=254,4 Ref=360,100
☐ DAD2 B, Sig=245,4 Ref=360,100
☐ DAD2 C, Sig=230,4 Ref=360,100
☐ DAD2 D, Sig=214,4 Ref=360,100
☒ MSD1 TIC, MS File
☐ MSD1 190, EIC=189.7:190.7
☐ MSD1 263, EIC=239.7:286.0
☐ MSD1 261, EIC=260.7:261.7
☒ MSD1 262, EIC=254.8:270.1
 Ok

Multiple signals can be managed using the signal selector. Use checkmarks for choosing signals to be displayed.



Select highlighted/all items

As a shortcut, multiple signals can quickly be highlighted by control- or shift-clicking signal names for the ²D chromatogram panel and clicking the selection buttons.



time/apex/average per range

Spectra can be displayed for a time or a peak apex by clicking the ²D chromatogram or dragging the mouse for an average spectrum per time range.




Pin spectra

Spectra can be pinned for storing them per associated cut for display and reporting:


MSD1 SPC,time=2.788 (Cut# 4 of 001-P1-B1-2DLC F25.D)
 MSD1 SPC,time=2.913 (Cut# 4 of 001-P1-B1-2DLC F25.D)
 Ok




Extract ion chromatogram



Ion chromatograms can be extracted for a mass or mass range selected in the spectrum or entered for the 2D chromatogram(s).



Ion chromatograms can also be extracted from spectra by clicking the EIC button and then selecting a mass.



Separate versus overlaid mode

Multiple (extracted ion) chromatograms or spectra can be displayed in the same window (overlay) or in separate windows.



same scale/full scale

Diagrams can be displayed using the same scale or the full scale for each diagram.



linked x and y axes

For zooming, the axes can optionally be linked for all displays. In case of linked axes, all displays will be zoomed to the same range.

Diverter Valve Configuration

Diverter Valve

Diverter Valve (G1170A:DE11998844) 6Port...



Waste: Port 1 -> 6 MSD: Port 1 -> 2

A diverter valve can be configured e.g. for removing buffers before they reach an MS detector. Any 2-position LC valve can be chosen in the configuration. Valves should be installed to external valve drives G1170A. Control is done via the CAN connection between 2D-LC modules and the diverter valve, which is not yet available for diverter valves in MS detectors. The user interface shows, which valve outlet ports go to waste vs. detector.

Diverter Valve Method Parameters


Diverter valve

☒ Switch valve to MSD after

1.00


min (²D time)

A two-position LC valve (e.g. a 2-pos/6-port valve) can be configured in the 2D-LC Configuration for diverting ²D eluent with buffers before they enter an MS detector. As a method parameter, a ²D time can be set, after which the flow path is directed to the detector.



Detectors			
Module name	Usage	Peak trigger	Transfer volume [ul]
DAD (G4212A PP00055025)	¹ D Detector	☐	13.00
DAD (G4212A DEBAF01751)	² D Detector		
G6150B MSD (G6150B)	² D Detector		
Transfer volume of reference signal 15.00 ul			

A previously acquired ¹D MSD signal can be used as a reference signal after configuring an MS detector in ²D without configuring a ¹D peak trigger. In this case, the transfer volume of the reference signal can be entered for that previously installed ¹D detector:



- Reference signal is acquired with MSD in first dimension.
- Capillary connections are changed and MSD is moved to ²D. Then, there is no peakdetector in ¹D. The transfer volume V_t needs to be calculated as difference $V_{t,2}-V_{t,1}$ and configured.
- Reference signal is used for defining cuts in sampling table.
- MSD is used for detection in second dimension.

Agilent 2D-LC Solution with high-end MS detection




Figure 117 Typical system setup

NOTE

OpenLab CDS ChemStation and MassHunter have passed the test successfully for co-execution regarding parallel data acquisition. For further information please check the recommended conditions in the technical note *2D Chromatogram Creator for MassHunter*.

The LC part is an LC system as described in the previous pages of this 2D-LC system manual. A typical system uses a UV detector in the first dimension and an optional UV detector in the second dimension. The 2D-LC instrument control and UV signal acquisition is still done with OpenLab CDS ChemStation (CS) Edition. Therefore the 2D-LC software requires configuring at least one UV detector.

The MS detector, a TOF, Q-TOF, or Triple Quadrupole, is used in the second dimension. It is controlled by MassHunter (MH) Acquisition software, which acquires the MS signal.

For synchronizing signals of LC and MS, an APG remote cable connects both parts of the system, which communicates start and stop signals.

NOTE

For synchronizing run times (that is starting and stopping measurements simultaneously), please use **External start** as a run parameter in the Masshunter Acquisition software. You will find the information in section **Additional information**, add parameter **run type** and choose option **external run**. Set stop time **no limit** in MassHunter.

NOTE

Only UV detectors can be configured in ¹D as peak triggering detectors.

For comprehensive 2D-LC MS measurements, we recommend using the GC Image LCxLC-HRMS Edition software. This software can display and analyze the 2D-LC MS data. For further information please check the following:

- Chapter “Data Analysis for Comprehensive 2D-LC (LCxLC)” on page 219,
- LC Video Tutorial on the GC Image DVD,
- Or in the help file of the GC Image SW.

Table 25 Software for using Agilent 2D-LC Solution with high-end MS detection

Type of 2D-LC Mode	Instrument Control	Data Acquisition	Data Analysis
Multiple Heart-Cutting ¹ , High Resolution Sampling	ChemStation plugin	ChemStation plugin MassHunter	LC Image ¹ or MassHunter in combination with 2D Chromatogram Creator for MassHunter software
Comprehensive	ChemStation plugin	ChemStation plugin MassHunter	LC Image

¹ MHC data analysis can be displayed as a false result because peaks from cuts that not belong together can be integrated as one bulb.

But the story is different if you are using Multiple Heart-Cutting (MHC) or High-Resolution to get the correct cut and identify the mass correctly with high-end MS detection. MassHunter has no knowledge about the size, order, and position of cuts created by OpenLab CDS ChemStation Edition. The raw MS TIC signal continuously measures MS data of the 2D eluent. Therefore, you need a tool that aligns the different cut operations from the OpenLab SW with the MS data of the eluent. For this purpose, you must use the 2D Chromatogram Creator for MassHunter.

User Interface/Features




Figure 118 Use single data or sequences/worklists and set delay between optional 2D UV detector and MS detector for aligning signals

²D Chromatogram Creator for MassHunter provides an easy access to 2D-LC measurements for Agilent MassHunter users. The ²D Chromatogram Creator for MassHunter is used for combining UV signals and metadata of cuts in OpenLab CDS ChemStation Edition with the MS signal acquired by MassHunter.

As output, it creates following MassHunter chromatograms:

- A series of ²D chromatograms, one per cut including MS and optional UV signals
- If applicable, a ¹D UV chromatogram (for a typical but optional ¹D UV detector) in MassHunter data format including cut data
- An optional raw ²D UV signal, and a
- ²D MS raw signal

You will find the 2D Chromatogram Creator for MassHunter on the 2D-LC add-on DVD.

NOTE

More information about Installation, Example Data, and Workflow is available in the Agilent technical note *²D Chromatogram Creator for MassHunter* (G2198-90101).

NOTE

Further information about the 2D-LC Diverter Valve Solution can be found in the technical note *Agilent InfinityLab 2D-LC Solution with mass spectrometric detection and diverter valve* (G4236-90100).

15

Theoretical Background

Theoretical basis of 2D-LC	323
Orthogonality	325
Resolution	325
Peak Capacity	328
2nd dimension as detector	331
Successful Mode Combinations	333
Solvent Elution Modes	334
Practical Issues	340

This chapter gives the theoretical background of 2D-LC and describes the system components (soft- and hardware) of the Agilent 1290 Infinity II 2D-LC Solution.

Theoretical basis of 2D-LC

In 2D-LC, fractions from a chromatographic system (1st dimension) are transferred to a second chromatographic separation system (2nd dimension). So 2D-LC bases on the application of two independent liquid phase separation systems to a sample. 2D-LC is mainly used to improve resolution and sensitivity or to decrease analysis time.




Figure 119 Peak capacity relationship between peak capacities of orthogonal first and second dimension

The most important benefit of 2D-LC over 1D-LC is the increase of resolving power, which is especially important if dealing with complex samples.

For an overview on the main differences between 1D- and 2D-LC, refer to the following topics:

- “Orthogonality” on page 325,
- “Resolution” on page 325, and
- “Peak Capacity” on page 328

The following different methods of 2D-LC exist:

- Heartcutting (LC-LC)
Only interesting portion of the first dimension effluent transferred to the second dimension.
- Comprehensive (LCxLC)
Entirety of first dimension effluent sequentially transferred to the second dimension.

Orthogonality

The 2D-LC separation power depends the fact that the two selectivity mechanisms in the different separation stages must be as different as possible. If the mechanisms are completely different and independent the two separations are called *orthogonal*.

Any correlation between the selectivity mechanisms degrades orthogonality and reduces the efficiency of the 2D-LC system.

For strategies to achieve maximum orthogonality, refer to [Table 27](#) on page 333 and [Table 28](#) on page 338.

Resolution

A chromatographic separation can be optimized based on physical parameters of the HPLC column such as particle size, pore size, morphology of the particles, the length and diameter of the column, the solvent velocity, and the temperature. In addition, the thermodynamics of a separation can be considered and the properties of the solute and the stationary and mobile phases (percentage of organic solvent, ion strength, and pH) can be manipulated to achieve the shortest possible retention and highest selectivity.

1D-LC Resolution (R_S) can be described as a function of three parameters:

- Column efficiency or theoretical plates (N),
- Selectivity (α),
- Retention factor (k).

$$R_s = \frac{\sqrt{N}}{4} \left[\frac{\alpha - 1}{\alpha} \right] \left[\frac{k'_2}{k'_2 + 1} \right]$$

Figure 120 Resolution equation

This means that the selection of appropriate mobile and stationary phase properties and temperature is critical in achieving a successful separation.

Resolution in a one-dimensional separation usually is measured with:

$$R = \frac{\Delta t}{4\sigma}$$

R = Resolution
Δt = Difference in retention time maxima of two components
σ = Average standard deviation of two Gaussian peaks

Following results of this formula are important in practice:

- $R > 1.5$
Peaks are completely baseline resolved
- $R > 1$
Difference in retention time is larger than peak broadening, and therefore peak spacing is adequate to observe distinct component zones
- $R < 0.5$
Peaks are completely fused

2D-LC In 2D-LC the separation behaviour is more complex and described below.




Figure 121 Diagram of 2D resolution measurement: Slice for resolution (top) and 2-dimensional contour plot (bottom)

The distance between two spots in the contour plot may be calculated by the Pythagorean expression:

$$dr = \sqrt{dx^2 + dy^2}$$

For the resolution along the axis of each dimension applies:

$$R_1 = \frac{dy}{4\sigma_1}$$

and

$$R_2 = \frac{dx}{4\sigma_2}$$

So for two dimensions the resolution may be calculated as follows:

$$R_{2D} = \frac{dr}{4\sigma} = \sqrt{\left(\frac{dx}{4\sigma}\right)^2 + \left(\frac{dy}{4\sigma}\right)^2}$$

Figure 122 2D-Resolution (Pythagorean relation)

or, σ approximated by the average of σ_1 and σ_2 , using the easy to measure peak to valley ratio ($P = f/g$) and assuming that peaks are Gaussian:

$$Rs = \sqrt{-\frac{1}{2} \ln \left(\frac{1-P}{2} \right)}$$

Figure 123 2D-Resolution (peak to valley ratio relation)

Table 26 Definitions

Symbol	Denotation
R	Resolution
Δt	Difference in retention time maxima of two components
σ	Average standard deviation of Gaussian peaks
dr	Distance between two spots in a plane
P	Peak to valley ratio
f	Difference between amplitude at the valley, h, and g
h	Valley
g	Average peak maximum

Peak Capacity

Peak capacity may be differently defined:

- As the maximum number of peaks that can be resolved in the available separation space (*Geometrical Definition*), or
- As the ratio of the total area of the chromatogram to the area required for the resolution of any zone (*General Definition*)

Geometrical Definition

The peak capacity may be defined as the maximum number of peaks that can be resolved in the available separation space. So peak capacity n_c is related to the number of theoretical plates N :

$$n_c = PN^{1/2}$$

(P depends on the retention time range)

In practice peaks are usually not distributed randomly over the chromatogram and often overlap. Or in other words: In practice peaks don't fill the available separation space evenly. This is the reason, why the number of detectable components of a sample in 1D-LC is relatively small.

2D-LC separation offers an alternative possibility for increasing n_c : Orthogonal retention mechanisms generate a separation plane. Thus, the peak capacity in 2D-LC is the product of the peak capacities of the individual columns. Due to peak broadening in 1st and 2nd dimension, components in 2D-LC are present as two-dimensional ellipses on the retention plane.

How to calculate n_c depends on the method:

- For comprehensive 2D-LC:

$$n_c = \frac{L_1 L_2}{ab} = n_{c1} n_{c2}$$

L = Separation space for dimension

ab = Area for rectangle circumscribing the ellipse on the separation plane

- For heart-cutting 2D-LC:

$$n_c = \sum_{i=1}^k n_{ci}$$

**General
Definition**

Alternatively peak capacity may be defined as the ratio of the total area A of the chromatogram to the area A_0 required for the resolution of any zone:

$$n_{c,alternat} = \frac{A}{A_0}$$

n_c defined that way is related to the geometrical definition by a factor:

$$n_c = \frac{\pi}{4} n_{c,alternat} \approx 0.79 n_{c,alternat}$$

**Limits of Peak
Capacity in
2D-LC**

Under ideal circumstances (*orthogonality*), the overall peak capacity ($n_{c,2D}$) should be equal to the product of the individual peak capacities of the first and second dimension separations (1n_c and 2n_c)

$$n_{c,2D} = ^1n_c \times ^2n_c$$

In practice the increase in peak capacity is not directly proportional to increase in ability to resolve peaks.

Probable reason for this:

- In 1D-LC, with a baseline width of a single component peak $x_0 = 6\sigma$, x_0 units of component free space on both sides of the maxima is necessary to ensure baseline resolved peaks.
- In 2D-LC the single component zone is $A_0 = 2\pi r^2$ and an area of component free space of $\pi(2r)^2$.
- As a result: For every two component free widths in one dimension, four component free areas are required in two dimensions.

**Conclusions
for 2D-LC**

1D-LC is inadequate for the separation of complex mixtures, as the number of observable peaks compared to number of peaks to observe is too low. One theoretical model (Statistical Model of Overlap = SMO), that correlates well with real world observations, predicts, that the maximal fraction of the total peak capacity that can be seen as chromatographic peaks is 37 % and even only 18 % as single peaks. This implicates that extremely high peak capacities are needed to separate complex samples with lots of components which is extremely difficult to achieve.

Compared to 1D-LC separations, it's complicated to predict the number of observable peaks in 2D-LC. For example, at a given peak capacity and a given number of components, the aspect ratio in the two axes of separation has impact on how effective the two separation are.

From the practical point of view the performance between 1D- and 2D-LC should be compared, considering the following aspects:

- Peak capacity
- Number of peaks observed in experimental chromatograms

Ideal 2D Peak Capacity

One major problem in 2D-LC is loss of 1st dimension resolution due to 2nd dimension sampling process. The determining factors are:

- Gradient time of the 2nd dimension separation cannot exceed the sampling interval of the 1st dimension separation
- Resolution of a pair of peaks in the two-dimensional space is related to the resolution on the first and second dimensions as the Pythagorean average (see [Figure 122](#) on page 327)

A 2D chromatogram is only a way of displaying a lengthy series of sequential chromatograms obtained on the second column and the second column and detector are just a unique type of chemically selective detector of what comes out of the first column (see, "[2nd dimension as detector](#)" on page 331). The peak width observed on the second column is independent of the sampling time used in the 1st dimension.

This leads to two extreme scenarios, on how mixtures of components may behave:

- Unresolved mixture is injected into second column and second column separates analytes perfectly

$R_{s,2D}$ is independent of 1st dimension sampling rate

- Partially resolved mixture is injected into second column and analytes co-elute on the second column

$R_{s,2D}$ strongly depends on first dimension sampling rate.

This indicates, that it's very important to respect, how often the 1st dimension effluent must be sampled to avoid loss of resolution.

NOTE

The theoretical limits for ideal 2D peak capacity are defined by the Murphy-Schure-Foley Criterion (M-S-F sampling criterion). According to this criterion, the effluent must be sampled at least 3 – 4 times over 8σ width of the first dimension peak.

2nd dimension as detector

Functionally the second dimension of 2D-LC operates like a chemically sensitive detector for the peaks that elute from the first dimension column. Thus, 2D-LC may be understood as a three step process:

- 1st dimension separation (1)
- Sampling of the 1st dimension (2)
- 2nd dimension separation and detection (3)




Figure 124 Diagram of instrumentation for 2D-LC




Figure 125 Principle of 2D-LC (example for LCxLC): Effluent of first column (1) is sampled (2) and injected to second column (3). Peaks of second column separation are detected and reconstructed.

- | | |
|-----|---------------------------------|
| (1) | First dimension separation |
| (2) | Sampling of the first dimension |
| (3) | Second dimension separation |

Successful Mode Combinations

2D-LC separations are the more effective, the more the selectivity mechanisms involved in the two stages differ. Completely different and independent mechanisms are said to be orthogonal. Any correlation between the selectivity mechanisms degrades orthogonality and reduces the efficiency of the 2D-LC system.

Thus, selecting the best combination of stationary and mobile phase is the major issue to improve 2D-LC methods. [Table 27](#) on page 333 summarizes the advantages and disadvantages of combinations of normal phase (NP), reverse phase (RP), ionexchange (IEC) and size exclusion chromatography (SEC) for 2D-LC operation.

Table 27 Mode combinations in 2D-LC (LCxLC)

Combination	Orthogonality	Peak capacity	Application	Comment
RP x RP	1	++ ²	Peptidomics, metabolomics, pharmaceuticals, foods, cosmetics	Miscible solvents, broadest application, fast speed, gradient elution on both dimensions
IEC and RP	+ ³	-	Proteomics, peptidomics	
SEC and RP	+	4 ₋	Polymers, proteomics	
NP and RP	+		Polymers, pharmaceuticals, oils	Solvent incompatibility, limited application
Affinity and RP	+	-	Proteomics	
SEC and NP	+	-	Polymers	
SEC and IEC	+	-	Proteomics	

¹ Orthogonality, depends on the column choice or mobile phase choice

² very good

³ good

⁴ not so good

Solvent Elution Modes

Table 28 on page 338 focuses on the effects of elution modes for second dimension separation.

The following elution modes for second dimension separation are commonly used:

- Gradient

A standard gradient of solvent A vs. solvent B for the second dimension separation will be repeated during the complete first dimension separation




Figure 126 Standard gradient mode

- Shifted Gradient

From each second dimension separation to the next the start-%B and end-%B values of the individual second dimension gradients will be increased in a defined way. Additionally, the gradient span can be increased from each second dimension gradient to the next.




Figure 127 Shifted gradient mode with increase of start-%B

- Isocratic

All second dimension separations will be carried out in an isocratic mode.




Figure 128 Isocratic mode

- Advancing isocratic

Nearly isocratic conditions are used in each second dimension separation, with slightly increasing solvent strength in each successive run.

The second dimension pumping system is fed with a shallow gradient in eluent composition over the course of the 2D-separation.




Figure 129 Advancing isocratic mode

Table 28 Different elution modes in the 2nd dimension (pros and cons)

Criterion	Gradient/Shifted gradient	Isocratic/Advancing isocratic
Peak capacity	Superior	Inferior
Diversity of samples (complex samples)	Superior	Inferior
Baseline performance (sensitivity)	Inferior (baseline drift caused by solvent gradient)	Superior
Pressure stress (column lifetime!)	Inferior (large changes within every 2nd dimension gradient)	Superior (no pressure changes with isocratic, gradually changing with advancing isocratic)

All modes are easily available with the Agilent 2D-LC Acquisition software.

Each mode has advantages and disadvantages. No single mode is superior in all applications of 2D-LC.


Effect of shifted gradient elution mode in the 2nd dimension

Figure 130 2nd dimension gradient mode compared to isocratic mode and its effect on resolution

α_{sg} as achieved in shifted gradient mode is larger than α_g achieved in standard gradient elution mode. This can lead to an improved peak detection and improved separation.

See D. Li and O. J. Schmitz "Use of Shift Gradient in the Second Dimension to Improve the Separation Space in Comprehensive Twodimensional Liquid Chromatography" *Anal. Bioanal. Chem.* 405, 6511-6517 (2013)

Practical Issues

The table below gives an overview, which practical issues have to be considered in 2D-LC.

Table 29 Practical issues in 2D-LC

Issue	Theoretical base	Comment
Choice of first dimension column diameter	Has impact on trade off between optimum first dimension flow rate and amount of sample injected into the second dimension column for each second column run	
Ratio of column diameter in the two dimensions	Causes significant analyte dilution effects	
Goals of the analytical method	Chosen parameter depend on what is important in analysis: • separate as many constituents as possible or • focused on resolution and quantitation of a specific constituent	True gradient elution in the second dimension separation provides better peak capacity than in isocratic elution. Gradient elution is the best available mechanism for achieving peak focusing.
Selection of the stationary phases and column formats	For RPLC in both dimensions the retentivity of the second dimension column must be much higher than that of the first dimension column required because: • a relatively large volume of the sample will be collected and injected into the second column • to minimize peak broadening the sample should be focused at the inlet of the second column	

Based on theory, in most cases following approaches to achieve best possible 2D-LC should be respected:

- Methodology

As in Comprehensive 2D-LC is no direct need¹ for UV-detection in the first dimension, other eluents than acetonitril or methanol are possible. This implies the possibility to use unconventional organic solvents in the first dimension.

NOTE

Take care when using any unconventional organic solvents that these are still compatible with the used instrumentation. In doubt, refer to the module documentation or call Agilent.

- Instrumentation

It is important to use very low delay-volume-gradient pumping systems that are able to produce high flow rates to achieve fast second dimension gradients with only little gradient delay - like the Agilent 1290 Infinity LC.

- Columns

Total orthogonality is difficult to achieve, as there are relatively few combinations sufficiently phase selective.

- Detection methods

Compared to mass spectrometry DAD based UV detection is faster, cheaper and offers higher reproducibility, thus mass spectrometry offers additional increase in peak capacity by expanding the separation space into the MS-domain. A high sensitivity UV-detector is recommended since a dilution of the first dimension peaks occurs in the second dimension separation – an Agilent 1260 or 1290 Infinity Diode-Array-Detector with 60 mm flow cell is ideal as second dimension detector.

- Data analysis

2D-LC-data are complex. Use of special software is advisable.

¹ In case the peak and time triggered operation of the second dimension separation, which is optionally available with the Agilent 1290 Infinity 2D-LC solution, an UV-detector is required between the first dimension column and the modulation valve.

General Safety Information	343
General Safety Information	343
Safety Standards	343
General	344
Before Applying Power	344
Ground the Instrument	344
Do Not Operate in an Explosive Atmosphere	345
Do Not Remove the Instrument Cover	345
Do Not Modify the Instrument	346
In Case of Damage	346
Solvents	347
Safety Symbols	348
Waste Electrical and Electronic Equipment (WEEE) Directive	350
Radio Interference	351
Sound Emission	352
Solvent Information	353
Further Information	354
Agilent Technologies on Internet	355

This chapter provides addition information on safety, legal and web.

General Safety Information

General Safety Information

The following general safety precautions must be observed during all phases of operation, service, and repair of this instrument. Failure to comply with these precautions or with specific warnings elsewhere in this manual violates safety standards of design, manufacture, and intended use of the instrument. Agilent Technologies assumes no liability for the customer's failure to comply with these requirements.

WARNING

Ensure the proper usage of the equipment.

The protection provided by the equipment may be impaired.

- ✓ **The operator of this instrument is advised to use the equipment in a manner as specified in this manual.**

Safety Standards

This is a Safety Class I instrument (provided with terminal for protective earthing) and has been manufactured and tested according to international safety standards.

General

Do not use this product in any manner not specified by the manufacturer. The protective features of this product may be impaired if it is used in a manner not specified in the operation instructions.

Before Applying Power

WARNING

Wrong voltage range, frequency or cabling

Personal injury or damage to the instrument

- ✓ Verify that the voltage range and frequency of your power distribution matches to the power specification of the individual instrument.
- ✓ Never use cables other than the ones supplied by Agilent Technologies to ensure proper functionality and compliance with safety or EMC regulations.
- ✓ Make all connections to the unit before applying power.

NOTE

Note the instrument's external markings described under ["Safety Symbols"](#) on page 348.

Ground the Instrument

WARNING

Missing electrical ground

Electrical shock

- ✓ If your product is provided with a grounding type power plug, the instrument chassis and cover must be connected to an electrical ground to minimize shock hazard.
- ✓ The ground pin must be firmly connected to an electrical ground (safety ground) terminal at the power outlet. Any interruption of the protective (grounding) conductor or disconnection of the protective earth terminal will cause a potential shock hazard that could result in personal injury.

Do Not Operate in an Explosive Atmosphere

WARNING

Presence of flammable gases or fumes

Explosion hazard

- ✓ Do not operate the instrument in the presence of flammable gases or fumes.
-

Do Not Remove the Instrument Cover

WARNING

Instrument covers removed

Electrical shock

- ✓ Do Not Remove the Instrument Cover
 - ✓ Only Agilent authorized personnel are allowed to remove instrument covers. Always disconnect the power cables and any external circuits before removing the instrument cover.
-

Do Not Modify the Instrument

Do not install substitute parts or perform any unauthorized modification to the product. Return the product to an Agilent Sales and Service Office for service and repair to ensure that safety features are maintained.

In Case of Damage

WARNING

Damage to the module

Personal injury (for example electrical shock, intoxication)

- ✓ Instruments that appear damaged or defective should be made inoperative and secured against unintended operation until they can be repaired by qualified service personnel.
-

Solvents

WARNING

Toxic, flammable and hazardous solvents, samples and reagents

The handling of solvents, samples and reagents can hold health and safety risks.

- ✓ When working with these substances observe appropriate safety procedures (for example by wearing goggles, safety gloves and protective clothing) as described in the material handling and safety data sheet supplied by the vendor, and follow good laboratory practice.
- ✓ Do not use solvents with an auto-ignition temperature below 200 °C (392 °F). Do not use solvents with a boiling point below 56 °C (133 °F).
- ✓ Avoid high vapor concentrations. Keep the solvent temperature at least 40 K below the boiling point of the solvent used. This includes the solvent temperature in the sample compartment. For the solvents methanol and ethanol keep the solvent temperature at least 25 K below the boiling point.
- ✓ Do not operate the instrument in an explosive atmosphere.
- ✓ Do not use solvents of ignition Class IIC according IEC 60079-20-1 (for example, carbon disulfide).
- ✓ Reduce the volume of substances to the minimum required for the analysis.
- ✓ Never exceed the maximum permissible volume of solvents (8 L) in the solvent cabinet. Do not use bottles that exceed the maximum permissible volume as specified in the usage guideline for solvent cabinet.
- ✓ Ground the waste container.
- ✓ Regularly check the filling level of the waste container. The residual free volume in the waste container must be large enough to collect the waste liquid.
- ✓ To achieve maximal safety, regularly check the tubing for correct installation.

NOTE

For details, see the usage guideline for the solvent cabinet. A printed copy of the guideline has been shipped with the solvent cabinet, electronic copies are available in the Agilent Information Center or via the Internet.

Safety Symbols

Table 30 Symbols

	The apparatus is marked with this symbol when the user shall refer to the instruction manual in order to protect risk of harm to the operator and to protect the apparatus against damage.
	Indicates dangerous voltages.
	Indicates a protected ground terminal.
	The apparatus is marked with this symbol when hot surfaces are available and the user should not touch it when heated up.
	Sample Cooler unit is designed as vapor-compression refrigeration system. Contains fluorinated greenhouse gas (refrigerant) according to the Kyoto protocol. For specifications of refrigerant, charge capacity, carbon dioxide equivalent (CDE), and global warming potential (GWP) see instrument label.
	Flammable Material For Sample Thermostat which uses flammable refrigerant consult Agilent Information Center / User Manual before attempting to install or service this equipment. All safety precautions must be followed.
	Confirms that a manufactured product complies with all applicable European Community directives. The European Declaration of Conformity is available at: http://regulations.corporate.agilent.com/DoC/search.htm
	Manufacturing date.
	Power symbol indicates On/Off. The apparatus is not completely disconnected from the mains supply when the power switch is in the Off position
	Pacemaker Magnets could affect the functioning of pacemakers and implanted heart defibrillators. A pacemaker could switch into test mode and cause illness. A heart defibrillator may stop working. If you wear these devices keep at least 55 mm distance to magnets. Warn others who wear these devices from getting too close to magnets.

Table 30 Symbols

	Magnetic field Magnets produce a far-reaching, strong magnetic field. They could damage TVs and laptops, computer hard drives, credit and ATM cards, data storage media, mechanical watches, hearing aids and speakers. Keep magnets at least 25 mm away from devices and objects that could be damaged by strong magnetic fields.
	Indicates a pinching or crushing hazard
	Indicates a piercing or cutting hazard.

WARNING

A WARNING

alerts you to situations that could cause physical injury or death.

- ✓ Do not proceed beyond a warning until you have fully understood and met the indicated conditions.

CAUTION

A CAUTION

alerts you to situations that could cause loss of data, or damage of equipment.

- ✓ Do not proceed beyond a caution until you have fully understood and met the indicated conditions.

Waste Electrical and Electronic Equipment (WEEE) Directive

This product complies with the European WEEE Directive marking requirements. The affixed label indicates that you must not discard this electrical/electronic product in domestic household waste.



NOTE

Do not dispose of in domestic household waste

To return unwanted products, contact your local Agilent office, or see <http://www.agilent.com> for more information.

Radio Interference

Never use cables other than the ones supplied by Agilent Technologies to ensure proper functionality and compliance with safety or EMC regulations.

Test and Measurement

If test and measurement equipment is operated with equipment unscreened cables and/or used for measurements on open set-ups, the user has to assure that under operating conditions the radio interference limits are still met within the premises.

Sound Emission

Manufacturer's Declaration

This statement is provided to comply with the requirements of the German Sound Emission Directive of 18 January 1991.

This product has a sound pressure emission (at the operator position) < 70 dB.

- Sound Pressure L_p < 70 dB (A)
- At Operator Position
- Normal Operation
- According to ISO 7779:1988/EN 27779/1991 (Type Test)

Solvent Information

Observe the following recommendations on the use of solvents.

- Brown glass ware can avoid growth of algae.
- Avoid the use of the following steel-corrosive solvents:
 - solutions of alkali halides and their respective acids (for example, lithium iodide, potassium chloride, and so on),
 - high concentrations of inorganic acids like sulfuric acid and nitric acid, especially at higher temperatures (if your chromatography method allows, replace by phosphoric acid or phosphate buffer which are less corrosive against stainless steel),
 - halogenated solvents or mixtures which form radicals and/or acids, for example:
$$2\text{CHCl}_3 + \text{O}_2 \rightarrow 2\text{COCl}_2 + 2\text{HCl}$$
This reaction, in which stainless steel probably acts as a catalyst, occurs quickly with dried chloroform if the drying process removes the stabilizing alcohol,
- chromatographic grade ethers, which can contain peroxides (for example, THF, dioxane, diisopropyl ether) should be filtered through dry aluminium oxide which adsorbs the peroxides,
- solvents containing strong complexing agents (e.g. EDTA),
- mixtures of carbon tetrachloride with 2-propanol or THF.
- Avoid the use of dimethyl formamide (DMF). Polyvinylidene fluoride (PVDF), which is used in leak sensors, is not resistant to DMF.

Further Information

Further information is available:

- Folder Documents on the software DVD:
 - Document Primer 2D-LC 5991-2359EN.pdf gives an introduction to principles, practical implementation and applications for Two-Dimensional Liquid Chromatography.
- Folder Documentation of the Agilent OpenLab 2D-LC Software CD:
 - Technical Note G2198-90101 describes software ²D Chromatogram Creator for MassHunter.
- Folder Documents on the Driver CD:
 - Software Status Bulletin (SSB)
The SSB is updated regularly. Please visit our Website for the latest version at
<http://www.agilent.com/cs/library/Support/Patches/SSBs/2D-LC.html>.
 - Software Release Bulletin (SRB)
The SRB is an excerpt from the SSB which lists issues which have been fixed with this revision.
- Firmware and firmware documentation are available for download from
http://www.chem.agilent.com/_layouts/agilent/downloadFirmware.aspx?wid=69761.
- Press **F1** in the software user interface for the Online Help with more information on specific software functions.
- For more information about applications, please visit InfinityLab Application Finder
<https://www.agilent.com/en/promotions/applicationfinder?s=learnmore>.
- For more information about Agilent hardware, and software, please visit the Agilent web site at <http://www.agilent.com>.

Agilent Technologies on Internet

For the latest information on products and services visit our worldwide web site on the Internet at:

<http://www.agilent.com>

In This Book

The manual describes the following:

- introduction,
- installing,
- configuring,
- using,
- data analysis,
- safety and related information.

www.agilent.com

© Agilent Technologies Inc. 2012-2020

Published in Germany
02/2020



Part No: G2198-90001 Rev. B
Document No: SD-29000383 Rev. B

