

Maximum throughput and lowest cost for population-scale whole-genome sequencing.

HiSeq X System Specifications

Parameter	Specification ^a
Instrument Configuration	Computer and touch screen display Installation setup and accessories Data collection and analysis software
Instrument Control Computer	Base Unit: 2× Intel Xeon E5-2697V2 64 BIT 2.7 GHz CPU 30 MB Cache Memory: 128 GB RAM Hard Drive: 8 × 1.0 TB SATA 7200 RPM 3.5" Constellation ES.3 (2 RAID 0 drives, 6 RAID 5 drives) Solid-State Drive: 5 × 400 GB (5 RAID 0 drives) Operating System: Microsoft Windows 7 Professional Note: Computer specifications will be regularly upgraded. Contact your local account manager for current configuration.
Operating Environment	Temperature: 22°C ± 3°C Humidity: Noncondensing 20–80% Altitude: Less than 2000 m (6500 ft) Air Quality: Pollution degree rating of II Ventilation: Maximum of 4000 BTU/h For Indoor Use Only.
Laser	532 nm, 660 nm, 650 nm (barcode reader)
Dimensions	W × D × H: 118.6 cm × 76.0 cm × 94.0 cm (46.7 in × 30.0 in × 37.0 in) Weight: 225.9 kg (498 lbs) Crated Weight: 316.6 kg (698 lbs)
Power Requirements	100–240V AC 50/60Hz, 20A, 1500 W Illumina provides a region-specific uninterruptible power supply for all HiSeq instruments
Product Safety	CE-marked and ETL-listed instrument
a. Specifications refer to an individual HiSeq X System.	

Ordering Information

Product	Catalog No.
HiSeq X Ten Products	
HiSeq X Ten System ^a	SY-412-1001
HiSeq X Ten Reagent Kit v2.5	FC-501-2501
HiSeq X Ten Reagent Kit v2.5 - 10 pack	FC-501-2521
HiSeq X Five Products	
HiSeq X Five System ^b	SY-412-1011
HiSeq X Five Reagent Kit v2.5	FC-502-2501
HiSeq X Five Reagent Kit v2.5 - 10 pack	FC-502-2521
a. Catalog No. refers to an individual system. Minimum order for the HiSeq X Ten is 10 systems. b. Catalog No. refers to an individual system. Minimum order for the HiSeq X Five is 5 systems.	

Illumina • 1.800.809.4566 toll-free (U.S.) • +1.858.202.4566 tel • techsupport@illumina.com • www.illumina.com

For Research Use Only. Not for use in diagnostic procedures.

© 2016 Illumina, Inc. All rights reserved.
Illumina, BaseSpace, HiSeq, HiSeq X, Isaac, TruSeq, and the pumpkin orange color are trademarks of Illumina, Inc. and/or its
affiliate(s) in the U.S. and/or other countries. All other names, logos, and trademarks are the property of their respective owners
Pub. No. 770-2014-056 Current as of 22 March 2016



AGAATGATAACAGTAACACACTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAGAGCTACCGTCTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAAC
TCAACGTACCGTAACGAACGTATCATTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAGAGCTACCGTGCAACGACGAAAAAGATGATAACAGTAACACACTTCTGTTAACCTT
CGACGAAAAGATGATAACAGTAACACACTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAGAGCTAC
CAACGTACCATTAAGAGCTACCGTGCAACGAGTAACACACTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAGAGCTACCGTGCAACGACGAAAAAGATGATAA
AGAATGATAACAGTAACACACTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAGAGCTACCGTCTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAAC
GATTACTTGATCCACTGATTCAACGTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACGTATCAATTGAGCTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACGTATCAATTGAGACTAGCAACGAC
CGTATCAATTGAGACTAAATATTAACGTACCATTAAGAGTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAGAGCTACCGTGCAACGAAAAAGATGATAACAG