



# chemiSOLO

Personal Western Blot Imager





# Highly Sensitive Imaging


Femtogram detection of  
chemiluminescent Westerns





# Unique Software


Capture and download images  
instantly from any device.






# Flexible Imaging Options


Chemiluminescent blots, colorimetric SDS PAGE gels, colorimetric blots, color marker, densitometry




Chemiluminescent Blot  
with Color Marker




Coomassie Stained  
Protein Gel



Silver Stained  
Protein Gel



Colorimetric  
Stained Blot




Densitometry



Plant  
Bioluminescence

# Quantitative Data Acquisition

chemiSOLO's Extended Dynamic Range linearly detects wider ranges without saturation.



Scan for more information on  
chemiSOLO's quantitative imaging





# Advanced Technology

## Embedded LEDs

Used for capturing marker and colorimetric gel images

## Two USB Ports

Used for Internet connection or direct data transfer

## One Ethernet Port

Used for direct connection via Ethernet cable

## Imaging Stage

Field of View: 10.0 x 15.0 cm (4.0 x 5.9 in)

## 6.29 MPx Camera

Powerful back-illuminated Peltier-cooled CMOS camera

## Browser-Based Software

Unique Wi-Fi signal allows the chemiSOLO to be controlled by any device with Wi-Fi connectivity and a browser





# Space-efficient Design

## Darkroom

At least 8 x 6ft



## Gel Tank and Power Supply

16.5 x 22 x 13in  
42 x 56 x 33cm

## chemiSOLO

11.5 x 17 x 8.75in  
29.2 x 43.2 x 22.2cm

## Microwave

22 x 15 x 18in  
56 x 38 x 45cm



| chemiSOLO             | Specifications                                                                                         |
|-----------------------|--------------------------------------------------------------------------------------------------------|
| Part Number           | CS1000                                                                                                 |
| Camera                | 6.29MPx back-illuminated CMOS, 16-bit, 65,536 grayscale                                                |
| Peltier Cooling       | -40C regulated cooling from ambient                                                                    |
| Applications          | Chemiluminescent Imaging with Visible Marker, Densitometry, Visible Gel Imaging in Color and Grayscale |
| Light Source          | LEDs                                                                                                   |
| Field of View         | 10 x 15 cm (4.0 x 5.9 in)                                                                              |
| Footprint (W x H x D) | 29.2 x 22.2 x 43.2 cm (11.5 x 8.75 x 17.0 in)                                                          |
| Weight                | 9.1kg (20lbs)                                                                                          |
| Certifications        | CE RED, cTUVus, CB scheme, RoHS                                                                        |

Scan for more information  
on the chemiSOLO



6747 Sierra Court, Suites A-B • Dublin, CA 94568 USA  
 Phone: 925.307.7127 • Fax: 925.905.1816  
[info@azurebiosystems.com](mailto:info@azurebiosystems.com) • [www.azurebiosystems.com](http://www.azurebiosystems.com)

Copyright © 2022 Azure Biosystems. All rights reserved. The Azure Biosystems logo and Azure Biosystems® are trademarks of the Company. All other trademarks, service marks and trade names appearing in this brochure are the property of their respective owners.

MM-0087 Rev2