

## Technical Data

### Power Supply

|                          |                 |
|--------------------------|-----------------|
| 50/60Hz (max. load 9.8A) | 1 x 100V / 120V |
| 50/60Hz (max. load 4.9A) | 1 x 220V / 240V |

### Output voltage/current

|                  |          |                      |
|------------------|----------|----------------------|
| Polishing        | 0 - 100V | (1V steps) / 6A      |
| Etching          | 0 - 25V  | (0.5 V steps) / 6A   |
| External Etching | 0 - 15V  | (0.5 V steps) / 1.5A |

### Software and Electronics

|          |                                                               |
|----------|---------------------------------------------------------------|
| Display  | 128 x 240 dots (16 x 40 characters)                           |
| Controls | Touch pad                                                     |
| Database | 10 Struers methods + 20 user definable methods (non volatile) |

## Dimensions and Weight

|        | Control unit   | Polishing unit  |
|--------|----------------|-----------------|
| Width  | 385 mm / 15.2" | 220 mm / 8.7"   |
| Depth  | 350 mm / 13.8" | 350 mm / 13.8"  |
| Height | 160 mm / 6.3"  | 160 mm / 6.3"   |
| Weight | 18 kg / 40 lbs | 6.3 kg / 14 lbs |

## Specifications

### LectroPol-5

Fully automatic, microprocessor controlled electrolytic polishing and etching apparatus. Complete with control unit, polishing unit (04936201) and external etching unit (02250121). One set of masks (04936915) and two electrolyte containers (04936902) included.

### Cat. no.

Control unit with external etching unit 04936233

Polishing unit with one set of masks (04936915) and two electrolyte containers (04936902) 04936201

Polishing unit for low temperature, with one set of masks (04936915), one cathode (04936908) and two electrolyte containers (04936902). External cooling unit and tubing for connection to the external cooling unit are not included and have to be sourced alternatively. Internal dia. of tubing to be used, 6 mm 04936202

Extension Arm, for the preparation of specimens up to 100 mm high. Can be mounted on pol. unit 04936202 and on pol.unit 04936201 from serial no. 4935250 04936909

1 set of masks (0.5, 1, 2, 5 cm<sup>2</sup> and 1 without aperture.) 04936915

1 set of 5 masks 0.5 cm<sup>2</sup> 04936916

1 set of 5 masks 1 cm<sup>2</sup> 04936917

1 set of 5 masks 2 cm<sup>2</sup> 04936918

1 set of 5 masks 5 cm<sup>2</sup> 04936919

1 set of 5 masks without aperture 04936920

Electrolyte container with lid 04936902

External etching unit 02250121

## List of Electrolytes

| A2              | A3                | D2     |
|-----------------|-------------------|--------|
| Aluminium       | Martensitic Steel | Brass  |
| Iron            | Molybdenum        | Copper |
| Nickel          | Stainless Steel   | Gold   |
| Silver          | Titanium          |        |
| Stainless Steel | Vanadium          |        |
| Steels          | Zirconium         |        |
| Tin             |                   |        |
| Titanium        |                   |        |

*Struers' equipment is in conformity with the provisions of the applicable International Directives and their appurtenant Standards. (Please contact your local supplier for details)*

*Struers' products are subject to constant product development. Therefore, we reserve ourselves the right to introduce changes in our products without notice.*



### Struers A/S

Pederstrupvej 84  
DK-2750 Ballerup, Denmark  
Phone +45 44 600 800  
Fax +45 44 600 801  
struers@struers.dk  
www.struers.com

### AUSTRALIA & NEW ZEALAND

**Struers Australia**  
27 Mayneview Street  
Milton QLD 4064  
Australia  
Phone +61 7 3512 9600  
Fax +61 7 3369 8200  
info.au@struers.dk

### BELGIQUE (Wallonie)

**Struers S.A.S.**  
370, rue du Marché Rollay  
F-94507 Champigny  
sur Marne Cedex  
Téléphone +33 1 5509 1430  
Télécopie +33 1 5509 1449  
struers@struers.fr

### BELGIUM (Flanders)

**Struers GmbH Nederland**  
Elektraweg 5  
3144 CB Maassluis  
Telefoon +31 (10) 599 7209  
Fax +31 (10) 599 7201  
netherlands@struers.de

### CANADA

**Struers Ltd.**  
7275 West Credit Avenue  
Mississauga, Ontario L5N 5M9  
Phone +1 905-814-8855  
Fax +1 905-814-1440  
info@struers.com

### CHINA

**Struers Ltd.**  
No. 1696 Zhang Heng Road  
Zhang Jiang Hi-Tech Park  
Shanghai 201203, P.R. China  
Phone +86 (21) 6035 3900  
Fax +86 (21) 6035 3999  
struers@struers.cn

### CZECH REPUBLIC

**Struers GmbH**  
Organizační složka  
Havlíčková 361  
CZ-252 63 Roztoky u Prahy  
Phone +420 233 312 625  
Fax +420 233 312 640  
czechrepublic@struers.de

### DEUTSCHLAND

**Struers GmbH**  
Carl-Friedrich-Benz-Straße 5  
D-47877 Willich  
Telefon +49 (0) 2154 486-0  
Fax +49 (0) 2154 486-222  
verkauf@struers.de

### FRANCE

**Struers S.A.S.**  
370, rue du Marché Rollay  
F-94507 Champigny  
sur Marne Cedex  
Téléphone +33 1 5509 1430  
Télécopie +33 1 5509 1449  
struers@struers.fr

### HUNGARY

**Struers GmbH**  
Magyarországi Fióktelep  
Tatai út 53  
2821 Gyermely  
Phone +36 (34) 880546  
Fax +36 (34) 880547  
hungary@struers.de

### IRELAND

**Struers Ltd.**  
Unit 11 Evolution@ AMP  
Whittle Way, Catcliffe  
Rotherham S60 5BL  
Tel. +44 0845 604 6664  
Fax +44 0845 604 6651  
info@struers.co.uk

### ITALY

**Struers Italia**  
Via Monte Grappa 80/4  
20020 Arese (MI)  
Tel. +39-02/38236281  
Fax +39-02/38236274  
struers.it@struers.it

### JAPAN

**Marumoto Struers K.K.**  
Takara 3rd Building  
18-6, Higashi Ueno 1-chome  
Taiko-ku, Tokyo 110-0015  
Phone +81 3 5688 2914  
Fax +81 3 5688 2927  
struers@struers.co.jp

### NETHERLANDS

**Struers GmbH Nederland**  
Elektraweg 5  
3144 CB Maassluis  
Telefoon +31 (10) 599 7209  
Fax +31 (10) 599 7201  
netherlands@struers.de

### NORWAY

**Struers A/S, Norge**  
Sjøskegveien 44C  
1407 Vinterbro  
Telefon +47 970 94 285  
info@struers.no

### ÖSTERREICH

**Struers GmbH**  
Zweigniederlassung Österreich  
Betriebsgebiet Puch Nord 8  
5412 Puch  
Telefon +43 6245 70567  
Fax +43 6245 70567-78  
austria@struers.de

### POLAND

**Struers Sp. z o.o.**  
Oddział w Polsce  
ul. Jasnowiecka 44  
31-358 Kraków  
Phone +48 12 661 20 60  
Fax +48 12 626 01 46  
poland@struers.de

### ROMANIA

**Struers GmbH**  
Sucursala Sibiu  
Str.Scoala de Inot, nr. 18  
RO-550005 Sibiu  
Phone +40 269 244 558  
Fax +40 269 244 559  
romania@struers.de

### SCHWEIZ

**Struers GmbH**  
Zweigniederlassung Schweiz  
Weissenbrunnstraße 41  
CH-8903 Birmensdorf  
Telefon +41 44 777 63 07  
Fax +41 44 777 63 09  
switzerland@struers.de

### SINGAPORE

**Struers Singapore**  
627A Aljunied Road,  
#07-08 BizTech Centre  
Singapore 38942  
Phone +65 6299 2268  
Fax +65 6299 2661  
struers.sg@struers.dk

### SPAIN

**Struers A/S, Spain**  
Camino Cerro de los Gamos 1  
Building 1 - Pozuelo de Alarcón  
CP 28224 Madrid  
Teléfono +34 917 901 204  
Fax +34 917 901 112  
struers.es@struers.es

### SUOMI

**Struers Suomi**  
Kalevankatu 43  
00180 Helsinki  
Puhelin +358 (0)207 919 430  
Faksi +358 (0)207 919 431  
finland@struers.fi

### SWEDEN

**Struers Sverige**  
Ekbacksvägen 22  
168 69 Bromma  
Telefon +46 (0)8 447 53 90  
Telefax +46 (0)8 447 53 99  
info@struers.se

### UNITED KINGDOM

**Struers Ltd.**  
Unit 11 Evolution @ AMP  
Whittle Way, Catcliffe  
Rotherham S60 5BL  
Tel. +44 0845 604 6664  
Fax +44 0845 604 6651  
info@struers.co.uk

### USA

**Struers Inc.**  
24766 Detroit Road  
Westlake, OH 44145-1598  
Phone +1 440 871 0071  
Fax +1 440 871 8188  
info@struers.com