


PCR setup STARlet

Automated PCR assay setup

for e.g. endpoint-, real-time-, RT- or multiplex PCR



- ▶ Highly flexible
- ▶ Easy to use
- ▶ Resource saving
- ▶ Reliable

Discover the possibilities

The Hamilton PCR setup STARlet is suitable for a wide range of applications such as gene expression, genotyping, sequencing or pathogen detection.

It has been developed for 96/384 PCR plates but different labware types such as strips or even rotor discs can be used.

The outstanding deck capacity of 25 SBS plate positions compared to the small system footprint is unique in this category; and with its available configurations from 4 to 8 independent 1 ml channels the PCR setup STARlet serves a broad range of needs and throughputs, from simple to highly complex assay setups.

Highly flexible

- ▶ **PCR setup software:** enables the flexible combination of PCR components, from simple two step- (e.g. template + primer MasterMix) to complex setups with up to 7 pipetting steps
- ▶ **Sample normalization:** varying sample input concentrations can be normalized to required target concentration
- ▶ **Multiplexing:** up to 4 primer mixes can be pipetted per reaction well
- ▶ **MasterMix preparation:** the multiplexing function can be used for simple "MasterMix preparation in reaction" with up to 6 components
- ▶ **Reagents/labware:** different PCR kits and labware formats can be used, there is no supplier limitation
- ▶ **LIMS compatibility:** the assay setup information is generated with a GUI but it can also be imported from a LIMS system

Easy to use

- ▶ **Visualized setup:** a Graphical User Interface (GUI) easily guides through the assay setup
- ▶ **Quick PCR setup function:** allows fast assay setup, higher throughput and decreased processing times
- ▶ **Advanced PCR setup function:** simplifies complex setups with up to 7 PCR components including multiplexing function and free plate design
- ▶ **Reusable protocols:** setups can be saved and reused again. Partial reuse is possible e.g. replacing samples while leaving remaining setup and plate design unchanged


Resource saving

- ▶ **Low volume PCR:** precise pipetting of small volumes for low volume reactions in e.g. 384 PCR plate format
- ▶ **Reagent saving pipetting:** with low dead volumes expensive reagents can be saved
- ▶ **Time saving:** reduced hands-on and processing times due to high deck capacity and optimized pipetting settings

Reliable

- ▶ **Accurate reaction setup:** highly precise liquid handling with Hamilton Monitored Air Displacement pipetting technology
- ▶ **High traceability:** optional tracking of all input components by barcode scanning; tracing complete setup information throughout the process
- ▶ **Verified configurations:** tested default configurations are included in delivery


For further information please visit www.hamiltonrobotics.com or contact your local distributor.




Easily choose function: quick- / advanced PCR assay setup or sample normalization



Either use predefined or your desired labware



Design PCR plates row-/column wise or completely free



Create your own GUI environment by renaming input terms e.g. "templates" instead of "samples"

© 2012 Hamilton Bonaduz AG. All rights reserved.
Lit. No. F-1210-02/00 QTY: 200, 10/12 Printed in Germany

HAMILTON

Web: www.hamiltonrobotics.com

USA: 800-648-5950

Email: infoservice@hamiltonrobotics.com

United States
Tel: +1-775-858-3000

United Kingdom & Ireland
Tel: +44 (0)121-717-0199

Brazil
Tel: +55 (11) 9677-4093

China
Tel: +86-21-6164-6567

France
Tel: +33 (01) 69751616

Italy
Tel: +39-39-689-33-93

Denmark, Norway,
Sweden, Finland
Tel: +46 (0) 8 410 273 73

Germany, Switzerland,
Austria, Benelux
Tel: +49 (0) 89 552649-0

To find a subsidiary or distributor in your area, please visit hamiltonrobotics.com/contacts.