

Guided Resource Map

Agilent 708-DS Dissolution Apparatus



One platform for all dissolution needs

The Agilent 708-DS Dissolution Apparatus provides a versatile and robust platform for standardizing dissolution testing in your lab. This instrument delivers consistent performance and results and complies with current pharmacopeia specifications.

The modular design of the 708-DS supports you in a wide range of dissolution tests, from simple manual tests to automated configurations for high-throughput sampling and online analysis. You can configure the 708-DS for 6- or 8-position testing and choose from a variety of vessel sizes and accessories based on the dosage form being tested. In addition, the 708-DS also offers you advanced usage options, including automated dosage delivery, in-vessel temperature monitoring, automated sampling, and in-situ fiber-optic analysis.



Paddles

- 14-3603 Lower interchangeable paddle, electropolished stainless steel, 1 L
- 14-3602 Lower interchangeable paddle, PTFE-coated, 1 L



Thermal Paper

- G7910-87002 Thermal printer paper for use with 708-DS, 10 rolls/pk



Filters

- 17-4000 Full Flow Filter UHMWPE, blue, 10 µm, 100/pk
- 17-4005 Full Flow Filter UHMWPE, blue, 10 µm, 1000/pk



Basket Shafts

- 14-3631 Basket shaft, lower interchangeable, 3-clip, USP



Vessels

- 12-5149 TruAlign vessel, 1 L, clear
- 12-5149V TruAlign vessel, 1 L, clear, including certificate
- 12-1450 TruAlign vessel, 250 mL conversion kit, for Chinese Pharmacopeia (includes vessel, adapter, two-piece mini paddle, evaporation cover, and manifold adapter), set of 6




Shafts

- 14-3613 Upper receptor shaft, 21 inch (total length when assembled)



Cannulas

- 12-3222 Bent cannula with adjustable gauge 7.75 inch 708-DS



Visit the following link to request your 850-DS Method Setup Guide and optimize your dissolution sampling today!
www.agilent.com/chem/850-ds-method-setup-guide



Basket

- 12-2100 Basket, USP, 40-mesh
- 12-2100V Basket, USP, 40-mesh with certificate



Sinkers

- 12-3000 Capsule weight wire, 316 stainless steel, 0.032 inch dia., 50 ft
- 12-3050 Capsule weight, 3-prong
- 12-3070 Sinkers basket, Japanese Pharmacopoeia (JP)

Supplies Quick Reference Guide
Download the complete list of supplies.
www.agilent.com/chem/708-ds-quick-reference-guide



Sampling Position Guidance
Optimize your 708-DS for precise sampling positioning.
<https://www.agilent.com/chem/708-ds-sampling-position-guide>



Are You Ready? Pretest Checklist
Use this quick checklist for an error-free start to every test!
www.agilent.com/chem/708-ds-pretest-checklist



Get the Answers You Need
Join the Agilent Dissolution Community for support.
<https://community.agilent.com/technical/dissolution>

