



ΔΙΑ•Μ  
судебная лаборатория

www.dia-m.ru  
analysis on-line

*Bringing power and flexibility down to size.*

**BIOMEK® NX<sup>P</sup>**

LABORATORY AUTOMATION WORKSTATION





Small  
footprint

Easy fit on  
standard  
lab benches

Many  
configurations  
to suit  
application  
needs

## Biomek® NX<sup>P</sup> workstation

By incorporating a wide range of next-generation features into a small-footprint design, Biomek NX<sup>P</sup> sets the standard for flexible laboratory solutions. It puts every aspect of liquid handling – including pipetting, dilution, dispensing and integration – into a single, automated system that is as powerful and flexible as it is efficient and economical. All while offering unprecedented flexibility for changing needs. Combined with robust Biomek Software and our ongoing applications development, Biomek NX<sup>P</sup> is bringing new possibilities to life.

## *The right features – right where you need them*

Knowing that the ability to adapt to changing situations is critical for today's labs, our systems are designed with the flexibility to accommodate whatever your needs may be. Build your Biomek system to handle today's requirements, make modifications as necessary, add accessories whenever and wherever extra capability is important, and use our wide variety of modules to manage tomorrow's needs with specifically targeted solutions. From standalone to industrial robotics solutions, it is no wonder our flexible liquid handling systems are the foundation of discovery solutions worldwide.

## *Eight-probe, 96- or 384-channel pipetting*

The NX<sup>P</sup> workstation offers multiple configurations, including multichannel, Span-8 and Span-8 with an optional gripper. The multichannel model is available with 96- and 384-channel pipetting heads for rapid plate-to-plate applications and assay miniaturization.

Interchangeable  
pipetting  
heads with  
Biomek FX<sup>®</sup>

Optional  
gripper  
mechanism  
on Span-8

Multiple  
syringe and  
probe options





### *360° rotating gripper for maximum flexibility*

With the ability to rotate a full 360°, the optional gripper handle provides access from the left, right or from integrated devices on the back of the instrument. It turns to pick up plates in different orientations, enhances on-the-fly barcode reading and saves deck space by moving titer lids while the Span-8 is pipetting.

### *Open frame design allows future integration of additional devices*

The flexible, open architecture of Biomek NX<sup>®</sup> simplifies integration, as plate readers, washers and storage devices can be located on the left, right or back of the instrument. The open frame design allows future vertical workspace integration, for applications where limited lab space is an issue.



Open base  
for vertical  
integration

Gripper  
allows device  
integration  
for Span-8  
configuration

Open platform  
for flexible  
integration

### *Pipetting performance*

Biomek NX<sup>®</sup> offers proven pipetting performance for low-volume reaction setup and assay miniaturization, as well as accurate and repeatable results extending into the submicroliter range. Different techniques can be used for various volume ranges, to enhance liquid handling capabilities. Improved X, Y and Z stability delivers reliable results for small volumes and, to ensure optimum pipetting performance, an integrated pump purges air from the lines.





Confirmatory  
barcode reads

No deck  
space needed  
for barcode  
reader

Error  
handling  
gives users  
recovery  
options

### *Built-in bar code reader delivers on-the-fly tracking*

With the ability to scan plates being held by the gripper, sample integrity is ensured. Sample IDs can be imported based on the plate barcode and tracked throughout the run. If barcode inconsistencies are detected, Biomek Software allows a wide range of user-controlled corrective actions, logging and error handling instructions.




### *Software for your needs*

First-time user of automated liquid handling technology? Experienced specialist? Either way, you'll find that Biomek® Software is packed with easy-to-use features to meet any need. And with a single, consolidated user interface driving all Biomek liquid handlers, it's simple to move between platforms. Since many functions are combined into a single step, setup is fast and effortless. And specialized algorithms speed programming even more by suggesting the best settings for the job at hand. All of which lets any user make the most of everything that Biomek NX® has to offer.

### *21 CFR Part 11 compliance*

To help achieve regulatory compliance with 21 CFR Part 11— and perform validation assays in a regulated environment—Biomek Software delivers the ability to set different user accounts and permissions, track revision control and enable electronic signatures. It also requires users to register any change to methods and stores all versions.





## Intuitive interface

The screenshot shows the 'Liquid Type Editor' dialog box. The 'Liquid Type' list on the left contains the following items: Solid, SPED, Flatten, Etched, Etch/Ink, Laser, Photo, Inkjet, Digital, Laser, Photo, Inkjet, Digital, Laser, Photo, Inkjet, Digital. The 'Solid' type is selected. The 'Solid' settings on the right include: 'Tracing for Color' (set to 'CMYK'), 'Color' (set to 'CMYK'), 'Speed' (set to '100'), 'Color Intensity' (set to '100'), 'Height' (set to '100'), 'Angle' (set to '0'), 'Speed' (set to '100'), 'Color' (set to 'CMYK'), 'Default Lines' (set to '100'), 'Default Values' (set to '100%'), 'Liquid Level Setting' (set to '100%'), 'Clear Selection' (set to '100%'), 'Clear' (set to '100%'), 'OK', 'Cancel', 'Apply'.

## Point-and-click methods development





## *The right applications – the right way*

The demands of the life science world are changing daily. So at Beckman Coulter, we're continuously delivering new ways to accelerate your valuable research. Our dedicated applications group, for instance, is focused exclusively on creating automated solutions for a variety of commercially available kits. eLabNotebook gives you quick access to a wide variety of in-house and customer-developed applications. And our renowned Biomek® Software makes it easy to develop methods to meet your individual needs. Put it all together with our robust new Biomek NX<sup>®</sup> Laboratory Automation Workstation and it is easy to see we are giving you the tools to do your best science, your way.



## *Septum piercing and clot detection*

To maintain the sample integrity so important to any application, Biomek NX<sup>®</sup> offers advanced septum piercing abilities. Clot detection circuitry ensures that samples are consistently and reliably dispensed to their proper location in a non-contaminated form – ideal for applications such as gDNA purification from whole blood.

## *Detection assays*

Cytokines have been shown to regulate immunologic responses and cell-to-cell communication, as well as dose responses to infectious agents. Their presence at elevated levels can indicate disease progress. Combined with our Biomek NX<sup>®</sup> workstation, the Immunotech Cytokine kits deliver fast and efficient methods to measure the different types of Cytokines present in a cell, while adding to the proteomics and cellular analysis solutions available from Beckman Coulter.

Use Span-8  
device to add  
reagents  
individually,  
precisely  
when needed

Optional  
integrated  
storage  
devices  
accommodate  
long  
incubation  
times





96 samples  
in less than  
one hour

DNA yield  
8-10 µg

Expandable  
system for  
higher  
throughput

### *Sample preparation*

The need for increased throughput for molecular techniques – particularly DNA/RNA purification methods – has led to a number of purification protocols on automated platforms. These automated methods, however, must generate high-quality results in order to be used by downstream applications such as sequencing, PCR\* and transfection. The Biomek® NX<sup>®</sup> provides a complete portfolio of nucleic acid sample preparation solutions for DNA and RNA with the magnetic bead-based Agencourt® SPRI® technology.

### *Automated PCR and sequencing reaction setup*

Biomek NX<sup>®</sup>'s reliable, repeatable pipetting, tube-to-plate capability and efficient automation combine to provide the ideal automation solution for this routine process, while fitting seamlessly into any lab's workflow. And its low-volume capability helps miniaturize reaction setups and save on costly reagents.



PCR &  
sequencing  
applications

Submicroliter  
transfers

Reaction setup  
in less than  
30 minutes

### *PCR and sequencing reaction cleanup*


Reaction cleanup is a simple process, to remove unwanted nucleotides and primers before running samples on a DNA analyzer, such as the CEQ™ Series from Beckman Coulter or other sequencing platforms from Applied Biosciences or GE Healthcare. To accommodate a wide variety of cleanup protocols, the Biomek FX<sup>®</sup> has been optimized for automation of paramagnetic technology. Configure a Biomek FX<sup>®</sup> with the Agencourt CleanSEQ® Dye Terminator Removal kit to process up to 17 plates per hour for your high-throughput sequencing needs.



Magbead  
technology

Validated  
methods  
for SPRI  
protocols

Up to  
17 plates  
per hour



The Biomek NX<sup>®</sup> and all our Systems Biology offerings are an important part of a broad continuum of Beckman Coulter products, including automated liquid handling, capillary electrophoresis, centrifugation, ultracentrifugation, DNA sequencing, electrochemistry, flow cytometry, fragment analysis, HPLC, integrated core systems, microarrays, particle characterization, scintillation counting, and spectrophotometry.

For information about our comprehensive line of Laboratory Automation Workstations, please contact your local Beckman Coulter representative or visit our web site at

[www.beckmancoulter.com/labautomation](http://www.beckmancoulter.com/labautomation)

\* The PCR process is covered by patents owned by Roche Molecular Systems and E. Hoffman-La Roche, Ltd. All other trademarks are the property of their respective owners.

For laboratory use only. Not for use in diagnostic procedures.

## Biomek<sup>®</sup> NX<sup>®</sup> Laboratory Automation Workstation Specifications

### Workstation

|                    |                       |
|--------------------|-----------------------|
| Weight             | 86.18 kg (190 lb.)    |
| Height             | 111.76 cm (44 in.)    |
| Width              | 91.44cm (36 in.)      |
| Depth              | 81.28 cm (32 in.)     |
| Power Requirements | 50/60 Hz, 100-240 VAC |

### Pipetting Options

96-Channel Disposable Tip Pipetting Head  
384-Channel Disposable Tip Pipetting Head  
Span-8 Fixed 100 Probes  
Span-8 Fixed 60 Probes  
Span-8 P250 Disposable Tips  
Span-8 P250 Liquid Level Sensing Tips  
Span-8 P20 Disposable Tips  
Span-8 P20 Liquid Level Sensing Tips  
Septum Piercing Probes

### Gripper Tool

Used to move labware and devices around the Biomek worksurface to automate fully filtration-based processes such as nucleic acid purification and provide access to integrated plate readers.

### Minimal Table Requirements

64 cm (25 in.) Width × 61 cm (24 in.) Depth  
Sufficient to support 136.2 kg (300 lb.)

### BIOMEK NX<sup>®</sup> PIPETTING PERFORMANCE SPECIFICATIONS

#### SPAN-8 SYSTEMS

| Transfer Volume | Span-8 Syringe Volume | Tip Types             | Accuracy ± % | Precision < % |
|-----------------|-----------------------|-----------------------|--------------|---------------|
| 0.5 µL          | 250 µL                | P20, Fixed 60 mm      | 5            | 10            |
| 1 µL            | 250, 500, 1000 µL     | P20, P50, Fixed 60 mm | 3            | 7             |
| 5 µL            | 250, 500, 1000 µL     | P20, Fixed 60 mm      | 3            | 5             |
| 10 µL           | 500 µL                | P50, P250             | 3            | 5             |
| 50 µL           | 500 µL                | P250                  | 3            | 5             |
| 100 µL          | 500 µL                | P250, P1000           | 3            | 5             |
| 250 µL          | 500 µL                | P1000                 | 2            | 2             |
| 950 µL          | 1000 µL               | P1000                 | 1            | 1             |

#### MULTICHANNEL SYSTEMS

| Transfer Volume | Multichannel Configuration | Tip Types     | Accuracy ± % | Precision < % |
|-----------------|----------------------------|---------------|--------------|---------------|
| 0.5 µL          | 384                        | P30           | 3            | 7             |
| 1 µL            | 96 & 384                   | P20, P50, P30 | 3            | 5             |
| 5 µL            | 96 & 384                   | P20 & P30     | 3            | 5             |



000 «Диаэм»

Москва  
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ [sales@dia-m.ru](mailto:sales@dia-m.ru)

[www.dia-m.ru](http://www.dia-m.ru)

С.-Петербург  
+7 (812) 372-6040  
[spb@dia-m.ru](mailto:spb@dia-m.ru)

Новосибирск  
+7 (383) 328-0048  
[nsk@dia-m.ru](mailto:nsk@dia-m.ru)

Воронеж  
+7 (473) 232-4412  
[vrn@dia-m.ru](mailto:vrn@dia-m.ru)

Йошкар-Ола  
+7 (927) 880-3676  
[nba@dia-m.ru](mailto:nba@dia-m.ru)

Красноярск  
+7 (923) 303-0152  
[krsk@dia-m.ru](mailto:krsk@dia-m.ru)

Казань  
+7 (843) 210-2080  
[kazan@dia-m.ru](mailto:kazan@dia-m.ru)

Ростов-на-Дону  
+7 (863) 303-5500  
[rnd@dia-m.ru](mailto:rnd@dia-m.ru)

Екатеринбург  
+7 (912) 658-7606  
[ekb@dia-m.ru](mailto:ekb@dia-m.ru)

Кемерово  
+7 (923) 158-6753  
[kemerovo@dia-m.ru](mailto:kemerovo@dia-m.ru)

Армения  
+7 (094) 01-0173  
[armenia@dia-m.ru](mailto:armenia@dia-m.ru)

