

COULTER® AC•T diff 2™ Hematology Analyzer

Method and Technologies

Coulter Principle: counting and sizing
Triplicate counting
Extended platelet counting
Cyanmethemoglobin

Parameters

WBC	MPV
RBC	Pct*
Hgb	PDW*
Hct	GR%
MCV	GR#
MCH	LY%
MCHC	LY#
RDW	MO%
Plt	MO#

*Not intended for diagnostic use.

Sample Analysis Characteristics

Closed or open vial sampling
Whole blood closed vial mode 18 µL
Whole blood open vial mode 18 µL
Predilute 20 µL
Automatic probe wipe

User Interface

Touch screen
Language independent icons
Single screen display of all results
Single screen for all input

Linearity

Parameter	Linearity Range	Difference
WBC	0.0-99.9 x 10 ³ cells/µL	±0.3 or ±5.0%
RBC	0.0-7.00 x 10 ⁶ cells/µL	±0.05 or ±5.0%
Hgb	0.0-25.0 g/dL	±0.2 or ±3.0%
Plt	0.0-999 x 10 ³ cells/µL	±10 or ±10.0%

Linearity, when tested using a stable sample having no interfering substances.

Precision

Parameter	Range	% CV
WBC	6.0-15.0 x 10 ³ cells/µL	≤3.0
RBC	3.00-6.00 x 10 ⁶ cells/µL	≤3.0
Hgb	12.0-18.0 g/dL	≤2.0
MCV	80.0-100.0 fL	≤3.0
Plt	200-500 x 10 ³ cells/µL	≤7.0
MPV	5.0-20 fL	≤3.0
RDW	12.0-15.0%	≤3.0
LY	20-50%	≤1.5 (S.D)
MO	2.0-10%	≤1.5 (S.D)
GR	30.0-70.0%	≤3.0 (S.D)

Imprecision is based on 31 replicate determinations of the same sample.

Operating range

Parameter	Range
WBC	0.0-150 x 10 ³ cells/µL
RBC	0.00-8.00 x 10 ⁶ cells/µL
Hgb	0.00-30.0 g/dL
MCV	50.0-130.0 fL
Plt	0.00-3000 x 10 ³ cells/µL
MPV	5.0-20.0 fL
LY%,	MO%, GR% 0-100%
LY#,	MO#, GR# 0-99.9 x10 ³ cells/µL

Throughput

Results in less than 60 seconds.
Minimum of 50 samples per hour.

Data Management

Patient ID auto-numbering
14 digit manual patient ID
Three user-definable patient ranges for reporting results
Linearity range for result limits reporting
Data storage for 250 patient results
QC storage for 3 controls and 279 results
Levy Jennings graphs

Printer

Graphics printer, roll printer or ticket printer

Host Communication

RS232C with ASTM standards

Dimensions/Weight

Height	49cm (19in)
Width	45cm (18in)
Depth	43cm (17in)
Weight	20kg (45lb)

Power

100v + 10%, 50/60Hz
120v + 10%, 50/60Hz
220v + 10%, 50/60Hz
240v + 10%, 50/60Hz

Consumption

Less than 120W

Installation

Category II per IEC 1010-1

Results You Can Trust with the Superior Support Your Lab Needs.

COULTER® AC•T diff 2™ Hematology Analyzer



- General Chemistry
- Immunodiagnosics
- Centrifugation
- Disease Management
- Clinical Research
- Hematology
- Hemostasis
- Lab Automation
- Data Management
- Flow Cytometry
- Primary Care



SIMPLIFY • AUTOMATE • INNOVATE



Africa/Middle East/Eastern Europe: Switzerland, Nyon (41) 22 994 0707. Australia, Gladesville (61) 2 9844 6000. Canada, Mississauga (1) 905 819 1234. China, Beijing (86) 10 6515 6028. Hong Kong (852) 2814 7431, 2814 0481. France, Villepinte (33) 1 49 90 90 00. Germany, Krefeld (49) 2151 33 35. Italy, Milan (39) 02 953921. Japan, Tokyo (81) 3 5404 8424. Latin America (1) (305) 380 4709. Mexico, Mexico City (525) 575 6805. Netherlands, Mijdrecht (31) 297 230630. Singapore (65) 339 3633. South Africa, Johannesburg (27) 11 805 2014. Sweden, Bromma (46) 8 564 85 900. Switzerland, Nyon 0800 850 810. Taiwan, Taipei (886) 2 2378 3456. Turkey, Istanbul (90) 216 309 1900. UK, High Wycombe (44) 01494 441181. USA, Brea, CA (1) 800 352 3433, (1) 714 993 5321.



The AcT diff 2 Delivers Reliable Results at an Incredible Value



Gain confidence in results with the AcT diff 2.

The COULTER® AC•T diff 2™ hematology analyzer is an affordable, easy-to-use system that is ideal for a physician's office lab or smaller laboratory. As a fully automated, 16 or 18 parameter system, the AC•T diff 2 offers a CBC

and Coulter histogram differential for consistently accurate results. Plus, this compact analyzer provides closed vial sampling to maximize safety and efficiency in your lab. It also offers a small sample volume in both open and closed analysis modes, a wide operating range and multiple output formats. With the reliable performance of the AC•T diff 2 and our superior customer support, you'll achieve maximum uptime in your laboratory.



Specially formulated reagents, controls and calibrators ensure accurate results.

Efficient Data Management Ensures Accurate Results

- User-definable high/low flagging of patient results and quality control simplifies data interpretation
- Storage of three levels of control and automated calibration procedures enhances productivity and assures accurate capture and calculation of data
- Meet your lab's individual needs with six different reporting styles including WBC, RBC & PLT histograms and tickets

Ease of Operation and Superior Sample Handling Increases Productivity

- CBC and differential results are easily obtained with sophisticated, yet user-friendly technology
- Wide operating ranges and extended platelet counting simplify handling of samples with abnormally high or low results
- Closed vial sampling for controls and calibrators provides operator safety
- Small 18 µl sample size enables processing of hard-to-obtain specimens

Service You Can Count On

We continually prove our dedication to our customers through our world-class service and support network. We provide fast, efficient technical support with a quick service response time. Regardless of which Beckman Coulter systems your lab relies on, you can count on us to ensure that they will continue to operate at peak efficiency for years to come. Contact your Beckman Coulter representative today for more information on how the AC•T diff 2 can deliver greater productivity and efficiency in your lab.



Running Samples is Easy

1. Select sample mode and verify ID via the intuitive touch screen



2. Place sample or control tube into the closed vial sample station



3. Review results in less than 60 seconds.



Superior service and a quick response time increase uptime in your lab.