

Stoelting

EVOLVING NEUROSCIENCE AND BEHAVIOR RESEARCH

- Stereotaxic Instruments 6
- Surgical Instruments50
- Anesthesia34
-  **ANY-maze**108
-  **Ugo Basile**136
- Behavior Testing Equipment
...and much more!



Company & Order Information

In need of assistance? We can help!

Contact Stoelting for all your scientific needs.



U.S. & Canada Order Information:

Web: www.StoeltingCo.com
 Email: Neuro@StoeltingCo.com
 Call: 1-800-860-9775 or 630-860-9700
 Fax: 630-860-9775
 Mail: Stoelting Co • 620 Wheat Lane
 Wood Dale, IL 60191 • USA



Europe Order Information:

Web: www.stoeltingeurope.com
 Email: info@stoeltingeurope.com
 Call: +353 1 524 2200
 Mail: Stoelting Europe • Landscape House • Landscape Rd
 Churchtown, Dublin • D14A6P3 • Ireland



India Order Information:

Web: www.StoeltingCo.com
 Email: gentechmd@yahoo.com
 Call: +91-11-47563575 or +91-11-47472393
 Mail: Gentech Marketing & Distribution Pvt Ltd
 B116 & B117, R G Complex, D B Gupta Road
 Motia Khan • New Delhi 110055



CONNECT WITH US

ALL INTERNATIONAL CUSTOMERS

(Outside of North America, Europe or India) Please contact Stoelting's North American office for assistance.

SALES INQUIRIES

Please contact our Sales Department at Stoelting with Neuroscience or ANY-maze related sales inquiries.

TECHNICAL ASSISTANCE

Neuroscience

Please contact our Neuroscience Technical Staff at Neuro@StoeltingCo.com for technical assistance.

ANY-maze

Please contact ANY-maze Support at Support@ANY-maze.com or TechSupport@ANY-maze.com for technical assistance.

ORDERS

Orders may be placed by phone, fax, e-mail (orders@StoeltingCo.com) or through our website at www.StoeltingCo.com. For orders placed in the U.S., Canada and Mexico, please address your order to Stoelting Headquarters located in Wood Dale, IL.

For orders placed in Europe, please address your order to Stoelting Europe.

For orders placed in India, please address your order to Stoelting India.

All International orders (outside of North America, Europe and India) please contact our North American office for further assistance. Payments can be made using a credit card (Visa, Mastercard, American Express), bank transfer or purchase order (PO; pre-approved accounts only).

RETURNS AND REPAIRS

Please contact our Returns Department at Stoelting for assistance with your return/repair. All returned items must be assigned a Return Materials Authorization (RMA) and be accompanied by a completed Equipment Decontamination form. Any merchandise that has not been assigned an RMA or does not include the Equipment Decontamination form will be returned to the sender. A restocking fee may be applied to all returned merchandise at the discretion of Stoelting Co. Customized equipment is not eligible for return. Following inspection, you will be contacted with a repair estimate prior to initiation of the repair. Repair costs can be paid by VISA, MASTERCARD, AMEX, bank transfer or purchase order (PO; pre-approved accounts only).

WARRANTY INFORMATION

Stoelting Co. warrants its products against defects in materials and workmanship for a period of one-year from the date of shipment. For warranty repairs, please include a copy of your original invoice.

**Prices and product availability subject to change without notice. Call or visit www.StoeltingCo.com for current prices.*

We thank you for your interest in Stoelting products, and look forward to serving you today, and in the years to come.

Rat and Mouse Behavior Mazes



60135
Morris Water Maze

Water Maze

The Morris Water Maze is one of the most widely used methods in behavioral neuroscience to assess spatial and learning memory. The Water Maze tank is filled with water and equipped with a platform (60030) which acts as a means of escape for animals. During testing, animals learn the location of the platform using spatial cues. The Water Maze tank is made of aluminum and is available in three different sizes for testing of Mice, Rats or both. Platform (60030) is included.

WATER MAZE

60135	Water Maze, Mouse, 3.5-4.5 ft Diameter (1.07-1.37m)
60136	Water Maze, Mouse/Rat, 4.5-5.5 ft Diameter (1.37-1.68m)
60235	Water Maze, Rat, 5.5-6.5 ft Diameter (1.68-2.0m)

AVAILABLE ISLANDS

60030	4 in (10.2 cm) Diameter/12 in (30.5cm) Height
60035	4 in (10.2 cm) Diameter/20 in (50.8cm) Height

Barnes Maze

The Barnes Maze is designed for testing spatial learning and memory. Within the maze, animals are motivated to escape from the brightly lit, open platform to a small dark recessed chamber (escape box). During testing, the animal learns the spatial location of the escape box. The maze also includes false escape boxes, which help to remove visual cues that might be observed from a distance or through an open hole.



60170 Barnes Maze

CATALOG NUMBER		TABLE DIA.	HOLE DIA.	HEIGHT	PRICE
60170	Barnes Maze, Mouse	91cm	5cm	90cm	
60270	Barnes Maze, Rat	122cm	10cm	90cm	

Save money when purchasing an ANY-maze bundle!

Maze bundle packages include: Laptop computer, USB Camera, Camera Accessories, ANY-maze software and a Stoelting maze of your choice! See page 131 for more details.

