

Stoelting

EVOLVING NEUROSCIENCE AND BEHAVIOR RESEARCH

- Stereotaxic Instruments 6
- Surgical Instruments50
- Anesthesia34
-  **ANY-maze**108
-  **Ugo Basile**136
- Behavior Testing Equipment
...and much more!



Company & Order Information

In need of assistance? We can help!

Contact Stoelting for all your scientific needs.



U.S. & Canada Order Information:

Web: www.StoeltingCo.com
 Email: Neuro@StoeltingCo.com
 Call: 1-800-860-9775 or 630-860-9700
 Fax: 630-860-9775
 Mail: Stoelting Co • 620 Wheat Lane
 Wood Dale, IL 60191 • USA



Europe Order Information:

Web: www.stoeltingeurope.com
 Email: info@stoeltingeurope.com
 Call: +353 1 524 2200
 Mail: Stoelting Europe • Landscape House • Landscape Rd
 Churchtown, Dublin • D14A6P3 • Ireland



India Order Information:

Web: www.StoeltingCo.com
 Email: gentechmd@yahoo.com
 Call: +91-11-47563575 or +91-11-47472393
 Mail: Gentech Marketing & Distribution Pvt Ltd
 B116 & B117, R G Complex, D B Gupta Road
 Motia Khan • New Delhi 110055



CONNECT WITH US

ALL INTERNATIONAL CUSTOMERS

(Outside of North America, Europe or India) Please contact Stoelting's North American office for assistance.

SALES INQUIRIES

Please contact our Sales Department at Stoelting with Neuroscience or ANY-maze related sales inquiries.

TECHNICAL ASSISTANCE

Neuroscience

Please contact our Neuroscience Technical Staff at Neuro@StoeltingCo.com for technical assistance.

ANY-maze

Please contact ANY-maze Support at Support@ANY-maze.com or TechSupport@ANY-maze.com for technical assistance.

ORDERS

Orders may be placed by phone, fax, e-mail (orders@StoeltingCo.com) or through our website at www.StoeltingCo.com. For orders placed in the U.S., Canada and Mexico, please address your order to Stoelting Headquarters located in Wood Dale, IL.

For orders placed in Europe, please address your order to Stoelting Europe.

For orders placed in India, please address your order to Stoelting India.

All International orders (outside of North America, Europe and India) please contact our North American office for further assistance. Payments can be made using a credit card (Visa, Mastercard, American Express), bank transfer or purchase order (PO; pre-approved accounts only).

RETURNS AND REPAIRS

Please contact our Returns Department at Stoelting for assistance with your return/repair. All returned items must be assigned a Return Materials Authorization (RMA) and be accompanied by a completed Equipment Decontamination form. Any merchandise that has not been assigned an RMA or does not include the Equipment Decontamination form will be returned to the sender. A restocking fee may be applied to all returned merchandise at the discretion of Stoelting Co. Customized equipment is not eligible for return. Following inspection, you will be contacted with a repair estimate prior to initiation of the repair. Repair costs can be paid by VISA, MASTERCARD, AMEX, bank transfer or purchase order (PO; pre-approved accounts only).

WARRANTY INFORMATION

Stoelting Co. warrants its products against defects in materials and workmanship for a period of one-year from the date of shipment. For warranty repairs, please include a copy of your original invoice.

**Prices and product availability subject to change without notice. Call or visit www.StoeltingCo.com for current prices.*

We thank you for your interest in Stoelting products, and look forward to serving you today, and in the years to come.

Rat and Mouse Behavior Mazes



60450 Sociability Cage

Sociability Apparatus

Ideal for the investigation of:

- Autism
- Maternal Behavior
- Dominance
- Social Memory
- Sociability

The Sociability Apparatus is used to assess social interaction and conditions thereof. The three- chambered apparatus includes two gray “stranger” enclosures with wide grid bars to promote interaction, and a non-reflective gray floor to maximize contrast, for video tracking purposes.

Rat sizes available.



SOCIABILITY APPARATUS

60450	Sociability Cage, Mouse
60449XL	Sociability Cage, Rat

66000 Repeated Acquisition and Performance Chamber



Repeated Acquisition and Performance Chamber (RAPC)

RAPC is a behavioral monitoring apparatus designed to completely automate and assess spatial learning and memory. RAPC enhances the accuracy of recording temporal measurements, and error scoring, during spatial learning and memory experiments, through the use of the ANY-maze data acquisition software and electronic sensors. In fact RAPC is an ANY-maze software driven product. The software can be used to track as many as 16 chambers simultaneously, and performs full statistical analysis of results immediately upon completion of the experiment.

- Automated data collection results in high throughput and data collection accuracy
- Acrylic construction allows for easy animal observation
- ANY-maze Software control provides ease of operation
- A single personal computer can collect data from up to 16 chambers simultaneously
- ANY-maze software automatically configures the experiments
- ANY-maze controls the experiment flow, data collection and statistical analysis
- The chamber is equipped with two Mylar speakers to generate white noise or provide auditory cues

**Automated
assessment of
spatial learning
and memory!**

STOELTING REPEATED ACQUISITION AND PERFORMANCE CHAMBER (RAPC)

66000	RAPC Repeated Acquisition and Performance Chamber
--------------	---