

Stoelting

EVOLVING NEUROSCIENCE AND BEHAVIOR RESEARCH

- Stereotaxic Instruments 6
- Surgical Instruments50
- Anesthesia34
-  **ANY-maze**108
-  **Ugo Basile**136
- Behavior Testing Equipment
...and much more!



Company & Order Information

In need of assistance? We can help!

Contact Stoelting for all your scientific needs.



U.S. & Canada Order Information:

Web: www.StoeltingCo.com
 Email: Neuro@StoeltingCo.com
 Call: 1-800-860-9775 or 630-860-9700
 Fax: 630-860-9775
 Mail: Stoelting Co • 620 Wheat Lane
 Wood Dale, IL 60191 • USA



Europe Order Information:

Web: www.stoeltingeurope.com
 Email: info@stoeltingeurope.com
 Call: +353 1 524 2200
 Mail: Stoelting Europe • Landscape House • Landscape Rd
 Churchtown, Dublin • D14A6P3 • Ireland



India Order Information:

Web: www.StoeltingCo.com
 Email: gentechmd@yahoo.com
 Call: +91-11-47563575 or +91-11-47472393
 Mail: Gentech Marketing & Distribution Pvt Ltd
 B116 & B117, R G Complex, D B Gupta Road
 Motia Khan • New Delhi 110055



CONNECT WITH US

ALL INTERNATIONAL CUSTOMERS

(Outside of North America, Europe or India) Please contact Stoelting's North American office for assistance.

SALES INQUIRIES

Please contact our Sales Department at Stoelting with Neuroscience or ANY-maze related sales inquiries.

TECHNICAL ASSISTANCE

Neuroscience

Please contact our Neuroscience Technical Staff at Neuro@StoeltingCo.com for technical assistance.

ANY-maze

Please contact ANY-maze Support at Support@ANY-maze.com or TechSupport@ANY-maze.com for technical assistance.

ORDERS

Orders may be placed by phone, fax, e-mail (orders@StoeltingCo.com) or through our website at www.StoeltingCo.com. For orders placed in the U.S., Canada and Mexico, please address your order to Stoelting Headquarters located in Wood Dale, IL.

For orders placed in Europe, please address your order to Stoelting Europe.

For orders placed in India, please address your order to Stoelting India.

All International orders (outside of North America, Europe and India) please contact our North American office for further assistance. Payments can be made using a credit card (Visa, Mastercard, American Express), bank transfer or purchase order (PO; pre-approved accounts only).

RETURNS AND REPAIRS

Please contact our Returns Department at Stoelting for assistance with your return/repair. All returned items must be assigned a Return Materials Authorization (RMA) and be accompanied by a completed Equipment Decontamination form. Any merchandise that has not been assigned an RMA or does not include the Equipment Decontamination form will be returned to the sender. A restocking fee may be applied to all returned merchandise at the discretion of Stoelting Co. Customized equipment is not eligible for return. Following inspection, you will be contacted with a repair estimate prior to initiation of the repair. Repair costs can be paid by VISA, MASTERCARD, AMEX, bank transfer or purchase order (PO; pre-approved accounts only).

WARRANTY INFORMATION

Stoelting Co. warrants its products against defects in materials and workmanship for a period of one-year from the date of shipment. For warranty repairs, please include a copy of your original invoice.

**Prices and product availability subject to change without notice. Call or visit www.StoeltingCo.com for current prices.*

We thank you for your interest in Stoelting products, and look forward to serving you today, and in the years to come.

Rat and Mouse Behavior Mazes

Open arms
are equipped
with height-
adjustable lip to
minimize falls.



60140
Elevated Plus-Maze



Close-Up of
Elevated Plus-Maze's
Adjustable Rim

Elevated Plus-Maze

The Elevated Plus-Maze is used to assess anxiety in rodents. The maze exploits the conflict between the innate fear that rodents have of open areas, versus their desire to explore novel environments. The Elevated Plus-Maze is designed with an elevated plus-shaped platform, with two closed arms and two open. Both open arms are equipped with height-adjustable lip to minimize falls.

CATALOG NUMBER		LANE WIDTH	ARM LENGTH	WALL HEIGHT	LEG HEIGHT	PRICE
60140	Plus-Maze, Mouse	5cm	35cm	15cm	40cm	
60240	Plus-Maze, Rat	10cm	50cm	40cm	40cm	

Elevated Zero-Maze

The Elevated Zero-Maze is a modification of the Elevated Plus-Maze, which incorporates both traditional and novel ethological measures, in the analysis of anxiety. The design comprises an elevated circular platform, with two opposite enclosed quadrants, and two open. This removes any ambiguity in interpretation of time spent on the central square of the traditional "plus" design, and allows uninterrupted exploration.



60190
Elevated
Zero-Maze

Includes
adjustable rim
similar to
the elevated
Plus-Maze.

CATALOG NUMBER		DIA.	LANE WIDTH	WALL HEIGHT	LEG HEIGHT	PRICE
60190	Zero-Maze, Mouse	50cm	5cm	15cm	40cm	
60290	Zero-Maze, Rat	100cm	10cm	30cm	40cm	