

Stoelting

EVOLVING NEUROSCIENCE AND BEHAVIOR RESEARCH

- Stereotaxic Instruments 6
- Surgical Instruments50
- Anesthesia34
-  **ANY-maze**108
-  **Ugo Basile**136
- Behavior Testing Equipment
...and much more!



Company & Order Information

In need of assistance? We can help!

Contact Stoelting for all your scientific needs.



U.S. & Canada Order Information:

Web: www.StoeltingCo.com
 Email: Neuro@StoeltingCo.com
 Call: 1-800-860-9775 or 630-860-9700
 Fax: 630-860-9775
 Mail: Stoelting Co • 620 Wheat Lane
 Wood Dale, IL 60191 • USA



Europe Order Information:

Web: www.stoeltingeurope.com
 Email: info@stoeltingeurope.com
 Call: +353 1 524 2200
 Mail: Stoelting Europe • Landscape House • Landscape Rd
 Churchtown, Dublin • D14A6P3 • Ireland



India Order Information:

Web: www.StoeltingCo.com
 Email: gentechmd@yahoo.com
 Call: +91-11-47563575 or +91-11-47472393
 Mail: Gentech Marketing & Distribution Pvt Ltd
 B116 & B117, R G Complex, D B Gupta Road
 Motia Khan • New Delhi 110055



CONNECT WITH US

ALL INTERNATIONAL CUSTOMERS

(Outside of North America, Europe or India) Please contact Stoelting's North American office for assistance.

SALES INQUIRIES

Please contact our Sales Department at Stoelting with Neuroscience or ANY-maze related sales inquiries.

TECHNICAL ASSISTANCE

Neuroscience

Please contact our Neuroscience Technical Staff at Neuro@StoeltingCo.com for technical assistance.

ANY-maze

Please contact ANY-maze Support at Support@ANY-maze.com or TechSupport@ANY-maze.com for technical assistance.

ORDERS

Orders may be placed by phone, fax, e-mail (orders@StoeltingCo.com) or through our website at www.StoeltingCo.com. For orders placed in the U.S., Canada and Mexico, please address your order to Stoelting Headquarters located in Wood Dale, IL.

For orders placed in Europe, please address your order to Stoelting Europe.

For orders placed in India, please address your order to Stoelting India.

All International orders (outside of North America, Europe and India) please contact our North American office for further assistance. Payments can be made using a credit card (Visa, Mastercard, American Express), bank transfer or purchase order (PO; pre-approved accounts only).

RETURNS AND REPAIRS

Please contact our Returns Department at Stoelting for assistance with your return/repair. All returned items must be assigned a Return Materials Authorization (RMA) and be accompanied by a completed Equipment Decontamination form. Any merchandise that has not been assigned an RMA or does not include the Equipment Decontamination form will be returned to the sender. A restocking fee may be applied to all returned merchandise at the discretion of Stoelting Co. Customized equipment is not eligible for return. Following inspection, you will be contacted with a repair estimate prior to initiation of the repair. Repair costs can be paid by VISA, MASTERCARD, AMEX, bank transfer or purchase order (PO; pre-approved accounts only).

WARRANTY INFORMATION

Stoelting Co. warrants its products against defects in materials and workmanship for a period of one-year from the date of shipment. For warranty repairs, please include a copy of your original invoice.

**Prices and product availability subject to change without notice. Call or visit www.StoeltingCo.com for current prices.*

We thank you for your interest in Stoelting products, and look forward to serving you today, and in the years to come.

Mouse Behavior Mazes



60100 Open Field

Open Field

The Open Field test is commonly used to assess locomotor activity, exploratory behavior and anxiety in rodents. Open field apparatuses are available with a 40cm square (clear, gray or IR black) or round (clear or gray) enclosure, and gray base plate to maximize contrast, for video tracking.



60101



60102

CATALOG NUMBER		WIDTH	DEPTH	HEIGHT	PRICE
60100	Open Field, Mouse, Clear Walls	40cm	40cm	35cm	
60101	Open Field, Mouse, Gray Walls	40cm	40cm	35cm	
60102	Open Field, Mouse, IR Black* *For use with photobeam arrays See page 116	40cm	40cm	35cm	

Visit our website
for more options!

Novel Object Recognition Test

The Novel Object Recognition Test is based on the premise that rodents will explore a novel object, more than a familiar one, if the animal is able to remember the familiar one. The flexibility of the Novel Object Recognition Test allows for testing of short- or long-term memory, and can be used to selectively test the effects of an acute drug treatment on a specific stage of memory formation.

CATALOG NUMBER	ENCLOSURE			PRICE
	WIDTH	DEPTH	HEIGHT	
62017 Novel Object Test, Mouse	40cm	40cm	35cm	
62007 Novel Object Accessories Only (see page 129)	Set of 12 (Includes 2 ea. square and round blocks in gray, black and white)			



62017 Novel Object Test